## ENGLISH

- Finish and polish.
- Try fit cured Protemp Crown and make necessary adjustments.






- Check occlusion.
  - Use a standard method to finish and polish.
  - Finish using a fine carbide bur and/or 3M ESPE Sof-Lex™ Contouring and Polishing Discs.
- Polish using a dry muslin rag wheel or polishing brush.

- Check all adjusted margins of the finished and polished crown prior to cementation, and make further adjustments if necessary.
- Check all adjusted margins of the finished and polished crown prior to cementation, and make further adjustments if necessary.

Protemp Crown temporization material can be used to model pre-existing tooth structure and anatomy for use in digital imaging applications.

**Note:**

- Proper adaptation of the margin is needed for optimal marginal integrity. Adaptation and smooth finishing of the crown margin (facial, lingual, interproximal) are needed to reduce the likelihood of leakage and gingival tissue irritation.*
- A professional assessment of oral health, chewing habits, and patient concerns should be factored into the decision.*

Size Selection Chart		
<b>Crown Name</b>	<b>Approximate Range Mesial Distal Width</b>	
<b>BICUSPID (Premolar) Upper</b>		
Bicuspid (Premolar) Upper S	6 - 7 mm	
Bicuspid (Premolar) Upper L	7 - 8 mm	
<b>MOLAR (1st &amp; 2nd) Upper</b>		
Molar Upper S	9 - 11 mm	
Molar Upper L	11 - 12 mm	
<b>CUSPID (Canine) Universal</b>		
Cuspid (Canine) Universal S	6 - 8 mm	
Cuspid (Canine) Universal L	8 - 9 mm	
<b>BICUSPID (Premolar) Lower</b>		
Bicuspid (Premolar) Lower	6 - 9 mm	
<b>MOLARS (1st &amp; 2nd) Lower</b>		
Molar Lower S	9 - 11 mm	
Molar Lower L	11 - 12 mm	

**Abutment Preparation**

Prepare abutment using conventional techniques. (Note: adequate abutment structure is necessary to retain crown). In case of small abutment dimensions (e.g. narrow cast and core), relining is recommended to obtain optimal fit and retention.

**Precautionary Information for Patients:**

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known acrylate allergies. If allergic reaction occurs, seek medical attention as needed, remove the product if necessary and discontinue future use of the product.

**Precautionary information for Dental Personnel:**

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If contact occurs directly with the skin, wash skin with soap and water. Use of protective gloves is recommended; however, acrylates may penetrate commonly used gloves. After product placement, remove and discard gloves, wash hands immediately with soap and water and then re-glove. If allergic reaction occurs, seek medical attention as needed.

3M ESPE MSDSs can be obtained from [www.3MESPE.com/MSDS](http://www.3MESPE.com/MSDS) or contact local subsidiary.

**Instructions For Use**

### 1) Size Selection

**Determine proper size crown.** After tooth is prepared, determine the correct Protemp Crown size.

Measure the mesial-distal width of the prepared tooth using the 3M ESPE Crown Size Tool or a caliper.

**3M ESPE Crown Size Tool.**

*The Crown Size Tool contains graduated widths from 7-12 millimeters. Insert the tool until contact is made with the adjacent teeth and in line with contact points and resting on the occlusal surface of the crown preparation. Slide tool between adjacent teeth (the tool will wedge). The wedging point, relative to the tool establishes the approximate mesial distal width of the Protemp Crown. The 3M ESPE Crown Size Tool is designed for single use.*

### Select crown.

- Select the Protemp Crown that coincides with the approximate mesial distal width range provided on the Size Selection Chart.

*When in doubt about the size, use the smaller sized crown.*

- Remove the Protemp Crown from the sealed crown case.

- Hold film between thumb and finger and carefully remove film from Protemp Crown.

**Note:**

- Measure carefully, once removed from sealed package, the crown cannot be placed back into the package for re-use at a later date.*

- Protemp Crowns are light sensitive and preparation should be taken to minimize ambient and operatory light exposure.*

### 2) Protemp Crown Preparation

- Gently roll crown in gloved fingers to pre-warm and straighten axial walls.

**Note:**

- If the Protemp Crowns are stored in a cool environment, gently rolling the axial wall of the crown will make it more malleable and less prone to fracture during placement.*

*Gently rolling eliminates the flare at the margin and reduces the circumference of the crown at the margin.*

*As the crown is placed onto the tooth preparation it will expand and hug the preparation.*

### Trim Protemp Crown.

• Determine amount of excess to be removed from the Protemp Crown by measuring the approximate height of the adjacent teeth. This can be done by using the 3M ESPE Crown Size Tool, perio probe, or by placing the Protemp Crown on the preparation.

- Follow the gingival contour when trimming excess material to establish proper crown length using crown scissors.
- Make sure to trim enough material. If too little material is trimmed before checking occlusion, the occlusal surface of the crown can be damaged which will lead to fracturing when the crown is cemented.*

### 3) Adaptation

#### Adapt Protemp Crown.

- Seat trimmed Protemp Crown to verify length.
- Adapt, shape and establish interproximal contacts and preparation margins for a snug fit onto the moist preparation.
- Have patient gently close to establish occlusion, open to adapt occlusal surface.

#### Adapt facial margin.

- While patient is in occlusion adapt the facially accessible margin.

**Note:**

- After the crown is adapted to the prepared tooth a tack-cure is required to reduce the possibility of distortion upon removal. Do Not Over Tack Cure.*
- Track cure the facial surface using a dental curing light for 2-3 seconds.

#### Adapt lingual margin.

- Have patient open, holding crown in place.
- Adapt lingual margin. While adapting the lingual margin apply finger pressure to the adapted facial surface to prevent dislodging.
- Track cure the lingual surface for 2-3 seconds.
- Track cure the occlusal surface for 2-3 seconds.

*After tack curing, remove the crown and reseat several times to check fit before final cure.*

**Note:**

- Dental curing light with a wavelength range 420-500 nm and a minimum light intensity of 400 mW/cm<sup>2</sup> is recommended.*
- Avoid looking directly at the light guide and do not stare at the light reflection from the crown surfaces while curing.*

#### 4) Final Light Cure

- Carefully remove the Protemp Crown.
- Light cure adapted Protemp Crown outside of the mouth for a total of 60 seconds.
- Fully light cure by exposing all surfaces of the Protemp Crown to the light.

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#### Place crown.

- Use gentle thumb pressure, slowly seat crown into place.
- Have patient gently close into occlusion and allow excess cement to flow out.

#### Remove excess cement.

- Follow cement manufacturer's instructions for proper set time.
- Carefully check the sulcus of the treated teeth and the surrounding areas: remove any cement residue that remains in the sulcus and along the Protemp Crown margin.

- Clean the interproximal areas by drawing a piece of knotted floss back and forth under the contact points.

**Note:** It is important to follow the instructions for use very closely. In addition, for optimal performance:

- When in doubt, select the smaller crown size*
- Make sure to trim enough material.*

*If too little material is trimmed before checking the occlusion, the occlusal surface of the crown may be damaged which can lead to fracturing upon cementation.*

*After tack curing, remove the crown and reseat several times to check the fit before final cure.*

### Protemp Crown for Digital Imaging

In its uncured state, Protemp Crown temporization material is soft and moldable. In this malleable state, the crown can be shaped to model pre-existing tooth structure and anatomy. Once the desired shape is achieved, the crown is fully cured for 60 seconds. This crown can then be placed and scanned to produce a digital image.

**Custom Characterization and Repair**
Use a flowable; for example 3M ESPE Filtek™ Supreme Plus Flowable, or ProTemp™ Plus for custom characterization. Custom characterization using either of these materials does not require the use of a bonding agent. If you intend to use other products, refer to the pertinent manufacturer's instructions for use.

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- Thoroughly dry with a stream of air.
- Apply the composite materials in increments to the Protemp Crown.

- Follow recommended light cure times for composite materials.
- Finishing and polish as recommended in the "Finish and Polish" instructions.

The finished crown should approximate, as nearly as possible, the original tooth in terms of space occupied in the arch, proper relationship with adjacent teeth and functional occlusion with the opposing dentition.

### Storage and Use

This product is designed to be stored and used at room temperature. Shelf life at room temperature is 24 months. Ambient temperatures routinely higher than 27°C/80°F or lower than 10°C/50°F may reduce shelf life. Do not store materials in proximity to eugenol containing products. See outer package for expiration date.

#### Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

#### Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for your's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

#### Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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