**D\_** Devices

Advanced systems for dental devices

**D\_SMALL STEAM STERILIZERS CAL23**

**EN\_** User manual

**ES\_** Manual de instrucciones

## Introduction

- Before operating the sterilizer, please read carefully all safetycautions and instructions of operation. This operator's manual will help you understand all functions of Small Steam Sterilizers as far as possible.
- Please follow carefully instructions in this operator's manual While servicing and maintaining the units.
- Please keep this manual safe for your future references.
- When error occurs during operation of the machine, please get in touch with your local distributor or us for our best qualified services and assistance.

## Contents

1. Important information on Safety.....	4
2. Description & Usage.....	4
3. Specifications .....	5
4. Installation .....	6
5. Control panel .....	7
6. Operation.....	8-10
7. Door Adjustment.....	11
8. Printer installation and use .....	12-13
9. System setup.....	14
10. Maintenance & Maintain .....	14-17
11. Transport & Storage.....	18
12. Warranty .....	18
13. Accessory.....	18
14. Appendixes.....	19-22
15. Troubleshooting .....	23
16. Important information .....	23

## 1. Important information on Safety

- Please plug the power in a correct voltage (according to nameplate) electrical socket. Don't use electrical outlets other than the correct voltage; otherwise it might cause fire or electric shock.
- Don't wet your hand before plugging or unplugging the power cord.
- Don't place your sterilizer on an unstable table top, moving desk, tilted and shaking surface
- Please don't cover or block the door, ventilation or radiation openings on the sterilizer.
- Please don't put heavy stuff or liquid containers that might spill on the top of the machine.
- Please unplug the power cord when the machine isn't used for a long while or when you travel etc.

## 2. Description & Usage

D\_Devices design and make pressurized steam Small Steam Sterilizer especially to meet the needs for sterilization of instruments and devices in clinics, departments in a hospital and laboratories. Operators include doctors or his/her assistance under doctors supervisions. The Small Steam Sterilizer is computerized, user friendly and easy to operate. Its real time display shows all data during the sterilization cycle. The software program provides auto diagnose, protection against over heating or over pressurizing, and thus ensures reliability of the sterilization results. Build- in condenser water collecting system prevents discharge of waste steam and keeps the machine clean and safe.

### OU-22L, OU-17L Steam sterilizers

OU uses the advanced technology of 3 times pre -vacuum and vacuum drying. It enables the hot steam to penetrate in thin tube and porous material. This product is suitable for dentistry, ophthalmology, surgical surgical instruments, dressing theatres and laboratories to sterilize dental and and other 135 °C autoclavable stuffs.

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<b>Chamber inner dimensions (Diameter X Depth)</b>	22 L: 249 mm x 450 mm
	17 L: 249 mm x 355 mm
<b>Power supply</b>	According to the nameplate
<b>Sterilizing temperature/pressure</b>	121°C/1,0~1,3 bar (under standard atmospheric pressure)
	134°C/2,0~2,4 bar (under standard atmospheric pressure)
 Caution: Temperature/Pressure here means the chamber temperature and pressure during the phase of sterilization. The value varies under different atmospheric pressure.	
<b>Power fuse</b>	F20AL 250 V/F12AL 250 V
<b>Water tank capacity</b>	3.5 L
<b>Environment temperature</b>	+5°C ~ +40°C
<b>Relative humidity range</b>	30% - 75%

#### Figures

Package Key

Test Cycle Key

Caution & Warning

Start/Stop key

Used Water Outlet

Type B Device

Select Key

Cycle Step Key

Grounding


Dry Key

Tip: Standard atmospheric pressure: It the atmospheric pressure at sea level. It refers to 1 Standard atmospheric pressure. The value is 0.1Mpa=1000Hpa=1.0 Bar. Atmospheric pressure decreases by higher altitude. An increase of altitude of 1000 meters will result in a decrease of ATM PR of approximately BAF decreases by about

#### 4. Installation

The The Small Steam Sterilizer should be installed in a drafty place allowing minimum of 10 cm space at all sides and 50 cm at the top. The ventilation openings at the side of the machine should not be blocked. The The Small Steam Sterilizer should be placed on a leveled counter top. The dimension of machine as the figure:


**⚠ Caution:** If the front of the machine is lower than the back the control will display "E3" (E3 means it does not dry completely).

#### Preparation before using

We installed the draining filter before the device leaving factory. Clean the draining filter once per month. You can see the instructions of how to clean it according to page 15.

Before using the device, you should make electrical connections by using the power cord that comes with the unit. Insert the hollowed end of the power cord in to the inlet of the machine and plug the other end on a power outlet. **Caution:** The power supply can afford 1800VA load. Ok, you can turn on the The Small Steam Sterilizer now.

**⚠ Caution:** The power supply must be can afford (22L) 2300VA or (17L) 1800VA.

**Tips:** Before starting the Small Steam Sterilizer, please check on the local atmospheric pressure. If the pressure is lower than 0.095 Mpa or the Altitude is higher than 500 meters, you will need to contact your distributor who should reconfigure the controller in the device.


## 5. Control panel

1. **Pressure:** It displays actual pressure of chamber, the unit is bar.
2. **Temperature:** It displays actual temperature of chamber.
3. **State displays / Error code displays:** It displays all the states of cycle (Look at the appendix 1). When it occurs error, it will show the error code, just according to error code mode.

 Caution: If you meet problems, please contact with dealer or manufacture for more details.

4. **«SELECT» key:** Cycle select. According to different things to be sterilized choose different temperature you need.
5. **«PACKAGE» key:** Package select. You can choose package or non package according to the instrument. Non package adopt 3 times pre vacuum, package adopt B class pre vacuum.
6. **«TEST» key:** It used for select the testing cycle. Press one time, B&D light goes on, means will run B&D cycle, press again, LEAK light goes one, means will run LEAK cycle. Press this key for 3s, two lights goes on at mean time, means to run cleaning cycle, for cleaning the steam generator.
7. **«DRY» key:** This key selects the Fast sterilization cycle. With unwrapped instruments, you may use the fast cycle by pressing the DRY key once, the indication light goes on and it indicates that machine is in fast cycle mode. Using the Small Steam Sterilizer made by us, you may still use this selection to shorten the cycle time.
8. **«STEP» key:** This is a post cycle selection key. When a sterilize cycle completes, you may press this key to keep on working or stop working. The indication light goes on, and it means that when the cycle is complete, the heating element will stop before the chamber door is opened. When the indication light is off, it means that the chamber will be kept warm after completing a cycle and heating up time is shorter during next cycle. If the user does not start any operation after completing the cycle, the system will enter into last cycle automatically an hour later.
9. **«START/STOP» key:** When you choose the right program, push this button then it starts a cycle. Keeping pressing the key 3 seconds before desiccation phase, the cycle of sterilization will be quitted. Then it will enter into 3 minutes dryness. In this program if you press the key again, the program of dryness will be over at any moment. If you press the key at end of dryness, the cycle will be quitted directly.
10. **«WATER» indication light:** When the water volume is low, the light is glinting.
11. **«READY» indication light:** When the light goes on, the machines gets ready for sterilization, and the operator can press "start" key to start a cycle.
12. **«RUNNING» indication light:** Working light. It means sterilize program is working.

 Caution: In case of faults, please contact your distributor or the manufacture immediately.

## 6. Operation

Please connect the power cord before you start the operation. Press the main power switch at right lower corner of the front panel. The power switch light turns on indicating that electric power is connected. The screen no display, automatically examine. If it is ok, the screen will show LD.

### 6.1. Fill water

If the water level is lower than the lowest water level. When the power switch is on, select the sterilize program, the "water" indication light may glinting while alarm Beeps. At this time, you need to refill the reservoir. Please open the water tank and pour water into it by container. Please pay attention to the water level during filling water, you should stop filling when the surface is close to the red water level mark.


Figure 6-1

#### Attention:

- 1 Distilled water must be good for prolonging life of the machine.
2. Don't tilt the Small Steam Sterilizer where there is water in the reservoir.
3. The used water tank should be emptied frequently. You should evacuate the used water reservoir usually when water reservoir is empty.

### 6.2. Start a cycle

When water indication light is off, the machine can be started.

- a. Press "SELECT" to choose one of the programmes. The key toggles between the cycles at 121°C, 134°C and B+ mode. One of the indication lights turns on to show the programme you selected (Fig 6-2).
- b. You may need to select instruments package mode. According to the instrument you put into the chamber whether it is bareness or sop up, press "package" to select package mode.
- c. You may select "DRY" key whether starts the fast sterilization cycle. Press "DRY" key, it will reduce the sterilization time. If you press the button "DRY" on the panel, there are three selections, 1min (with 1min LED on), 5min (with 5min LED on), and 10min (with two LEDs off).

Figure 6-2

- d. After selecting above setup, put the instruments into the Steam Sterilizer which need sterilize.


Figure 6-3

**⚠ Caution:** Instruments that are ready to sterilize should be placed in the tray to touching one another allowed space for steam to pass through. It is suggest that the tray holder be used to put the tray with instruments inside the chamber to avoid burning injuries.

- e. After the instruments are loaded, you may close and lock the chamber door by tuning the door knob clockwise to it is locking position. The "ready" indication light goes on. (Fig 6-4, 6-5, 6-6).

Figure 6-4


Figure 6-5

Figure 6-6

**⚠ Caution:** The door knob must be tuned to its lock position, otherwise, there might be steam leak during a cycle and the alarm goes on and displays error code "E6". If the chamber is in keep warm mode, it is sometimes difficult to close the chamber door because there is still steam inside. You can keep the door open for a while until steam and moisture is vaporized, and the door is easy to close. You may also push the door inwards with some force while tuning the door knob to its locking position. If there is still steam leak after you did the above, you may try to readjust the door mechanism by a spanner as per item 7. Door adjustment.

- f. Press "START/STOP" key, the "ready" light goes off and display "HE" the machine starts a sterilization cycle. It takes about 30-60 minutes to complete an automatic cycle that includes heating up, sterilizing and drying. The cycle time depends on the volume of the instruments loaded, the chamber temperature at the beginning of the cycle and the cycle program you selected.

- g. When a cycle completes, the display "Ed" the indication light "RUNNING" Turns off and the indication light "READY" turns on with sounds. At this time you may then open the chamber door and unload the instruments which are sterilized. (Fig 6-6) At mean time, temperature indicate "prt", it means that it is performing print, whether you have install printer, the process can not be omitted. But you needn't care about it. It will not effect the process of opening the door and unloading the instruments. The indication will disappear when printing is over.

Figure 6-7

**⚠ Attention:** To ensure small steam sterilizer's on the safe side, When water level low allann sounds, Please be advised add some water first. When you add the water, please don't forget to let condenser water out.

**⚠ Caution:** Always use the tray holder to load or unload the tray with instruments in it, in order to avoid burning injuries.


Figure 6-8

- h. If you decide to stop using the small steam sterilizer for a while, please turn the main power switch off, or unplug the power cord should you not want to use the machine for a long time.

**⚠ Caution:** Never try to open the chamber door if and when the pressure display doesn't show "0.0". It is suggested that chemical indicator or bio-indicator be used for each and every cycle to ensure the reliability of each sterilization cycle run.

**6.3. The instrument rack should be put inside the boiler as the following diagrammatic sketch shows.**

The first way (fig. 6-9, fig. 6-10)

Figure 6-9


Figure 6-10

The second way: clockwise rotate the instrument rack by 90°, based on the first way (fig. 6-11, fig. 6-12).

Figure 6-11


Figure 6-12

**⚠ Notes:** the circle stamped angle is upalong.

## 7. Door adjustment

Under normal circumstance the chamber door lock needs no readjustment. Once steam leaking occurs (the seal fails), you may use the spanner to turn the door seal thightener.

- Open the door first.
- Insert the spanner in the gap beneath the plastic cover; use the spanner to lock on the adjusting nut (Fig. 7-2). Turn the nut counter clockwise as the figure below (Fig. 7-1). This will tighten the sealing plate.
- Turn the nut until the sealing plate is tight. If the door knob is too tight, you may also turn the nut clockwise to loosen it as long as the seal does not leak.

Caution: Never try to readjust the chamber door while the door is locked.



Figure 7-1


Figure 7-2

## 8. Printer installation and use

Due to the requirement of customer, D\_Devices sterilizer can adopt prescriptive printer. Connect the printer to the sterilizer with cables see fig 8-1. Switch on the power, two lights on the printer panel will goes on, means printer has been installed correctly. According to system setup you can setup whether you need print the process or not.

Figure 8-1

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<b>S/N:</b>	Serial No (default is close the function)
<b>Date:</b>	Date
<b>Cycle:</b>	Cycle Times
<b>C.N:</b>	Cycle name
<b>C.S.T:</b>	Cycle start time
<b>Vp1, Vp2, Vp31:</b>	Vacuum
<b>Pp1, Pp2, Pp31:</b>	Pressure strenghen
<b>H.S.T:</b>	Sterilization start time
<b>Max, Min, Ave:</b>	Max. Min and average pressure or temperature of chamber
<b>H.E.T:</b>	Sterilization end time
<b>D.S.T:</b>	Drying cycle start time
<b>C.E.T:</b>	Cycle end time
<b>Result:</b>	Result
<b>Operator:</b>	Operator, you can sign your name if it needed.

Sterilization Records	
Date: Y11M05D11	
Cycle: 00001	
C.N.= 134 2.10	
-B-SP=6mDP=10m	
-----	
C.C.S.T.: 13:05:30	
Vp1: 13:11:27	
0.00 048.9	
Pp 1: 13:12:57	
-0.81 050.2	
Vp2: 13:16:33	
0.29 107.5	
Pp2: 13:18:01	
-0.81 072.1	
Vp3: 13:20:36	
0.30 112.4	
Pp3: 13:22:14	
-0.81 081.4	
H.S.T.: 13:27:55	
Max: 2.14 134.6	
Min: 2.02 134.5	
Ave: 2.08 134.5	
H.E.T.: 1 3:34:06	
D.D.S.T.:13:34:43	
0.30 115.5	
C.E.T.: 13:45:32	
-----	
Result: Passed	
Operator:	
-----	
Thanks	

Figura 8-2

**⚠ Caution:** D\_Devices sterilizer only can select printer offered by D\_Devices at present. Please print when one cycle stops and another doesn't start, and the power switch is still on. Press "START" key for 3 seconds to print.

## 8.2. USB module function

- a. Main function:
  - Storage of printed data of sterilization Built-in FLASH memorizer
  - U-disk. Download
- b. Basic parameters:
  - Storage Capacity: 2M bytes
  - Quantity of printed data could be stored: 2048
  - Support FAT files
- c. Operation instruction::
  1. Turn on the power, the mode of USB will start automatically. After self-test, it goes into alert with the indicator flickering-1 time.
  2. After the sterilization cycle ends, the data of sterilizer is sent to the USB mode (The printer program should be turned on). Indicator will flicker-1 time after the data-transfer finish.
  3. During the alert, insert U-disk. The mode indicator will be on, with the U-disk indicator flickering. And now, it is in the process of data-download. Don't turn the power off or draw out the USB-disk, which will cause the data lose.
  4. After the download finish, USB mode indicator will flicker twice, reminding the user to draw out the U-disk.
  5. During the sterilization cycle, avoid data download, or else, the data will lose.
  6. Connect U-disk to computer. You will see files like \_2000\_, where the downloaded data is kept, following turns of Year, Month and Day. Easy to print it out.
  7. USB mode error (According to the flickering times of USB mode indicator) and solution:

Times	Definition	Solution
1	Self-test passed	Normal
2	Operation finish	Normal
3	No definition now	No definition now
4	U-disk needs formation	Format U-disk in FAT files
5	Not enough capacity	Delete the files in U-disk
6	U-disk needs formation	Format U-disk in FAT files
7	USB mode error	Contact us

Caution: This function is optional to the demand of customer.

## 9. System setup

- a. Enter function setup mode: Press "STEP" key for 20s when on non-work state (time display "Ld") pressure display "1", enter function setup mode.
  - b. Function mode setup: With 7 options, press "SELECT" key, you can select 1-7-1 cycle display. The definition and display of 7 options as following:
    - Option 1: year, display "yer", press "STEP" or « ↑ » key increase, press "DRY" or « ↓ » key decrease, setup numerical value.
    - Option 2: month, display "nth", press "STEP" or « ↑ » key increase, press "DRY" or « ↓ » key decrease, setup numerical value.
    - Option 3: day, display "day", press "STEP" or « ↑ » key increase, press "DRY" or « ↓ » key decrease, setup numerical value.
    - Option 4: hour, display "hor", press "STEP" or « ↑ » key increase, press "DRY" or « ↓ » key decrease, setup numerical value.
    - Option 5: minute, display "nin", press "STEP" or « ↑ » key increase, press "DRY" or « ↓ » key decrease, setup numerical value.
    - Option 6: Printing, it display "Pvt", press "STEP" key then you can choose the cycle between "on" and "off". "on" stand for "print", "off" stand for "no printing".
    - Option 7: Exit, when "Eit" displays, please press "START/STOP" to save and exit setup mode.
    - Option 8: Sterilization time adjust, display "S-T" when you setup this option, the sterilization time (SP) will increase at former value. The adjustable range are following 3 kinds. One is increased by 1 minute, second is increased by 2 minutes. Third is increased by 3 minutes.
    - Option 9: temperature data adjust automatically option, display D-P. When you setup this option, you can adjust the temperature more or less than 2C-3C in order to meet your request. But we suggest user should setup under guidance of technician or not setup.
- Press "SELECT" key after setup, you will come back to option 1. At this time, you must press STEP key 3 seconds in option 7. Then you can start the setup of option 8 and option. Otherwise it only can start cycle between option and option 7.

 Caution: If some function setup was changed by pressing the key without intention, please switch off the power, then restart it. It will exit from the wrong setup. In this case, the function won't be affected by this.

## 10. Maintenance & Maintain

### 10.1. How to clean the water tank

Please clean the sterilizer water tank and waste water tank regularly. The process is following:

- a. Empty the tanks water, please refer to item 10.2 draining the water.
- b. Open the tank cover board, you can use screw driver to loosen the five screws just as fig 10-1 indication.
- c. Pull up the water tank cover board by one hand, just, as fig 10-2 indication.
- d. Pull up the tank cover completely and you can disassemble the tank cover. As fig 10-3 indication.
- e. Now you can start to clean the tank, use cotton to dip in some alcohol or medical disinfectant to scour the tank wall, wash it by distilled water, then empty the tank and dry the tank. After it you should take out the filter net to clean.
- f. After cleaning the tank, assemble the tank cover and tighten the screws.

Figure 10-1


Figure 10-2

Figure 10-3

### 10.2. Draining of water from the machine

Drain the condenser water collector and water reservoir. Connect the water outlet using the one end of the tube without a connector and the other end of the tube to a drainage sink. Then turn the draining knob anti-clockwise to drain the water from the machine (Fig 10-4).


Figure 10-4

**⚠ Caution:** Before any kind of maintenance turn off the power to the small steam sterilizer and allow it to cool down. Use only distilled water for sterilizing and cleaning. The instruments being sterilized must be cleaned. Instruments with oil or impurities must be packed in a sterilization bag.

### 10.3. Clean and install the drain filter

Over time small impurities will deposit on the drain filter, and this may affect vacuum and drying. In order to prolong the life of the drain filter please use the following cleaning procedure:

1. Chamber water filter.

2. Back-out the water filter.

3. Take out the filter net, clean it by soft brush or replace a new one.

**⚠ Note:** Please check the O ring of Draining Filter during cleaning. If this part is broken please replace it with O ring in the accessories supplied.

#### 10.4. Change fuse (Fig 10-5)

- a. Unplug the small steam sterilizer before changing the fuse
- b. Unscrew the fuse holder with a screw driver counter clockwise.
- c. After screwing the fuse holder, pull the fuse out.
- d. Make sure to replace the faulty fuse with a correct fuse.
- e. Put the correct fuse back to the fuse holder, plug into the holder lightly.
- f. Put back the fuse holder and screw it on clockwise with a screw driver, make sure it will not pop-up.

Figure 10-5

#### 10.5. Clean the seal plate periodically

After using some time, the seal plate might have some scale build-up that might cause leaking. You should clean the seal plate periodically, using a soft cloth soaked with distilled water to wipe the rubber seal and the seal plate clean (Fig 10-6 & 10-7). If there is still leaking, you may need to take the rubber seal out for cleaning. Should the seal be broken, you need to change with a new one.

Figure 10-6

Figure 10-7

Remark: providing with seal packing collar for life.

#### 10.6. Change of the seal

You will need a flat screw driver without a sharp tip. Pay attention about head of the flat screw driver. The install step according to Fig 10-8.

Figure 10-8

- a. Hold the edge of the seal with one hand and insert the screw driver beneath the seal with your other hand and then pull out the seal slowly.
- b. Once a section of the seal is out, you may pull the whole sea! out slowly with both of your hands. You can now clean the seal groove and the seal. If the seal is broken or damaged, you must replace it with a new one.
- c. After the groove and the seal are cleaned, you may remount the seal back in place.

**⚠ Caution:** The seal should be mounted evenly inside the groove with first 4 positioning points in place. You may use your hand to press the seal in place evenly.

- d. Attention: When inserting the seal, the inner edge might evaginate. You may use a screw driver to press the edge down in the groove.

**⚠ Caution:** Power the machine off and let it cool down before you change the seal and avoid any burning injuries.

**10.7. Cleaning Steam generator**

- a. Aim: Clean the steam generator and steam pipe, to avoid the furring block the pipe.
- b. Description: 20 min cleaning + 10 min drying.
- c. Purpose: to clean the steam generator furring, take it out of the machine through pre vacuum.
- d. Operation: when the machine is inactivity, press "SELECT" for 5 seconds, the B&D light and LEAK light will go on at mean time, it means you chose the cleaning cycle. Close the door and press "START", it will start to clean.
- e. Cleaning remind: after cleaning, every 100 cycle, B&D light and LEAK light will going on at the mean time and glitter. Remind you to clean the steam generator.
- f. Please keep regular cleaning of steam generator to prolong the service life of the steam generator.
- g. Before executing clean steam generator program, we suggest to mat a wet towel on the bottom of the chamber (as shown in below picture), prevent excessive impurities absorption on the draining filter, which may influence pre-vacuum and dry effect.

Friendly remind: Cleaning function will not affect the small steam sterilizer cycle, only remind message.



Figure 10-9

**10.8. If the machine need to be maintained or the electric spare parts need to be changed, please cut off the power first. And ask your dealer or manufacturer to change it.**

**10.9. Maintenance operation frequency**

Maintenance operation	Frequency
a. Clean the water filter	every week
b. Clean the water tank and the water tank filter	every week
c. Clean the tray screen	every week
d. Clean the sealing ring	every month

For the replacement of spare parts, we recommend:	
a. Replace the air filter	every year or 1000 cycles
b. Replace the main silicone hose	every 2 years or 2000 cycles
c. Replace the sealing ring/door rubber	every year or 200 cycles
d. Replace the solenoid valve	every 4 years or 4000 cycles

## 11. Transport & Storage

**11.1. Preparation for transportation and storage:** Power the machine off, unplug the power cord and let the chamber cool down.

**11.2. The small steam sterilizer should be transported and stored in the follow conditions:**

Temperature: -40 ~ +55 °C

Humidity: ≤ 85 %

Atmospheric pressure: 500 hPa ~ 1060 hPa


 Caution: Don't drag the machine while moving.

## 12. Warranty

The manufacturer provides a warranty for repair at not change to the customer for two years after the invoice date with presentation of the invoice. Following are the disclaimers:

- a. Damage due to improper moving, installing and operating the machine;
- b. Damage due to dismantling the machine by unauthorized personal;
- c. Don't have guarantee of end user.
- d. Inappropriate operation on the machine without following the operator's manual
- e. Damage caused by power surge, fire or other external matters.
- f. Not having performed the annual maintenance or every 1000 cycles. Whichever comes first.

## 13. Accessories

- |  |   |
|--|---|
| 1. Instrument tray rack 1 pc and Instrument tray 3 pcs | 6. Filtrable net 1 pc (water tank)      |
| 2. Power cord 1 pc                                     | 7. Fuse (F12AL 250V) 2 pcs              |
| 3. Mainboard Fuse (F3AL 250V) 1 pc                     | 8. O ring (water draining filter) 2 pcs |
| 4. Water outlet tube 1 pc                              | 9. Door adjustment spanner 1 pc         |
| 5. Tray holder 1 pc                                    |   |

### Note:

1. The picture of spare parts in the attachment may be a little different to the parts supplied. D\_Devices keep the right of change the accessories. Please check the spare parts supplied with small steam sterilizer.
2. During disposing the obsolete device, please contact with professional recycle company. Don't throw them away.

## 14. Appendixes

### 14.1. Work mode chart

	134°C package	134°C non package	121°C package	121°C non package	134B+	BD test	Vacuum test cycle	CLEAN test
<b>Sterilizing temperature</b>	134	134	121	121	134	134	—	—
<b>Sterilizing pressure</b>	2.10	2.10	1.10	1.10	2.10	2.10	—	—
<b>Classes</b>	B	S	B	S	B	B	—	—
<b>Pre-heating (min)</b>	0-7	0-7	0-7	0-7	0-7	0-7	—	—
<b>Vacuumize time (min)</b>	10	4	10	4	10	10	4	Clean 20
<b>Admission time (min)</b>	5	5	5	5	5	5	—	Drying
<b>Sterilizing time (min)</b>	6	4	20	16	18	3.5	Keep pressure 10	—
<b>Drying time (min)</b>	10	10	10	10	10	1	—	10
<b>Quick drying time</b>	1/5	1/5	1/5	1/5	1/5	—	—	—
<b>Total time (min)</b>	38(29/33)	30(21/25)	52(43/47)	42(33/37)	50(41/45)	26.5	14	30

### 14.2. Error code chart

Item	Enor Code	Error	Description
1	E1	Seren fault	Check sensor
2	E2	Over pressure (>2.6Mba)	Check inner temperature sensor pressure sensor
3	E3	Chamer over temperatura (>145 degree)	Check outer temperature sensor
4	E4	Keep warm keep pressure faild	Rerun parameter autoadjust program
5	E5	Exhaust failed (20S<0.3Mba)	Check air relery valve vacuum valve
6	E6	Door opened when working	Check door switch electrical lock
7	E7	Steam generator over temperature (>230 degree)	Check steam generator temperature sensor
8	E8	Pre-heating fault>20min	Cheack heating ring or heating pipes
9	E9	Fill stern fault>30min	Check fin steam pipes
10	A01	Unsuitable water quality.	The use of water with a conductivity of more than 15µS/cm can affect the sterilization process and damage the sterilizer.
11	A02	Lack of water in the clean water tank	Refill tank with clean water
12	A03	Dirty water tank full	Empty the dirty water tank

### 14.3. D\_Devices working chart

134°C sterilization cycle working chart

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
## 14.4. Circuit diagram

VACUUM SMALL STEAM STERILIZER USE (LED)

14.5. D\_Devices piping

## 15. Troubleshooting

- a. During venting or drying, the thump sounds. It means that the condensed water collector is full. You must drain the water from The collector.
- b. Hissing from the chamber door (Solving this problem according to item 7 Door Adjustment).
- c. During a cycle, a banging sounds and large steam leaks from the chamber door and displays "E6" code. Please verify whether the door knob is properly turned and locked. If the steam leaks at the beginning of a cycle, you need to change the seal according to instructions at page 18.
- d. When a cycle completes, there are water drops on the instruments. You may correct the problem by:
  1. Please check the drain filter.
  2. Reduce the load of the instrument in the chamber, especially the porous and dressing and other water absorbent materials.
  3. You may also wait for a while at the end of a cycle "Ed" for better drying result. The small steam sterilizer will continue heating the chamber until you open the door.

## 16. Important information

- a. The small steam sterilizer must be installed on a leveled counter top.
- b. When operated at a higher elevation of more than 500 meters, The machine must be reprogrammed. Please get in touch with your distributor for instructions.
- c. Distilled water must be used to prolong the life of the machine.
- d. The radiation openings in the machine must not be blocked or covered.
- e. Instruments should be placed on the tray allowing space between one another and thus allowing steam to pass through.
- f. Condensed water collector must be emptied from time to time. It should be emptied when refilling the water reservoir (Draining of the collector can be followed as per item 10.2).
- g. Chamber door knob must be turned to its locking position.
- h. Never try to open the chamber door before pressure displays "0.0".
- i. To avoid burning injuries, power off the machine and let it cool down before replacing the seal.
- j. No dragging when the machine is moving.
- k. Power connection must be grounded properly.
- l. Must be provided enough power (1800VA).
- m. While the temperature of outside is lower, please pre-heating the machine without instruments for 30 minutes before sterilize.
- n. It is normal if the machine alarms when it starts for the first time. If it alarms again, please contact your distributor.

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**PROCLINIC S.A.U.**
**Ningbo Runyes Medical Instrument Co., Ltd.**