FICHA TÉCNICA

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Almacenar a temperatura ambiente de 24°C (75°F).

TECHNICAL DATA SHEET

PORCELAIN ETCH PROCLINIC EXPERT		
REFERENCE	78577: 9% Hydrofluoric Acid.	
MANUFACTURER	B.J.M. Laboratories Ltd. 12 Hassadna St. Industrial Park, Or-Yehuda 6022011 - ISRAEL Tel. 972-74-7000111 Fax. 972-3-7353020	
CLASSIFICATION OF THE PRODUCT	Medical Device class IIa, according to Directive 93/42/EEC.	
DESCRIPTION	Viscous hydrofluoric acid for etching any type of ceramic which creates micropores in the treated surface, thus achieving a strong bond to composite.	
COMPOSITION	Hydrofluoric acid, Water, Xanthan gum.	
PACKAGING	Package containing 2 x 1,2 ml syringes and 4 delivery tips.	
INTENDED USE	Cementation of ceramics for: <ul style="list-style-type: none"> • Preparation of porcelain crowns, veneers and inlays before definitive cementation. • Repair of ceramic pieces or fractured metal-ceramic. 	
SHELF LIFE	2 years.	
INSTRUCTIONS FOR USE	<ol style="list-style-type: none"> 1. Isolate the working field intra-orally with rubber dam. 2. For fractured ceramic repairs, roughen the surface of the ceramic to be etched using a medium grain bur. Remove the shiny surface and create enough space for the composite. Rinse and dry the ceramic surface. 3. Apply Proclinic Expert Porcelain Etch to the ceramic with a brush and leave to act for 60 seconds. 4. Rinse the ceramic with copious amounts of water to remove the acid and dry. The ceramic should have a matt, frosted appearance. If this is not the case, repeat the etching process until this appearance is achieved. 5. Apply Proclinic Expert Silane onto the etched ceramic and dry well with a stream of oil-free and water-free air. 	
SAFETY INFORMATION	Use with caution. Avoid contact with mucosa, eyes and skin. In case of accidental contact, wash immediately with copious amounts of water and consult a specialist. It is recommended to isolate the working field intra-orally with a rubber dam.	

HANDLING AND STORAGE	Store at room temperature of 24°C (75°F).
---------------------------------	-------------------------------------------