FICHA TÉCNICA

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TECHNICAL DATA SHEET

SILANE PROCLINIC EXPERT		
REFERENCE	78578	
MANUFACTURER	B.J.M. Laboratories Ltd. 12 Hassadna St - Industrial Park Or-Yehuda 6022011-ISRAEL Tel. 972-74-7000111	
CLASSIFICATION	Medical Device class IIa, according to Directive 93/42/EEC concerning medical devices.	
DESCRIPTION	Porcelain Etch is a viscous, hydrofluoric acid for etching any type of ceramic, which creates micropores in the treated surface thereby achieving a strong bond to composite. When applied to the etched ceramic, Proclinic Expert Silane improves the restoration's bonding strength and longevity.	
COMPOSITION	Glycidoxypropyltrimethoxysilane 3%, Ethyl alcohol 97%.	
PACKAGING	Package containing 2 units of Silane 2 ml + 2 empty syringes of 1.2 ml + 4 applicator tips.	
INTENDED USE	<u>Cementation of ceramics for:</u> <ul style="list-style-type: none"> • Preparation of porcelain crowns, veneers and inlays before definitive cementation. • Repair of ceramic pieces or fractured metal-ceramic. 	
SHELF LIFE	2 years.	
INSTRUCTIONS FOR USE	<ol style="list-style-type: none"> 1. Isolate the working field intra-orally with rubber dam. 2. For fractured ceramic repairs, roughen the surface of the ceramic to be etched using a medium grain bur. Remove the shiny surface and create enough space for the composite. Rinse and dry the ceramic surface. 3. Apply Proclinic Expert Porcelain Etch to the ceramic with a brush and leave to act for 60 seconds. 4. Rinse the ceramic with copious water to remove the acid, and dry. The ceramic should have a matt, frosted appearance. If this is not the case, repeat the etching process until this appearance is achieved. 5. Apply Proclinic Expert Silane on the etched ceramic and dry well with a stream of oil-free and water-free air. 	
SAFETY INFORMATION	Use with caution. Avoid contact with mucosa, eyes and skin. In case of accidental contact, wash immediately with copious amounts of water and consult a specialist. It is recommended to isolate the working field intra-orally with a rubber dam.	

HANDLING AND STORAGE	Store at room temperature of 24°C (75°F).
---------------------------------	-------------------------------------------