Circuiti e maschere per
Circuit and mask for
Schlauch mit maske für
Circuitos y mascarillas para
Circuitos e máscaras para

Master Flux

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1. INTENDED USE
2. CONNECTION
3. CLEANING

1 INTENDED USE

The Tecno-gaz circuits and masks are accessories of device Master Flux and should only be used with the device Master Flux for the supply of medical gases to the patient.


TECNO-GAZ is not responsible if the circuits and masks are used in a manner different from that stated above.

TECNO-GAZ is not liable for misuse of circuits and masks.

The masks are sold separately from the MASTER FLUX ONLY AS PARTS OF THE DEVICE.

2 | **CIRCUIT AND MASK**


Components CIRCUIT

MF820ZMF	CIRCUIT COMPLETE	
MF835ZMF	MASK size 1	
MF836ZMF	MASK size 2	
MF837ZMF	MASK size 3	
MF838ZMF	MASK size 4	
SMFA197	CIRCUIT EXHAUST PASSIVE	
MF824ZMF	CIRCUIT CONNECTION OF ASPIRATOR MASTER FLUX PLUS	
MF800IMF	CIRCUIT CONNECTION OF ASPIRATOR MASTER FLUX SMART	



The circuit **MF824ZMF / MF800IMF** is composed of:

<p>TERMINAL ASPIRATION</p>	
<p>REDUCTION FOR SUCTION SALIVA</p>	
<p>ONE-WAY VALVE</p>	
<p>MASCHIO PORTAGOMMA Only in MF800IMF for the MASTER FLUX SMART version</p>	

the circuit **MF820ZMF** is composed of:

<p>Fitting connection to flowmeter box (*) ▶ Connection 22 M/F 90°</p>	
<p>Two tubes length 1.8 m (material PVC/PE) ▶ Transparent and flexible tube. ▶ Corrugated tube (rigid spiral) externally (anti-crushing) smooth inside for low compliance and noise. ▶ Flexible connectors for easy insertion.</p>	
<p>Exhaust fitting connection (material: stirolox) (*) ▶ Reduction 22-15</p>	
<p>Reducing exhaust connection (material: TPE, PP) (*) ▶ Elastomeric connector 9-11 mm</p>	
<p>(*) DO NOT USE IN THE MASTER FLUX SMART VERSION</p>	
<p>PIPE CLIP SUPPORT AND VELCRO BAND ▶ Allows correct positioning of the tubes, the adhesive clip can be secured behind the headboard or via the Velcro strip can be positioned on the headboard without a secure fixing. ▶ It allows a positioning with fixed or rotary movement, according to the needs.</p>	

The mask is composed of:


<p>Kit Silhouette nasal and Support front Code 2ZMFA0050 size 1 Code 2ZMFA0051 size 2 Code 2ZMFA0052 size 3 Code 2ZMFA0054 size 4 (material: silicon) ▶ It reduces N2O environment emissions ▶ Create a mask placement comfortable</p>	
<p>Shell (material: PC) ▶ Transparent</p>	
<p>Small support front (material: PC) ▶ Transparent For mask size 1 and size 2</p>	
<p>Large support front (material: PC) ▶ Transparent For mask size 3 and size 4</p>	
<p>Fitting Y (PC) with valves (PC/silicone) ▶ Allows the correct flow of the gas mixture by preventing the exhaled gas can flow back into the supply circuit. ▶ The small room eliminates the risk of permanence CO2</p>	
<p>HEAD RIBBON WITH FASTENER FOR HEAD RIBBON ▶ The mask stabilizing and regulating via a velcro the length to adapt it to the patient.</p>	

CONNECTION INSTRUCTIONS

Check the integrity of the components.

Choose the most suitable mask for the patient.

Place the mask and adjust the head ribbon. The mask must be stable and comfortable for the patient. Connect the pipes, paying attention to the arrow on the Y fitting. The arrow indicates the flow direction. The non-return valves, placed in the fitting, push the gas flow in the direction indicated by the arrow.

Gas inlet pipe: Connect a hose to the special corrugated tube (Item G) to the flow control box by inserting the connector connection flow control box ".

Connect the free end to the Y connector placed on the mask, WITHOUT SIDE ARROW.


Gas inlet pipe: Connect the other corrugated pipe to the appropriate hose (Pos P) through the elastomeric connector. Connect the free end to the Y connector placed on the mask, SIDE WITH ARROW.

Fix the clip to stabilize the pipes behind the headboard of the chair.

And 'possible fixing with adhesive (fixed) or through the appropriate supplied Velcro (detachable).

Delivery knob
Flush button
Dosing knob


G Gas inlet hose fitting

F Anaesthesia bag

P Gas outlet hose fitting

Delivery knob

A

C

Dosing knob

B

Flush button

P

Gas outlet hose fitting

G

Gas inlet hose fitting

F


Anaesthesia bag

\

EVACUATION EXHALED GASES

The operator can choose the best solution for processing the exhaled gases, using the following connections:

PASSIVE EXHAUST	
UNIT-MOUNTED MASTER FLUX	WALL-MOUNTED MASTER FLUX
<p>Connect the hose, exhaust gas circuit mask to the hose in the center at the center of the unit.</p> <p>Connect the circuit SMA197, to the hose connection at the bottom of the cabinet.</p> <p>Connect the end of the spiral tube to the outside and make sure there is a final filter.</p> <p>The discharge must take place in an open, well-ventilated.</p>	<p>Connect the exhaust gas pipe of the appropriate circuit hose (SMFA257) after the spiral tube (SMFA197).</p> <p>Connect the end of the spiral tube to the outside and make sure there is a final filter.</p> <p>The discharge must take place in an open, well-ventilated.</p>
	

ACTIVE EXHAUST	
UNIT-MOUNTED MASTER FLUX	WALL-MOUNTED MASTER FLUX
<p>Connect the hose, exhaust gas circuit mask to the hose in the center of the cabinet.</p> <p>Connect the circuit MF824ZMF/MF800IMF to the hose on mobile.</p> <p>Connect the terminal to the suction.</p>	<p>MF824ZMF/MF800IMF connect the circuit to the discharge circuit.</p> <p>Connect the terminal to the suction.</p>
	
<p>SUCTION MF824ZMF/MF800IMF</p> 	
<p>SETTING THE VALVE</p> <p>Turn the ring perforated valve (SMFA236) until closing holes.</p> <p>Start dispensing OXYGEN setting ~ 10 NI / min.</p> <p>Start the surgical aspirator, set the minimum intake.</p> <p>Observe THE BAG ANESTHESIA.</p> <p>Turn the ring perforated valve gradually opening the holes until you get the ideal state of the BAG ANESTHESIA "HALF FULL HALF EMPTY".</p> <p>CAUTION: By varying FLOW, it may be necessary, adjust the correct intake flows of exhaled gases.</p>	

4 CLEANING

Use detergents that are not aggressive to plastic

Clean the surfaces of the device with a dry cloth.

Before and after use, clean, clean or sterilize the parts in contact with the patient or with his breath.

Cleaning table for MASTER FLUX parts

Attention:

Use a disinfectant that is compatible with plastics (avoid cream products).

Do not spray the disinfectant directly onto the component to be cleaned, but onto a cloth.

Use disinfectants in the dilution indicated by the manufacturer.

ITEMS:	DESCRIPTION	IMAGES				
MF835ZMF MF836ZMF MF837ZMF MF838ZMF	SHELL COMPLET FITTING VALVE ENF FRONT SUPPORT					
		MATERIAL	THERMAL DISINFECTION	N ° CYCLES	STEAM STERILIZZAZIONE	N ° CYCLES
		POLICARBONATO - SILICONE	SI	10*	134°	10*

MF835ZMF MF836ZMF MF837ZMF MF838ZMF MF824ZMF MF800IMF	SILHOUETTE NASAL AND SUPPORT FRONT - RINGS - EXAUST CIRCUIT PASSIVE					
		MATERIAL	THERMAL DISINFECTION	N ° CYCLES	STEAM STERILIZZAZIONE	N ° CYCLES
		SILICONE - ALUMINUM - PLASTIC	SI	10*	134°	10*

MF835ZMF MF836ZMF MF837ZMF MF838ZMF MF820ZMF 3MEDN0002 SMFA197	HEAD RIBBON WITH FASTENER - CIRCUIT COMPLETE WITH CONNECTORS - RESERVOIR BAG 3LT - EXAUST CIRCUIT PASSIVE					
		MATERIAL	THERMAL DISINFECTION	N ° CYCLES	STEAM STERILIZZAZIONE	N ° CYCLES
		NEOPRENE - PE - PVC SBS - TPE - VELCRO	NO	-	NO	-

OPTIONAL		MATERIAL	THERMAL DISINFECTION	N ° CYCLES	STEAM STERILIZZAZIONE	N ° CYCLES
3MEDN0003	RESERVOIR BAG 3LT -	SILICONE	134°	10*	134°	10*
3MEDT0018	TUBES LENGTH 1.5 MT	HYTREL - POLIESTERE	134°	10*	134°	10*
3MEDT0019	TUBES LENGTH SINGLE USO 50MT	LDPE	NO	-	NO	-

(*) Maximum recommended number of cycles.

Manual cleaning and manual disinfection should only be considered where no facilities exist for automated washing and disinfection, or automated cleaning is contraindicated.

CLEANING – MANUAL

Hand wash using a mild detergent such as a neutral pH liquid soap or tri-sodium phosphate. Rinse thoroughly with clean water, and wipe with an appropriate soft absorbent material that will not leave fibers on the device. Hard brushes or other materials that might damage the products should not be used. Do not use a washing machine. Dry at room temperature. Reservoir Bags and Rubber Masks may be dried by hanging in a warm cabinet at a maximum temperature of 50° centigrade. Do not dry on a hot surface.

DISINFECTION – MANUAL

Use a proprietary disinfecting agent. Concentration and immersion times stated by the manufacturer must be strictly adhered to. Compatibility with polychloroprene, stainless steel, polysulfone, polycarbonate and polypropylene materials used in the manufacture of the product must be established prior to conducting any disinfection cycle or routine. Rinse thoroughly with water of a suitable quality, and wipe with an appropriate soft absorbent material that will not leave fibers on the device. Dry at room temperature.


STERILISATION

Products can be sterilized by moist heat (autoclaving) up to temperatures of 137° centigrade in accordance with methods specified in ISO17665. Do not sterilize with dry oven sterilizers.


MASK ASSEMBLY

After cleaning the device, assemble the mask according to the diagram below.

**ATTENTION TO THE DIRECTION OF FLOW
ARE THESE VALVES
BINDING THE FLOW IN THE DIRECTION SHOWN**


Circuit INLET

INLET


Circuit OUTLET

OUTLET

Exhaled Gas Circuit identification

For easy identification of discharge circuit, during are provided template-circuit, assembly, two blue rings (Cod. 1ZMFA0082).

Figures A B C D show how to use the rings. The indicator is autoclavable.

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TECNO-GAZ S.p.A