Coltolux Comfort
LED Curing Light
Owner's Guide


Table of Contents

EN	Owner's Guide.....	4
DE	Benutzerhandbuch.....	20
FR	Guide de l'utilisateur.....	36
ES	Guía de uso.....	52
IT	Guida per l'utente.....	68
NL	Gebruikershandleiding.....	84
SV	Bruksanvisning.....	100
DA	Brugsanvisning.....	116
FI	Käyttöohjeet.....	132
NO	Brukerhåndbok.....	148
PT	Manual do utilizador.....	164
EL	Εγχειρίδιο χρήσης.....	180

CAUTION: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

Prior to installation and operation of the device, please read all instructions provided in this user manual carefully.

Product Description

The Coltolux® Comfort is a light emitting diode (LED) curing light used for the polymerization of light-cured materials by dental professionals.

It is characterized by:

- Simple, one-button operation
- Lightweight and ergonomic design
- Flat tip (90°) for easier posterior access
- Homogenous beam profile
- Inductive charging
- User replaceable battery
- Durable and easy to disinfect material construction

Indications for Use

The Coltolux® Comfort LED curing light is indicated for the polymerization of intra- and extraoral dental restorations fabricated from materials (commonly referred to as visible light cured dental restoratives and/or composites and/or dental adhesives) that contain Camphorquinone as the photoinitiator.

Contraindications

Coltolux® Comfort should not be used by any persons who are taking photosensitizing drugs, or who have a known sensitivity to light.

Intended Users

The intended user population are dentists and registered dental assistants. The Coltolux® Comfort curing light is not intended to be used by the general population.

Clinical Benefits and Performance Characteristics


The Coltolux® Comfort does not directly treat any medical conditions. It is intended to work in conjunction with restorative dental materials to seal and protect decayed tooth enamel and dentin, to bond dental braces for treatment of malocclusions, and to repair/restore enamel and dentin damaged by physical trauma.

Compatible Materials

The Coltolux Comfort is designed to cure dental polymer-based restorative materials that contain Camphorquinone as the photoinitiator in the 440-470 nm range. See polymer-based restorative material manufacturer's complete directions for use for specific product compatibility and curing recommendations.

Light Output Characteristics		
Powered Polymerization Activator: Coltolux® Comfort Curing Light		
Class: 2 Type: 2		
Light output description: 11mm / 1.1cm diameter (circular cross-section)		
Light output cross-sectional area: 380 mm² / 3.80 cm²		
Irradiation time used: 20 seconds		
Mode used: Standard		
Peak wavelength: 449 nm		
Maximum light intensity: 1,440 mW/cm²		
Wavelength region nm	Radiant flux mW	Radiant exitance mW/cm ²
<380	<1	<1
380 to 515	1,140	1,200 (+/- 240)
>515	<1	<1

Coltolux® Comfort Wavelength Output

Handpiece

- Power Supply: Lithium-ion battery
Nominal Voltage: 3.7V
Capacity: 850 mAh

- Number of 20-second cycles with a new, fully-charged battery: 120 (minimum)

- Dimensions: 8.46 in / 215 mm length
0.78 in / 20 mm diameter

- Weight: 2.7 oz / 77 g

- Duty Cycle: 2 min on, 5 min off (at ambient temperature 77°F / 25°C)
1 min on, 10 min off (at ambient temperature >77°F / 25°C)

- Operating Temperature: 50°F - 100°F / 10°C - 38°C

- Relative Humidity: 80% at 77°F / 25°C

- Atmospheric Pressure: 500 hPa – 1200 hPa

- Applied Part: Type B

- Ingress Protection Rating: IPX3

Charger Base / Power Supply

- Operating Voltage / Frequency: 100-240V
50/60 Hz

- Current Consumption: 2A max

- Time to full charge: Approximately 2 hours

- Classification: Class II Equipment

Transportation and Storage Conditions

Ambient Temperature Range:	-4°F to 100°F / -20°C to 38°C
Relative Humidity	15% – 85% RH

Warnings and Precautions

 Warnings	
 WARNING	Failure to follow the instructions provided may lead to patient or user harm.
	All curing lights output high levels of light energy and thermal radiation. Failure to follow the instructions provided may lead to soft tissue burns, pulpal damage, ocular damage, or patient discomfort.
	Failure to follow the instructions provided may lead to the device becoming inoperable and cause delay of procedure.
	Failure to properly use barrier sleeves can expose patients or users to cross contamination and risk of infection. Never use this device on a patient without a new barrier sleeve.
	Failure to properly inspect barrier sleeve for damage may lead to an increased risk of cross-contamination.
	Do not activate curing light over soft tissue such as lips, gums, or tongue. Failure to follow this instruction may lead to soft tissue burns.
	Always ensure the provided light shields and proper blue light-filtering safety glasses are utilized during activation of the curing light. Failure to shield patients or users from direct or reflected blue light energy may lead to unsafe levels of viewed light energy, which may be harmful to the eyes.
	Repeated operation of the device will cause the case temperature to increase. Always allow for a cool-down period of 5 minutes if the device is operated for longer than 2 minutes repeatedly. Failure to allow for device cooldown may lead to case temperatures up to 55 °C.
	Always test the selected time and composite material combination using the provided test block prior to utilizing on a patient to verify that the cure time is sufficient to fully cure the material. Failure to test the combination of cure time and material may lead to an increased risk of uncured material and restoration failure.
	Always inspect the output window, including lens cap and barrier sleeve, for debris that may be blocking light output. Any debris blocking light output may significantly reduce the effective irradiance of the device and may lead to restoration failure.
Do not modify this equipment without authorization of the manufacturer.	

 Risk of Electric Shock	
 WARNING	Inspect power cord prior to use for fraying or breakage.
	Examine all contents for damage prior to assembly and use.
	Only handle battery while removed from handpiece with clean, dry hands.
	Do not submerge handpiece, charger base, or power supply in liquid.


 RISK OF EXPLOSION	
 WARNING	Inspect battery for swelling or damage prior to use. Immediately discontinue use of any battery that appears swollen or deformed.
	Do not attempt to charge battery with any device other than the provided charger base or listed replacement component.
	Do not attempt to disassemble battery assembly.
	Do not utilize this device in the presence of flammable anesthetics.
	Do not use the device in an oxygen rich environment.

EMC WARNINGS	
 WARNING	Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
	Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
	Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Coltolux Comfort, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
	This equipment is intended for use by healthcare professionals only. This equipment may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures such as re-orienting or relocating the Coltolux Comfort unit or shielding the location.

Precautions	
CAUTION	Only utilize the accessories provided with the device or listed replacement components.
CAUTION	Do not remove the battery during surface disinfection. Exposure of battery terminals to cleaning agents may lead to device becoming inoperable and delay of procedure.

Unpacking the Device

 WARNING	Examine all contents for damage prior to assembly and use.
---	--

Contents:

QTY	CONTENT
1	Handpiece
1	Battery
1	Charger Base
1	USB-C Power Supply with Plug Adapters
3	Extraoral Light Shield
3	Intraoral Light Shield
20	Lens Caps
100	Barrier Sleeves

Setup and Assembly

 WARNING	<p>Only handle battery while removed from handpiece with clean, dry hands.</p>
	<p>Inspect battery for swelling or damage prior to use. Immediately discontinue use of any battery that appears swollen or deformed.</p>
	<p>Do not attempt to charge battery with any device other than the provided charger base or listed replacement component.</p>
	<p>Failure to recharge battery between uses will degrade the battery functionality and may lead to deep discharge.</p>

After unpacking the device, attach the appropriate plug adapter onto the power supply. Plug the USB-C power supply fully into the wall, and then plug into the charger base. The charger base indicator light ring should briefly light to indicate the charger base is properly connected. Slowly thread battery onto handpiece in clockwise direction until threads fully engage, then tighten until fully seated. Place the handpiece with battery into the charger base.

Ensure battery is fully charged prior to first use. Always place battery on charger base when not in use to ensure optimal battery functionality.

Function Check

 WARNING	<p>Never look directly into the light emitted from the light lens or at reflected light from tooth or other surfaces. Retinal damage may occur.</p>
	<p>Always test the selected time and composite material combination using the provided test block prior to utilizing on a patient to verify that the cure time is sufficient to fully cure the material. Failure to test the combination of cure time and material may lead to an increased risk of uncured material and restoration failure.</p>

The Coltlux® Comfort includes a 20-second cycle timer with audible indication after 10 seconds. The device will automatically deactivate after 20 seconds of operation. The device can be manually deactivated at any time prior to 20 seconds by pressing the “Activation/Deactivation” button.

It is the responsibility of the personnel performing the restorative procedure to verify that the cure time and light output is sufficient to cure the selected material and to the depth desired for the clinical scenario. Prior to utilizing the Coltlux Comfort in a clinical procedure, always utilize the provided test block to verify that the selected composite and cure time will fully cure the dental material to the desired depth. To utilize the test block, place the dental material in the well that corresponds with the desired cure depth for the clinical application. Cure the material for the desired cure time, then remove the composite from the well with a dental spatula. Using the spatula, test the underside of the dental material to verify that the material has fully cured. If the material has not fully cured, ensure the lens cap, barrier sleeve, and LED lens are free from debris that may be blocking the output window and reducing light output. If these components are free of debris, adjust the cure time until the selected dental material fully cures to the depth desired.

Always follow dental material manufacturer’s instructions on proper curing technique of specific materials.

Indicators

Handpiece

1. **Ready Mode** – Indicator LED is solid green.
2. **Sleep Mode** – Indicator LED is off.
3. **Thermal Protection Mode** – Indicator LED is solid red. Alternating high and low audible tones are emitted. During this period, the handpiece will be inoperable, and pressing the handpiece button will not activate the device LED until the handpiece cools.
4. **Low Battery Mode** – Indicator LED flashes green and an audible tone is heard.

NOTE: If low-battery state is reached, place battery in charger base to allow to fully recharge. The Coltlux® Comfort battery may be charged without the handpiece attached.

Charger base

1. **Charging status** – Indicator LEDs are pulsing blue.
2. **Fully charged status** – Indicator LEDs are solid blue.
3. **Charge Error status** – Indicator LEDs are rapidly flashing.

Operating Instructions

 WARNING	<p>All curing lights output high levels of light energy and thermal radiation. Failure to follow the instructions provided may lead to soft tissue burns, pulpal damage, or patient discomfort.</p>
	<p>Do not activate curing light over soft tissue such as lips, gums, or tongue. Failure to follow this instruction may lead to soft tissue burns.</p>
	<p>Always ensure the provided light shields and proper blue light-filtering safety glasses are utilized during activation of the curing light. Failure to shield patient or user from direct or reflected blue light energy may lead to unsafe levels of viewed light energy, which may be harmful to the eyes.</p>
	<p>Repeated operation of the device will cause the case temperature to increase. Always allow for a cool-down period of 5 minutes if the device is operated for longer than 2 minutes repeatedly. Failure to allow for device cooldown may lead to case temperatures up to 55 °C.</p>

1. Ensure the lens cap is fully installed by pressing the lens cap until it snaps into place.
2. Install the intraoral light shield by first snapping onto the neck of the device under the reflector, then sliding up to snap into place around the reflector (see figure below).

3. The intraoral light shield is specially designed to allow placement of the lens flush with the curing plane. The light shield will block reflected blue light from the tooth surface, provided the lens is placed close to the curing plane. Always ensure the intraoral light shield is attached to mitigate the effects of reflected blue light. The intraoral light shield is also designed to fit beneath the single-use barrier sleeve, to provide increased infection control. Always ensure the intraoral light shield is installed beneath the barrier sleeve to provide a proper fit.
4. Examine the lens cap for any debris on the surface. Any debris on the lens cap may reduce the effective light output of the device. If any debris is observed, scrape all debris off or replace the lens cap prior to use.

⚠ WARNING: Always inspect the output window, including lens cap and barrier sleeve, for debris that may be blocking light output. Any debris blocking light output may significantly reduce the effective irradiance of the device and may lead to restoration failure.

CAUTION: Only utilize lens caps provided with this device or listed in these instructions as accessory components.

5. If any damage is observed to the lens cap or light shield, replace prior to use.
6. Install single-use barrier sleeve by sliding the open end of the sleeve over the handpiece. Ensure the seam of the barrier sleeve is not placed over the output window. Install the extraoral light shield by snapping into place in the desired location.
7. Examine the barrier sleeve for damage prior to use on patient. If any damage is observed to the barrier sleeve, immediately replace.

⚠ WARNING: Failure to properly inspect barrier sleeve for damage may lead to an increased risk of cross-contamination.

8. If the indicator LED on the handpiece is turned off, pressing the activation button once will bring the device into ready state. Once the indicator LED turns green, the handpiece is ready for use. Pressing the button will activate the curing LED. The handpiece will run for a 20-second cycle, with an audible tone at 10 seconds and two tones at the end of the cycle.
9. When ready to cure material, ensure the output window is flush with the curing plane. The distance and angle between the output window and the cure plane is critical for successful light curing; therefore, always ensure the output window is placed as close to the material being cured as possible and held steady for the entire desired cure time.

Hygiene tips:


- Always use single-use barrier sleeves and replace after each patient.
- Place intraoral light shield underneath barrier sleeve to minimize risks of cross-contamination.

Cleaning and Disinfection

⚠ WARNING: Failure to properly utilize single-use barrier sleeves can expose patients or users to cross contamination and risk of infection. Never use this device on a patient without a new barrier sleeve.

CAUTION: Do not remove the battery during surface disinfection. Exposure of battery terminals to cleaning and disinfecting agents may lead to device becoming inoperable and delay of procedure.

1. After use, remove the snap-on extraoral light shield from the handpiece.
2. Immediately remove and discard the single-use barrier sleeve. To remove barrier sleeve, grab the sleeve with a gloved hand and pull upwards toward the tip of the device. Promptly dispose of the used barrier sleeve to minimize risk of infectious material transfer.
3. Remove the intraoral light shield from the tip of the device by pulling down from the reflector and away from the neck of the device (see figure below).

4. To remove the lens cap, firmly grab the edges of the lens cap and pull off the face of the reflector. Wipe down the lens cap with 70% isopropyl alcohol and a soft cloth. If clouding of the lens cap is observed, replace prior to subsequent use.
5. Wipe the intraoral and extraoral light shields with an OPTIM-1 wipe until wet. Wipe the surface of the handpiece housing down with an OPTIM-1 wipe until wet. NOTE: Do not allow disinfectant solution to enter the reflector housing and make direct contact with the LED. Always follow the cleaning/disinfectant manufacturer's instructions for use.
6. Allow handpiece assembly and all components to dry fully.
7. Reassemble device following steps 1-8 in Operating Instructions Section.

Manufacturer: Dent4You AG Method: Low-level disinfection Symbol: N/A Device(s): Intraoral light shield, extraoral light shield, lens cap, handpiece assembly	
WARNINGS	Do not apply any liquids directly to the LED inside the reflector, as this may damage the component and lead to device inoperability.
Limitations on Reprocessing	Always inspect handpiece assembly, light shields, and lens caps prior to subsequent use. If damage is observed, do not continue use.

Instructions	
Initial treatment at the point of use	After procedure, remove the extraoral light shield, intraoral light shield, and lens cap from the handpiece assembly.
Cleaning/Disinfection	Wipe handpiece assembly, intraoral light shield, and extraoral light shield with an OPTIM-1 wipe until wet. Wipe lens cap with 70% isopropyl alcohol and a soft cloth. Refer to cleaning/disinfectant manufacturer's Instructions for Use. Thoroughly clean the device.
Drying	Allow all disinfected surfaces to air dry.
Maintenance, Inspection and Testing	Inspect for any contamination on handpiece after cleaning and disinfection. If any contamination is observed after cleaning and disinfection, repeat cleaning and disinfection procedures until no visual contamination is observed. Inspect handpiece assembly and light shields following disinfection and prior to subsequent use. If damage is observed, do not continue use.

Maintenance and Component Replacement

⚠ Warning: Do not modify this equipment without authorization of the manufacturer

All components not explicitly stated in this section as being repairable and replaceable by the end-user require that the device be sent back to the manufacturer for repair or replacement.

1. Battery removal and replacement

- Lithium-ion batteries decrease in charge capacity as they are put through repeated charge/discharge cycles. If you notice your curing light does not maintain the number of full cure cycles per charge expected with a full battery, it may be that the battery is beyond its useful service life and requires replacement.
- If the device is not being used for a long period of time, consider removing the battery from the handpiece to reduce current draw on the battery.
- To remove the battery, firmly grasp the battery and handpiece and twist the battery in a counter-clockwise motion. To attach new battery, slowly thread battery onto handpiece in clockwise direction until threads fully engage, then tighten until fully seated. CAUTION: Do not use excessive force during battery installation or removal; use of excessive force may lead to thread damage and inability to utilize device for subsequent procedures.

2. Lens cap replacement

- To remove the lens cap, firmly grab the edges of the lens cap and pull off the face of the reflector. Press the new lens cap so it is fully seated on the reflector. Replace the lens cap any time it is visibly cloudy or covered in debris.

3. Extraoral light shield replacement

- To remove the extraoral light shield, pull away from the handpiece. To attach the new extraoral light shield, press the shield perpendicular to the handpiece until it snaps into place.

4. Intraoral light shield replacement

- To remove the intraoral light shield, slide down the neck until it detaches from the reflector and then pull the shield away from the neck.

To replace, first snap the shield onto the neck below the lens, then slide upwards to snap into place around the reflector.

5. Power supply replacement

- To remove the power supply, unplug from the charger base. Plug in the new power supply fully. Ensure used power supply is discarded according to applicable local regulations.

Problem	Issue/Root Cause	Solution
The curing light indicator LED does not light green when the button is pressed.	Handpiece is on the charger base.	Ensure the handpiece is removed from the charger before attempting to turn on.
	Battery is depleted.	Place the battery in the charger base and allow to charge.
	Battery is not fully seated on the handpiece.	Twist the battery onto the handpiece clockwise fully.
The curing light indicator LED is lighting red when the activation button is pressed and unit is beeping.	Handpiece is in thermal protection mode.	Allow the unit to cool down for a period of 5 minutes prior to continuing use. Handpiece indicator LED will turn green to indicate the device is ready to continue use.
	Handpiece has deactivated due to a component error.	Contact COLTENE Customer Service.
The curing light indicator LED is blinking green and an audible tone is heard.	Battery is low.	Place the battery in the charger base and allow to charge fully.
The charger base LED indicator is blinking rapidly.	There is an error with the wireless charging system.	Remove the battery from the charger base. Ensure there are no foreign objects on the outside of the battery or in the charger base that may be interrupting the charging system. If problem persists, contact COLTENE Customer Service.
The curing light is emitting lower light output than expected.	Debris is blocking the output window of the lens.	Remove foreign debris or replace lens cap.
The curing light does not output the rated number of curing cycles in a full battery charge.	Battery has exceeded its operative life.	Replace battery.
The charger base LED indicator does not light up when plugged into the wall.	The power supply is not fully plugged in.	Ensure the power supply is plugged in properly to the wall and the charger base.
	Charger base is deactivated due to a component error.	Contact COLTENE customer service.
The charger base LED indicator does not light up when battery is placed inside the base.	The charger base is not plugged into the power supply.	Plug the power supply into the charger base.
	Debris is between charger base and battery.	Remove debris from the charger base.
	The battery is deactivated due to a component error.	Contact COLTENE customer service for a replacement battery.

NOTICE: Report all incidents related to the use of this device to the manufacturer and the member state in which they occurred.

Instructions for Safe Disposal

	<p>This product must not be disposed of with other waste. Instead, it is the user's responsibility to dispose of their waste equipment by handing it over to a designated collection point for the recycling of waste electrical and electronic equipment. The separate collection and recycling of your waste equipment at the time of disposal will help to conserve natural resources and ensure that it is recycled in a manner that protects human health and the environment. For more information about where you can drop off your waste equipment for recycling, please contact your local city office, your waste disposal service or where you purchased this product.</p>
--	---

Explanation of Symbols

	Caution/Warning		Do not look directly into light emitted from LED
	Type B Applied part		Temperature Limitation
	CE Symbol		Humidity Limitation
	EU Authorized Representative		Keep Dry

	Refer to Instruction Manual		Recycle
	Electric Shock Hazard		Consult the instructions for use
	Explosion Hazard		Fragile – handle with care
MD	Medical Device Indicator: this product meets the definition of a medical device	REF	Catalog Number
Rx ONLY	Sale Restricted: The purchase of the device is restricted to licensed healthcare practitioners only.		The entity that manufactures the medical device
	Single-use only (barrier sleeves)		Class II Equipment
	Date of manufacture		Non-ionizing radiation
SN	Serial Number	UK CA	UK Conformity Assessed
LOT	Lot Number (accessories only)		Do not dispose of in household waste

Order Information

Catalog Number	Description
60034390	Coltolux® Comfort Kit, International
60034391	Coltolux® Comfort Spare Battery
60034392	Coltolux® Comfort Universal Power Supply
60034393	Coltolux® Comfort Barrier Sleeves (Pack of 1,000)
60034394	Coltolux® Comfort Intraoral Light Shields (Pack of 5)
60034395	Coltolux® Comfort Extraoral Light Shields (Pack of 5)
60034396	Coltolux® Comfort Lens Caps (Pack of 20)
60034446	Coltolux® Comfort Composite Test Block

Warranty Information

For a period of 2 years from date of purchase, the legal manufacturer, Dent4You AG, guarantees that the Coltolux Comfort LED curing light, when manufactured in new and unused condition, will operate as intended and within the specifications outlined in this user manual. The legal manufacturer also guarantees working condition of the battery supplied with the device for a period of 1 year from date of purchase.

Electromagnetic Compatibility Information, Warnings, and Considerations

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Coltolux Comfort, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

NOTE

The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

EN/IEC 60601 - 1 -2 2014		
Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
The Coltolux Comfort is intended for use in the electromagnetic environment specified below. The customer or the end user of the Coltolux Comfort should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11:ed5.0 (with A1:2010)	Group 1	The Coltolux Comfort uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11: ed5.0 (with A1:2010)	Class A	The Coltolux Comfort unit is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings for domestic purposes, provided the following warning is heeded: WARNING: This equipment is intended for use by healthcare professionals only. This equipment may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures such as re-orienting or relocating the Coltolux Comfort unit or shielding the location.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/Flicker emissions IEC 61000-3-3	Complies	

EN/IEC 60601 - 1 -2 2014

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The Coltolux Comfort unit is intended for use in the electromagnetic environment specified below. The customer or the end user of the Coltolux Comfort unit should assure it is used only in such an environment.

Immunity Test	IEC60601 test level	Compliance Level	Intended Electromagnetic Environment
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV Contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	±8kV Contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines 100 kHz repetition frequency	±2kV for power supply lines 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV differential mode (line-line) ±2kV common mode (line-earth)	±1kV differential mode (line-line) ±2kV common mode (line-earth)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Coltolux Comfort unit requires continued operation during power mains interruptions, it is recommended that the Coltolux Comfort unit be powered from an uninterruptible power supply with sufficient capacity to run the unit for the maximum required time of interruption.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30A/m 50/60 Hz	30A/m 50/60 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note UT is the a.c. mains voltage prior to application of the test level.

EN/IEC 60601 - 1 -2:2014

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The Coltolux Comfort unit is intended for use in the electromagnetic environment specified below. The customer or the end user of the Coltolux Comfort unit should ensure it is used in such an environment.

Immunity Test	IEC60601 test level	Compliance Level	Intended Electromagnetic Environment
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms</p> <p>6Vrms (In ISM Bands) 150 kHz to 80 MHz</p> <p>80MHz to 2.7GHz</p>	<p>3Vrms</p> <p>6Vrms (In ISM Bands) 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Coltolux Comfort Curing Light, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p>$d = 1.2\sqrt{P}$</p> <p>$d = 1.2\sqrt{P}$ 80MHz to 800 MHz</p> <p>$d = 2.3\sqrt{P}$ 800MHz to 2.5GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended minimum separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> <div style="text-align: center;">  <p>Non-Ionizing Radiation</p> </div>

NOTE 1: At 80MHz and 800MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from objects, structures and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Coltolux Comfort Curing Light is used exceeds the applicable RF compliance level above, the Coltolux Comfort Curing Light should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Coltolux Comfort Curing Light.

b. Over the frequency range 150kHz to 80MHz, field strengths should be less than 3V/m.

EN/IEC 60601 - 1-2:2014			
Recommended separation distances between portable and mobile RF communications equipment and the Coltolux Comfort unit			
The Coltolux Comfort unit is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Coltolux Comfort unit can help prevent electromagnetic interference by maintaining a minimum distance between the portable and mobile RF communications equipment (transmitters) and the Coltolux Comfort unit as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter in watts (W)	Separation distance according to frequency of transmitter in meters (m)		
	150kHz to 80MHz $d = 1.2\sqrt{P}$	80MHz to 800MHz $d = 1.2\sqrt{P}$	800MHz to 2.5GHz $d = 2.3\sqrt{P}$
0.01	.12	.12	.23
0.1	.38	.38	.73
1.0	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Table 9 – Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment						
Test frequency (MHz)	Band ^{a)}	Service ^{a)}	Modulation ^{b)}	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380 – 390	TETRA 400	Pulse modulation ^{b)} 18 Hz	1.8	0.3	27
450	430 – 470	GMRS 460, FRS 460	FM ^{c)} ± 5 kHz deviation 1 kHz sine	2	0.3	28
710	704 – 787	LTE Band 13, 17	Pulse Modulation ^{b)} 217 Hz	0.2	0.3	9
745						
780						
810	800 – 960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse Modulation ^{b)} 18 Hz	2	0.3	28
870						
930						
1 720	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse Modulation ^{b)} 217 Hz	2	0.3	28
1 845						
1 970						
2 450	2 400 – 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse Modulation ^{b)} 217 Hz	2	0.3	28
5 240	5 100 – 5 800	WLAN 802.11 a/n	Pulse Modulation ^{b)} 217 Hz	0.2	0.3	9
5 500						
5 785						

NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

a) For some services, only the uplink frequencies are included.

b) The carrier shall be modulated using a 50 % duty cycle square wave signal.

c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

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Coltene/Whaledent Inc.
Coltène/Whaledent GmbH + Co. KG
Coltène/Whaledent Ltd.