4DESIGN 4DISKS

DIGITAL PRECISION SOLUTIONS ZIRCONIA

- DE** Zirkonoxid-Disc. Gebrauchsanweisung.
- EN** Zirconia disc. Instructions for use.
- ES** Disco de zirconio. Instrucciones de uso.
- FR** Disque de zirconio. Mode d'emploi.
- IT** Disco di zirconio. Istruzioni per l'uso.
- PT** Disco de zircónio. Instruções de uso.

DE ZIRKONOXID-DISC

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HT WHITE	ST COLOR	ST MULTILAYER	TT MULTILAYER
ZrO ₂ +HfO ₂ +Y ₂ O ₃ 99% Y ₂ O ₃ 4.5-6% Al ₂ O ₃ <0.5% Andere Oxide <0.5%	Nanometer Zirkonium Pulver>98% Fe ₂ O ₃ <0.3% Pr ₂ O ₃ <0.2% Er ₂ O ₃ <1.0% Andere Oxide <0.5%	Nanometer Zirkonium Pulver>98% Fe ₂ O ₃ <0.3% Pr ₂ O ₃ <0.2% Er ₂ O ₃ <1.0% Andere Oxide <0.5%	ZrO ₂ +HfO ₂ 86.3-94.2% Y ₂ O ₃ 5.8-9.7% Er ₂ O ₃ <2% Fe ₂ O ₃ <0.5% Al ₂ O ₃ <0.5% Andere Oxide <0.5%
WAK (25-500°C)			
(10.5±1.0) x 10 ⁻⁶ K ⁻¹	(10.5±1.0) x 10 ⁻⁶ K ⁻¹	(10.5±1.0) x 10 ⁻⁶ K ⁻¹	(10.5±1.0) x 10 ⁻⁶ K ⁻¹

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Instruction for use

This product must only be handled, fitted and adjusted by dental professionals.

Used for dental restoration using different CAD/CAM or manual milling machines. All Blanks must be processed in dental laboratories by dental professionals.

1. Fix the zirconia blank in the specific clamping device.
2. Start the milling process.
3. Remove milled framework from blank and pre-trim with diamond disc.
4. Remove dust from the framework with an air blaster.
5. Load into the furnace.
6. Select sintering program and start sintering (see sinter instructions).
7. Remove the sintered framework from the furnace after the furnace temperature has dropped to room temperature.
8. Check the frame on the working model for sitting and margin fit. Trim with grinder and polisher if required.
9. Blast using a 50 micron Alumina at 2-2.5 bars of pressure. Clean the frame with a brush and water. Do not steam-clean.
10. The frame is now ready for veneering.
11. Now the restoration is ready for trying on the patient.

NOTE: The product must only be used by trained personnel. The manufacturer is not liable for any undesired results, since we have no control if this instruction is properly followed or if the product is appropriately used. In the event of a technical query which has not previously been described, please contact our technical service: info@4desingdigital.com

Medical Device Hazards: There are no medical device warnings for this product.

Expiry: The product should be replaced every 5 years. The product is not expected to degrade during this time period.

Storage and Handling: Keep the container dry and at a stable room temperature.

Safety instructions:

1. Avoid collision with hard material. Avoid squeezing or shaking.
2. Do not apply to the patient before the material is fully sintered.
3. This product should only be handled by dental professionals.
4. Wear a suitable face mask to prevent fine dust inhalation. Avoid fine dust air pollution by using the air blaster in short bursts only rather than a continuous flow.

Disposal: The material is not hazardous to the environment.

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