- EN INSTRUCTIONS FOR USE
- FR MODE D'EMPLOI
- IT ISTRUZIONI PER L'USO
- DE GEBRAUCHSANWEISUNG
- ES INSTRUCCIONES DE USO
- PT INSTRUÇÕES DE UTILIZAÇÃO
- NL GEBRUIKSAANWIJZING
- SV BRUKSANVISNING
- DA BRUGSANVISNING
- NO BRUKSANVISNING
- FI KÄYTTÖOPAS
- PL INSTRUKCJA UŻYCIA
- LT NAUDOJIMO INSTRUKCIJA
- LV LIETOŠANAS INSTRUKCIJA
- ET KASUTUSJUHEND
- CZ NÁVOD K POUŽITÍ
- SK NÁVOD NA POUŽITIE
- HU HASZNÁLATI UTASÍTÁS
- RO INSTRUCȚIUNI DE UTILIZARE
- SR UPUTSTVO ZA UPOTREBU
- CNR UPUTSTVA ZA UPOTREBU
- EL ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
- RU ИНСТРУКЦИИ ПО ПРИМЕНЕНИЮ
- TU KULLANIM TALIMATLARI
- AR تعليمات الاستخدام

## Instructions for Use – ENGLISH

**CAUTION:** This is a medical device. For dental use only.

**USA:** Rx only.

### 1. PRODUCT DESCRIPTION

CEREC® Zirconia+ CAD/CAM block for CEREC® and inLab® is a material that can be used for the fabrication of fully anatomic crowns and bridges in the anterior and posterior region using a CAD/CAM procedure. Indirect restorations are fabricated by milling CEREC® Zirconia+ CAD/CAM blocks using a Dentsply Sirona CAD/CAM system.

Note that CEREC® Version 5.1.3 with Material Pack or higher or inLab® CAD 20.0.3 Service Pack and inLab® CAM 20.0.1 Material Pack or higher is required. After milling, the restorations are sintered in the Dentsply Sirona CEREC® SpeedFire furnace or another conventional sintering furnace, such as inFire HTC speed or inLab® ProFire. The sintered restorations are finalized by polishing or polishing and glazing. CEREC® Zirconia+ dental ceramic material type II, class 5 pursuant to ISO standard 6872:2015 + Amd.1:2018.

#### 1.1 Intended purpose

Ceramics for fixed dental prosthetic restorations.

#### 1.2 Indications

CEREC® Zirconia+ CAD/CAM blocks are indicated for all ceramic restorations for anterior and posterior locations:

- Fully anatomic crowns and bridges
- Bridges with maximum two pontics

#### 1.3 Intended User

Dental professionals (dental technicians and dentists).

#### 1.4 Intended patient population and medical conditions

CEREC® Zirconia+ CAD/CAM blocks are intended for patients in need of long-term restorative or prosthodontic dental therapy, or esthetic dental corrections. The use of ceramics is not limited to a specific patient population.

#### 1.5 Contraindications

- Insufficient oral hygiene
- Insufficient tooth structure
- Insufficient preparation results
- Insufficient space available

#### 1.6 Composition

Oxides	Unit	Concentration
Y <sub>2</sub> O <sub>3</sub>	wt.-%	< 8
HfO <sub>2</sub>	wt.-%	< 3
Al <sub>2</sub> O <sub>3</sub> , SiO <sub>2</sub> + shading pigments	wt.-%	< 2
The content of ZrO <sub>2</sub> is calculated value = 100 % - (HfO <sub>2</sub> + Y <sub>2</sub> O <sub>3</sub> + Al <sub>2</sub> O <sub>3</sub> + SiO <sub>2</sub> + shading pigments)		

#### 1.7 Compatible Stains and Glazes and veneering ceramics

The use of spray glaze or paint-on stains or glaze is optional with CEREC® Zirconia+ CAD/CAM blocks. CEREC® Zirconia+ CAD/CAM blocks are compatible with Dentsply Sirona Universal Stain and Glaze System, Dentsply Sirona Universal Spray Glaze and Dentsply Sirona Universal Spray Glaze Fluo.

CEREC® Zirconia+ can be veneered with Zirconia veneering ceramics e.g. Celtra® Ceram and Cercon® ceram veneering ceramic. See corresponding IFUs for processing instructions.

#### 1.8 Compatible Cements

CEREC® Zirconia+ CAD/CAM block restorations are compatible with Universal/Self-Adhesive Cement, Adhesive Resin Cement, and Conventional Cement systems, including all Dentsply Sirona cement systems (Calibra® Ceram, Calibra® Universal and Calibra® Bio) designed for Zirconia ceramic cementation (see complete Instructions for Use of selected cement).

#### 1.9 Technical data

The following specifications apply to CEREC® Zirconia+ restorations that have been sintered in a CEREC® SpeedFire, a inFire HTC speed, or a inLab® ProFire sintering furnace.  
Coefficient of thermal expansion (20 – 500 °C): 10.3 · 10<sup>-6</sup> K<sup>-1</sup>  
Flexural strength (3-point flexural strength): > 1000 MPa

## 2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these instructions for use.

### Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

### 2.1 Warnings

- If properly processed and used, adverse effects from this medical device are highly unlikely. However, reactions of the immune system (such as allergies) or localized paresthesia (such as an irritating taste or irritation of the oral mucosa) cannot be completely ruled out. In case of skin sensitization or rash, discontinue use and seek medical attention.
- CEREC® Zirconia+ CAD/CAM block restorations are not suitable for patients with clinical symptoms of parafunctional habits or bruxism.
- Do not inhale dust particles during milling or finishing. Wear a suitable protective mask.
- Do not use competitive materials except as indicated in 1.4 as it may impact the performance of CEREC® Zirconia+.
- Choosing the wrong glazing program may lead to damage of the restoration or CEREC® SpeedFire.
- Do not apply glaze to the surface of sintered CEREC® Zirconia+ restorations in an unpolished state. Excessive wear of opposing surfaces may occur. The restoration must be polished prior to glaze application (see step-by-step Instructions).
- In patients with hypersensitivity to any of the ingredients, this medical device may not be used at all or only under the advice of the dentist or physician in charge.

### 2.2 Precautions

- This product is intended to be used only as specifically outlined in these Instructions for Use. Any use of this product inconsistent with the Instructions for Use is at the discretion and is the sole responsibility of the practitioner.
- Wear suitable protective eyewear, clothing and gloves. Protective eyewear is recommended for patients.
- Contamination of the preparation or margin area with saliva, blood, water, or hemostatic agents during adhesive cementation may lead to an adhesive failure. Ensure adequate isolation and tissue management techniques during adhesive cementation.
- Devices marked “single use” on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- CEREC® Zirconia+ CAD/CAM block restorations require adequate preparation reduction and restoration thickness. Insufficient wall thickness may lead to premature failure.
- CEREC® Zirconia+ CAD/CAM block restorations must be sintered and must be polished or polished and glazed before insertion. It is recommended to polish the occlusal areas before glazing. Direct insertion without sintering lead to failure.
- Use only in well ventilated areas.
- Do not cement using provisional cements. Use of provisional/temporary cements could cause fracture of CEREC® Zirconia+ restorations.
- Insufficient data exists to support use of CEREC® Zirconia+ CAD/CAM blocks to fabricate partial coverage restorations (veneers, inlays, onlays), resin-bonded-retainer bridges (“Maryland” bridges), endodontic post and cores, or implant abutments.
- Avoid adjusting sintered restorations with milling tools, especially in the connector area. Flexural strength may be compromised (see step-by-step Instructions).
- Use of Dentsply Sirona Universal Spray Glaze Fluo with CEREC® Zirconia+ restorations may result in a brighter shade appearance. Use of Dentsply Sirona Universal Spray Glaze is recommended if this effect is not desired.
- Conventional sintering may lead to higher chroma.

### 2.3 Interactions

CEREC® Zirconia+ CAD/CAM blocks are designed to be fabricated using a Dentsply Sirona CAD/CAM system. Milling blocks using non-compatible CAD/CAM systems may lead to inadequate or unacceptable restorations.

### 2.4 Adverse reactions

No adverse reactions have been reported for CEREC® Zirconia+ CAD/CAM blocks. Should you hear or receive information about any adverse effects, notify Dentsply Sirona.

### 2.5 Storage conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product. Store in a dry place and protect from moisture. Do not use after the expiration date.

### 3. STEP-BY-STEP INSTRUCTIONS

#### 3.1 Preparation

Proper reduction of the tooth during preparation is essential for maximizing the strength, shade and retention of the finished restoration.

When preparing anterior or posterior teeth, the anatomical form must be reduced as shown.

It is recommended to prepare a taper which is between 4° and 8°. The preparation must be performed with either a chamfer or a shoulder with rounded internal angle.


**Preparation guidelines for crowns:** All internal line angles of a preparation should be rounded.

**Minimum wall thickness for crowns:** The following pictures show the specified minimum wall thickness for crowns.


The minimum wall thickness must still be ensured after all manual adjustments have been made.

**NOTE:** The standard min wall thickness in the Software is 0.7 mm. The minimal wall thickness can be reduced to 0.5 mm for crowns except for the Super Fast Milling Mode.

Min wall thickness for Super Fast Milling Mode.

Minimum wall Thickness for all milling modes except of the Super Fast Milling Mode.


#### Preparation guidelines for abutment teeth (anterior and posterior bridges)

The preparation guidelines for abutment teeth correspond to the guidelines for the above described crowns.

The reduction should be at least 0.6 mm at the lowest point of the main fissure, cusps and axial wall.

Bridge design can be conventional, single abutment crowns at each end, with no more than 2 pontics.

The cross sections of the connectors can be taken from the table below:

Connector cross section for	Minimum connector cross section in mm <sup>2</sup>
Anterior bridge restoration	9
Posterior bridge restoration	12

#### 3.2 CAD/CAM processing

CEREC® Zirconia+ restorations are produced with CEREC® or inLab® CAD/CAM systems by Dentsply Sirona. If you have any questions about these systems, please contact Dentsply Sirona.

##### 3.2.1 Software Requirements

CEREC® Zirconia+ CAD/CAM blocks are supported by CEREC® version 5.1.3 with Material Pack or higher or inLab® CAD 20.0.3 Service Pack and inLab® CAM 20.0.1 Material Pack or higher.

##### 3.2.2 Milling the restoration

In the CEREC® software, select CEREC® Zirconia+ CAD/CAM Block material. The Dentsply Sirona CAD/CAM milling unit will prompt you to insert a CEREC® Zirconia+ CAD/CAM Block. For detailed processing, please consult the Instructions for Use and technical manuals of the appropriate CAD/CAM systems. Make sure to follow the manufacturers' recommendations.

For CEREC® Zirconia+ CAD/CAM Blocks dry milling is recommended, but, if desired, wet milling is also possible.

#### Wet Milling - To avoid contamination and reduced translucency

Using wet milling, you have to be aware that contaminated cooling water (e.g. due to residuals of glass-ceramic particles) can diminish the translucency of the final restoration. Therefore, renewal of the cooling water and cleaning of the milling chamber and filter is strictly recommended before wet milling. Alternatively, three separate water tanks for glass ceramics, rinsing in-between, and zirconia wet milling could be used.

#### 3.3 Sintering

Classic and speed sintering is possible.

##### 3.3.1 Preparation for sintering

After the milling process (and prior to sintering) a tungsten carbide tool shall be used to separate the restoration from the blank. Remove the sprue prior sintering.

If further adjustments are needed, it is recommended to do the adjustments before sintering.

#### Airborne particles - To reduce inhalation risk

- Do not inhale abrasive dusts.
- Use a vacuum system and wear a mask.

##### 3.3.2 Sintering the restoration

A Dentsply Sirona Sintering furnace (CEREC® SpeedFire, inFire HTC speed, or inLab® ProFire) is recommended for the sintering process. Before sintering the restoration, it is recommended to free the restorations from dust using compressed air or a ceramic brush.

Do not inhale abrasive dusts. Use a vacuum system and wear a mask.

**NOTE:** Restorations made from CEREC® Zirconia+ must be sintered in dry conditions. Therefore, wet milled restorations have pre-drying steps included in their sintering programs, prolonging overall sintering time.

#### Sintering in the CEREC® SpeedFire

When sintering CEREC® Zirconia+ restorations in the CEREC® SpeedFire, the CEREC® software automatically transfers the job to the CEREC® SpeedFire if the milling machine and the CEREC® SpeedFire are connected.

With CEREC® SW 5.2.3 or inLab CAM 22.0.0 Multi Job sintering of up to 3 single restorations is possible. Select the jobs of the single restorations and the furnace will adjust the sintering time automatically.

**NOTE:** Maximum restoration size

Observe the maximum furnace chamber size when loading the Furnace:

- Diameter: 38 mm
- Height: 20 mm

The restoration (including the Glazing Support Single/Multi Unit) must not exceed the size of the chamber, length (38 mm) and height (20 mm), otherwise the chamber may be damaged.

For further information please refer to the Operating Instructions of the CEREC® SpeedFire

-  Place the restoration with the occlusal surface facing down directly on the top door insulation.
- Start the process by touching the start icon. The furnace shuts automatically once the process starts.
- The furnace opens automatically following successful heat treatment. The process is not yet complete, as a cooling-off phase occurs in an open condition. There will be a signal tone once the cooling process has completed. When the unit's LED status display is green, the furnace can be unloaded after approximately 30 seconds.

#### Risk of injury

The restoration and parts of the door insulation may still be hot when the unit's LED status display is green. Always use tweezers to remove the restoration from the furnace. Allow the restoration to cool down for another five minutes before picking it up with hands.

- Only unload the furnace using metallic or ceramic tweezers. Plastic tweezers are not suitable for removing restorations, as the restoration is still very hot in this phase. Place the restoration on the fireproof tray for further cooling.

**NOTE:** Sintering times of dry-milled restorations can be reduced by placing the restoration in a 400°C preheated CEREC® Speedfire. For this, use the preheating function of CEREC® Speedfire before starting the sintering process.

#### inFire HTC speed:

##### Dry milled restorations

When sintering in the inFire HTC speed use the pre-programmed program "inCoris TZI/ZI speed".

	Heating rate °C/min	Holding temperature/°C	Holding time/min
S4	99	750	0
S3	99	1510	0
S2	50	1510	30
S1	99	800	5

### Wet milled restorations

When sintering in the inFire HTC speed use the pre-programmed program "inCoris TZI/ZI speed wet".

	Heating rate °C/min	Holding temperature/°C	Holding time/min
S4	99	750	0
S3	50	1510	30
S2	99	800	5
S1	15	80	30

### inLab® ProFire:

#### Dry milled restorations

When sintering in the inLab® Profire use the pre-programmed program "CEREC® Zirconia+ speed".

	Heating rate °C/min	Holding temperature/°C	Holding time/min
S1	99	800	5
S2	50	1510	30
S3	99	750	0

### Wet milled restorations

When sintering in the inLab® Profire use the pre-programmed program "CEREC® Zirconia+ speed wet".

	Heating rate °C/min	Holding temperature/°C	Holding time/min
S1	15	80	30
S2	99	800	5
S3	50	1510	30
S4	99	750	0

### Sintering CEREC® Zirconia+ with other furnaces

For sintering the restorations in furnaces from other manufacturers, follow the firing parameters outlined in Tables below or visit [www.dentsplysirona.com](http://www.dentsplysirona.com)

### Dry milled restorations

Start	Heating rate 1	Temp. 1	Holding time 1	Heating rate 2	Temp. 2	Holding time 2	Cooling Rate	Temp. 3	Action
°C	°C/min	°C	min	°C/min	°C	min	°C/min	°C	Open door
RT	99	800	5	50	1510	30	99	750	

### Wet milled restorations

Start	Heating rate 1	Temp. 1	Holding time 1	Heating rate 2	Temp. 2	Holding time 2	Heating rate 3	Temp. 3	Holding time 3	Cooling Rate	Temp. 4	Action
°C	°C/min	°C	min	°C/min	°C	min	°C/min	°C	min	°C/min	°C	Open door
RT	15	80	30	99	800	5	50	1510	30	99	750	

### 3.3.3 Adjustments of sintered restorations

#### Surface Condition – To reduce risk of compromised flexural strength

The surface condition of ceramic materials is critical for their flexural strength. Adjustment of the sintered restorations with milling tools, especially in the connector area, must be avoided.

However, if adjustment is necessary, then follow these basic rules:

- Adjustment in the sintered state should be performed with fine diamonds in a highspeed handpiece with water cooling and with low pressure. Adjustments made with diamonds instruments must be followed by polishing.
- As an alternative it is possible to do adjustments with soft diamond rubber polishers and a handpiece at low speed and low pressure. The tool must be applied flat to minimize the chatter.
- Areas that are under tension in clinical use, i.e. primarily the connectors in bridge structures, should not be adjusted.

### 3.4 Polishing, Try-in and optional Staining and Glazing

CEREC® Zirconia+ restorations can either be polished, or polished and glazed. It is recommended to polish the occlusal areas before glazing.

#### 3.4.1 Polishing

- CEREC® Zirconia+ can be polished with standard polishing agents for zirconia ceramics.
- Subsequent heat treatment (depressurization fire) is not necessary/recommended.
- The MEISINGER LUSTER® Kits for Zirconia are recommended for polishing the restoration.

1. 9735H: Trimming the contact points, Recommended rotary Speed: 8.000 – 12.000 rpm
2. DCA06: Smoothing of the external shape, Recommended rotary Speed: 7.000 – 12.000 rpm

3. 9771M\*: Polishing of the occlusal surfaces, Recommended rotary Speed: 7.000 – 12.000 rpm
4. DCA12: High-shine polishing of the external shapes, Recommended rotary Speed: 7.000 – 12.000 rpm
5. 9771C\*: High-shine polishing of the occlusal surfaces, Recommended rotary Speed: 7.000 – 12.000 rpm

\*Note: To be used with little contact pressure!

For further information please refer Instructions for Use of the MEISINGER LUSTER® Kits. Meisinger is not part of Dentsply Sirona.

#### 3.4.2 Optional Try-in

- Try-in the restoration for marginal and proximal fit. Make any necessary adjustments as outlined above.

#### Contamination – To reduce risk of infection

Restorations should be polished, cleaned and disinfected before and after optional try-in. See Hygiene section below.

- Adjustments made during try-in must be re-polished as outlined below prior to optional stain and glaze application and final delivery.

#### 3.4.3 Optional Staining and Glazing

#### Surface Condition – To reduce risk of excessive wear

- Prior to glazing, ensure restorations are properly sintered and at least the occlusal areas are polished.
- Restorations must be clean and dry before application.

CEREC® Zirconia+ restorations can be glazed with: Dentsply Sirona Universal Stain and Glaze, Dentsply Sirona Universal Spray Glaze, or Dentsply Sirona Universal Spray Glaze Fluo.

First, follow the polishing instructions from Section 3.4.1. After polishing, clean the restoration using either an ultrasonic cleaner or steam cleaner, prior to Spray Glaze or paint-on application. Ensure the restoration is free of contamination and completely dried with oil-free air prior to paint-on or Spray Glaze application.

#### Preparation of the restoration for Staining and/or Glazing

Use the Glazing Support Single/Multi Unit holders and CEREC® SpeedPaste. Add a small portion of CEREC® SpeedPaste to the restoration to ensure a safe hold on the pin. Do not fill the entire restoration with the firing paste. Place the Glazing Support Single/Multi Unit holders in the paste and pick up the tweezers, if required. Observe the operating instructions supplied with the CEREC® SpeedPaste.

#### Glazing with Dentsply Sirona Universal Spray Glaze or Dentsply Sirona Universal Spray Glaze Fluo

1. Pick up the support with the restoration using tweezers.
2. Shake the spray can vigorously immediately before applying.
3. Maintain a distance of 6 – 10 cm (2.5 – 4.0 in) between the nozzle outlet and the restoration surface.
4. Spray on an even coat of Dentsply Sirona Universal Spray Glaze or Dentsply Sirona Universal Spray Glaze Fluo from all sides directly onto the (stained) restoration. Shake the spray can vigorously between individual sprays.
5. Apply a uniform layer glaze to the surface of the restoration. The applied glaze should be thin and uniform.
6. Ensure the intaglio surface of the restoration is free of glaze. If the glazing material is present on the intaglio surface of the restoration, remove it with a dry, stiff, short-bristled brush.
7. Keep the can as upright as possible during the spraying process.
8. Wait for a few seconds until the glaze dried and displays an even, whitish coating. Apply more spray as required.
9. If another application of spray is necessary after the firing, the spraying and firing can be performed once more in the same way.

Refer to Instructions for Use (IFU) of Dentsply Sirona Universal Spray Glaze/Dentsply Sirona Universal Spray Glaze Fluo for further details.

#### NOTE:

- If characterization is desired, the stains must be applied prior to application of Spray Glaze.
- The Use of Dentsply Sirona Universal Spray Glaze Fluo may lead to a brighter shade. Use Dentsply Sirona Universal Spray Glaze if this effect is not desired (see complete Instructions for Use).

#### Glazing with DS Universal Stain and Glaze

1. Withdraw desired amount of DS Universal Glaze from the jar and place it on the mixing palette.

**NOTE:** If the stain or glaze inside the jar has separated, mix thoroughly with a glass, plastic or zirconia spatula.

2. If a thinner consistency is desired, dilute the material with the Dentsply Sirona Stain and Glaze Liquid (REF 601315/15 ml, 601350/50 ml).
3. Apply a thin layer of the glaze material all over the crown surface in the usual manner using a brush. Make sure not to apply the glaze material too thick or too thin. Too thin of a layer will result in a more matte finish. Too thick layers lead to bubbles on the surface.

4. If a more intensive shade effect is desired, DS Universal Stains can be applied on the glazed surface of the restoration. Apply a thin layer of the stain material on the crown surface in the usual manner using a brush. Make sure not to apply the glaze material too thick or too thin.

Refer to Instructions for Use (IFU) of Universal Stain & Glaze for further details.

#### 3.4.4 Glazing the Restoration in the furnace

**NOTE:** A restoration can be glazed a max of two times in the furnace.

#### Glazing CEREC® Zirconia+ in the CEREC® SpeedFire

Glaze restorations individually and do not glaze multiple restorations at the same time. Position the Support with the restoration centrally on the top door insulation and ensure that the Glazing Support Single/Multi Unit or restoration do not protrude out from the door insulation; otherwise, these may collide with the furnace chamber.

Two different glazing programs are available. Select the "GLAZING (Spray)" program for use of DS Universal Spray Glaze and DS Universal Spray Glaze Fluo or "GLAZING (Stain & Glaze Paste)" program for use of DS Universal Stain and Glaze on the CEREC® SpeedFire control panel and allow the unit to cycle.

**NOTE:** Using DS Universal Stain & Glaze for crowns and bridges preheating to 400 °C is recommended.

#### Choose the right Glazing program

- Choosing the wrong glazing program may lead to damage of the restoration or CEREC® SpeedFire.

For further information please refer to the Operating Instructions of the CEREC® SpeedFire.

#### Glazing CEREC® Zirconia+ with other furnaces

For firing the restoration in furnaces from other manufacturers, follow the firing parameters outlined in Table below.

Drying	Closing	Pre-heat- ing tem- perature	Pre- heating	Heating rate	Final tempera- ture	Vacuum	Holding time	Cooling
min	min	°C	min	°C/min	°C	min	min	min
0	2	400	0	55	760	0	2:00	0

#### 3.5 Pre-Cementation Surface Preparation

- Sand-blast the internal surface of the restoration using 50 µm aluminum oxide at a max pressure of 2.5 bar (35 PSI).
- Clean and disinfect the restoration as outlined in Hygiene section below.

#### Contamination – To reduce risk of infection

Restorations should be polished, stained and glazed (optional, if desired) and sandblasted internally before cleaning and disinfection.

Restorations should be cleaned and disinfected immediately prior to delivery. Follow instructions in Hygiene section below before delivery.

#### 3.6 Cementation

- Strictly follow preparation guidelines described above regarding preparation angle of 4 – 8° for conventional and self-adhesive cementation.
- Reduction of at least 0.5 – 0.7 mm in the central fossa and along the axial walls is mandatory for all types of cementation.

**NOTE:** If contact points need to be adjusted after cementation, the adjusted areas need to be polished afterwards.

#### 3.6.1 Conventional Cementation (Full Coverage Crowns and Bridges)

Dentsply Sirona Calibra® Bio Bioceramic Luting Cement is recommended for conventional cementation of full coverage crowns and bridges fabricated from CEREC® Zirconia+ (see full instruction for use). If using conventional cements (resin modified glass ionomer (RMGI) or glass ionomer (GI) type cements) from other manufacturers, follow their respective instructions for use.

#### 3.6.2 Self-Adhesive Cementation (Full Coverage Crowns and Bridges)

Dentsply Sirona Calibra® Universal Self-Adhesive Resin Cement is recommended for self-adhesive cementation of full coverage crowns and bridges fabricated from CEREC® Zirconia+ (see full instructions for use). Universal / self-adhesive resin type cements indicated for zirconia ceramic cementation from other manufacturers can be used following their respective instructions for use.

#### 3.6.3 Adhesive cementation (Full Coverage Crowns and Bridges)

Dentsply Sirona Calibra® Ceram Adhesive Resin Cement is recommended for adhesive cementation of full coverage crowns and bridges fabricated from CEREC® Zirconia+ (see full instructions for use). Adhesive resin cement systems indicated for zirconia ceramic cementation from other manufacturers can be used following their respective instructions for use.

## 4. HYGIENE AND DISPOSAL

#### Cross-contamination

Do not reuse single use products. Dispose of in accordance with local regulations. Finished device should be disinfected per manufacturers' recommendation of disinfection material.

Following materials are considered compatible with CEREC® Zirconia+ restorations:

- 80% ethanol
- 70% 2-propanol

#### 4.1 Disposal

Dispose of in accordance with local regulations.

## 5. LOT NUMBER, EXPIRATION DATE AND CORRESPONDENCE

**5.1 Do not use after expiration date. ISO standard uses: "YYYY-MM-DD"**

**5.2 The following numbers should be quoted in all correspondences:**

- Reorder number
- Lot number
- Expiration date

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

A summary of the safety and clinical performance (SSCP) for this product can be found (upon activation) at <https://ec.europa.eu/tools/eudamed> by searching using the Basic UDI-DI number ++EDDISMCERAMICSZ2 and at <https://dentsplysirona.com/ifu> using the reference number (REF).

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DeguDent GmbH