CEREC MTL® Zirconia – Block for CEREC® Multi Transitional Layer Zirconia

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Block for CEREC®

INSTRUCTIONS FOR USE - ENGLISH

CAUTION: This is a medical device. For dental use only.
USA: Rx only.

1. PRODUCT DESCRIPTION

CEREC MTL® Zirconia blocks for CEREC® are material comprised of zirconia ceramics for the fabrication of individually designed restorations using a CAD/CAM procedure. The aesthetic features of CEREC MTL® Zirconia enable application as fully anatomical crowns and bridges. Indirect restorations are fabricated by milling CEREC MTL® Zirconia blocks using a Dentsply Sirona CAD/CAM system.

The CEREC MTL® Zirconia blocks are provided in a partially sintered state, then milled enlarged by the CEREC CAD/CAM system. Restorations are individually processed to specification, and finally, densely sintered in the CEREC SpeedFire Sintering Furnace. The software determines the sintering program according to the indication.

Note that CEREC® software 5.1.3 with material pack or higher is required.

CEREC MTL® Zirconia dental ceramic material is type II, class 5 pursuant to ISO standard 6872.

Patient target group: No restrictions.

Intended User: Dental professionals only: dentists and dental technicians (Rx only).

1.1 Indications

CEREC MTL® Zirconia are blanks used for fabricating dental restorations.

- CEREC MTL® Zirconia is indicated for:
 - fully anatomical anterior and posterior crowns
 - fully anatomical 3-unit anterior and posterior bridges
 - Onlays
 - Inlays
 - Veneers

1.2 Contraindications

CEREC MTL® Zirconia is contraindicated for:

- more than one bridge pontics
- more than one cantilever bridge unit
- patients with parafunctions in particular for bruxism
- insufficient oral hygiene
- insufficient preparation results
- insufficient hard tooth substance
- patients who have known allergies or sensitivities to the chemical ingredients of the material
- conventional or self-adhesive insertion of inlays, onlays, veneers

Intended Purpose

CEREC MTL® Zirconia products are ceramic materials for dental treatments.

1.3 Delivery Forms (Some delivery forms may not be available in all countries)

- CEREC MTL® Zirconia blocks are available in blended shades, to match the VITA classical A-D shade guide tab
- CEREC MTL® Zirconia blocks are available in 2 sizes: mono (20 mm) for most single-unit crowns and medi (39 mm) for most 3-unit bridges

1.4 Technical/physical data

Property	Unit	Value
3-point flexural strength (ISO 6872)	MPa	> 850
Coefficient of thermal expansion (20-500 °C)	10 ⁻⁶ · K ⁻¹	approx. 10.3
Chemical solubility (ISO 6872)	µg/cm ²	< 20
Density after sinter firing	g/cm ³	approx. 6.05

1.5 Chemical composition

Component	Wt%
ZrO ₂	88-93
Y ₂ O ₃	6-8
HfO ₂	1-3
Al ₂ O ₃	0-1
Pigments	0-1

Note:

- The technical/physical values given are typical measurement results and refer to in-house manufactured samples and measuring instruments in the company.
- If samples are prepared using different methods and measurement equipment, other measuring results may be obtained.

1.6 Compatible Stains and Glazes

The use of spray or paint on glaze or stains is optional with CEREC MTL® Zirconia CAD/CAM blocks. CEREC MTL® Zirconia CAD/CAM blocks are compatible with Dentsply Sirona Universal Stains and Glaze System, DS Universal Spray Glaze Fluo, DS Universal Spray Glaze and VITA AKZENT® Plus GLAZE LT Spray.

1.7 Compatible Cements

CEREC MTL® Zirconia CAD/CAM block restorations can be temporarily cemented and are compatible with traditional eugenol and non-eugenol temporary cements. For final cementation, CEREC MTL® Zirconia CAD/CAM block restorations are compatible with Universal/Self-Adhesive, Adhesive Resin Cement and Conventional Cement systems designed for Zirconia ceramic cementation (e.g. Dentsply Sirona cement systems Calibra® Bio, Calibra® Universal and Calibra® Ceram). See complete Instructions for Use of selected cement system.

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these instructions for use.

Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards.

Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Product reliability

Information on reporting serious incidents in connection with medical devices, general risks associated with dental treatments, residual risks and (if applicable) short clinical safety and performance reports (SSCPs) can be found at https://www.vita-zahnfabrik.com/product_safety
The safety data sheets can be downloaded at <https://www.vita-zahnfabrik.com> or requested by fax at (+49) 7761-562-233.

2.2 Precautions

- This product is intended to be used only as specifically outlined in these Instructions for Use. Any use of this product inconsistent with the Instructions for Use is at the discretion and is the sole responsibility of the practitioner.
- Wear suitable protective eyewear, clothing and gloves. Protective eyewear is recommended for patients.

- Contamination of the preparation or margin area with saliva, blood, water, or hemostatic agents during adhesive cementation may lead to an adhesive failure. Ensure adequate isolation and tissue management techniques during adhesive cementation.
- Devices marked “single use” on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- CEREC MTL® Zirconia CAD/CAM block restorations require adequate preparation reduction and restoration thickness. Insufficient wall thickness may lead to premature failure.
- CEREC MTL® Zirconia blocks are delivered in a pre-sintered state. In this state, the material can be processed very well, but does not yet have the properties it has after sintering. Careful handling in this state is necessary for this reason.
- CEREC MTL® Zirconia block restorations must be sintered and must be polished, or polished and glazed before insertion. Direct insertion without sintering and polishing or polishing and glaze may lead to failure.
- Final shading of CEREC MTL® Zirconia block restorations depends on correct sintering time. Therefore, only restorations of the same shade can be grouped into one sintering cycle. Do not group restorations of different shades/sintering times into one sintering cycle.
- Use only in well ventilated areas.
- Insufficient data exists to support use of CEREC MTL® Zirconia CAD/CAM blocks to fabricate resin-bonded-retainer bridges (“Maryland” bridges), endodontic post and cores, or implant abutments.
- Avoid adjusting sintered restorations with milling tools, especially in the connector area. Flexural strength may be compromised (See step-by-step Instructions).
- Please check the packaging and the material immediately upon receipt for an intact condition with regard to:
 - integrity of the packaging
 - integrity of the product (no disruptions, cracks or shade irregularities must be noticeable).
 - The manufacturer’s name, VITA Zahnfabrik and the CE marking must be present on the packaging.
- The success of any restoration depends on its fit onto the underlying tooth structure with minimal discrepancies.
- The ability to produce a routinely smooth, sound and well-fitting restoration requires strict adherence to certain fundamentals. Sharp edges of preparations or abutments may cause fracture of the respective zirconia structures. These sharp edges must be avoided (round off with wax before the scan).

2.3 Interactions

- CEREC MTL® Zirconia CAD/CAM blocks are designed to be fabricated using a Dentsply Sirona CAD/CAM system. Milling blocks using non-compatible CAD/CAM systems may lead to inadequate or unacceptable restorations.
- CEREC MTL® Zirconia CAD/CAM blocks are designed to be dry milled only. Do not mill in a wet-milling CAD/CAM system. Pre-sintered blocks and restorations absorb moisture. If restorations are milled wet or otherwise become wet in pre-sintered state, prolonged drying is required prior to sintering. Re-milling dry is recommended.

2.4 Adverse Reactions

- No adverse reactions have been reported for CEREC MTL® Zirconia CAD/CAM blocks. Should you hear or receive information about any adverse effects, notify Dentsply Sirona.

2.5 Storage Conditions

- Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product. Store in a dry place and protect from moisture.
- Store the CEREC MTL® Zirconia blocks in the original packaging and in a dry place.
 - Take care when handling these, that the CEREC MTL® Zirconia blocks are not exposed to any blows or vibrations.
 - Take care that the materials are not allowed to be touched with wet hands. Handle with dry gloves only.
 - The materials must not be contaminated with foreign substances (e.g., glass-ceramic or metal grinding dust).

- They can be disposed of with household waste. The products labelled with a pictogram for hazardous substances are to be disposed of as hazardous waste. Recyclable waste (such as attachments, paper and plastics) must be disposed of using appropriate recycling systems. If necessary, contaminated product residues should be pretreated in accordance with regional regulations and disposed of separately.

3. STEP-BY-STEP INSTRUCTIONS

Please adhere to the following workflow when working with this product:

3.1 Tooth shade determination

Tips on determining tooth shade:

- Whenever possible, determine the tooth shade using a standardized daylight lamp.
- Make sure that the surroundings are neutral in terms of color. Lipstick and cosmetics should be removed, and bright clothing covered by a grey cape.
- Make your selection quickly. In case of doubt, trust your initial decision as the eye already tires after approx. 5-7 seconds.
- For restorations that will be final polished only (no stain or glaze) the shade effect/shade intensity increases depending on the degree of polishing. It is recommended, to select a block shade brighter than the target shade.
- Disinfect the shade guide every time after performing shade determination according to manufacturer’s instructions.

Note:

- For digital shade determination, use VITA Easyshade® V. For visual shade determination, the VITA classical A1-D4® shade guide.

3.2 Preparation

3.2.1 Basic preparation guidelines for all-ceramic restorations

- No angles or edges.
- The ideal preparation is a shoulder preparation with rounded inner edges and/or a chamfer preparation.
- The indicated dimensions reflect the minimum thicknesses for CEREC MTL® Zirconia restorations.

3.2.2 General preparation guidelines

- Evenly reduce the anatomical shape while observing the given minimum wall thicknesses below.
- For conventional and/or self-adhesive cementation, the preparation must demonstrate retentive surfaces.
- Preparation angles: 4-8° for conventionally and self-adhesive cementation.


3.2.3 Design parameters

All information refers to sintered restorations	Minimum wall thickness in mm/ Connector cross-sections		
Anterior and posterior crowns, inlays, onlays and veneers	incisal	0.6	-
	occlusal	0.6	
	axial	0.5	
Fully anatomical 3-unit anterior bridges	incisal	0.6	9 mm ²
	axial	0.6	
Fully anatomical 3-unit posterior bridges	occlusal	0.7	12 mm ²
	axial	0.6	

Note:

- The minimum wall thicknesses refers to fully sintered restorations.

3.2.4 Design of the connectors

Note:

- The height of the connector surfaces is the highest possible that can be selected (Fig. 1).
- The height should be at least as high as the width (Figs. 1 and 2).
- Sharp corners and edges are to be avoided.
- Failure to observe the given minimum wall thicknesses and connector dimensions may result in clinical failures, such as cracks or fracture of the restoration.

3.2.5 Design criteria**Improper Design – To reduce restoration fracture risk**

The applicable minimum wall thicknesses must be observed.

- The aim is to obtain a uniform wall thickness.
- For implant-supported restorations, depending on the fabrication process, sharp edges may exist, which may cause fracture of the respective zirconia structures during the period of wearing. These sharp edges must be avoided in general (e.g. can be rounded off with wax before the scan).

3.3 CAD/CAM Process (scanning, designing, milling)

CEREC MTL® Zirconia restorations are produced with CEREC CAD/CAM systems by Dentsply Sirona. If you have any questions about these systems, please contact Dentsply Sirona.

3.3.1 Scanning the restoration with CEREC Primescan/Omniscam**3.3.2 Use CEREC SW 5.1.3 with material pack or higher.**

- Select CEREC MTL® Zirconia in the material dialog. The Dentsply Sirona CAD/CAM milling unit will prompt you to insert a CEREC MTL® Zirconia CAD/CAM Block. For detailed processing, please consult the Instructions for Use and technical manuals of the appropriate CAD/CAM systems. Make sure to follow the manufacturers' recommendations.
- Make sure that the chamber of the milling unit is clean and dry.

**Dry mill only – To reduce the risk of moisture contamination**

CEREC MTL® Zirconia CAD/CAM Blocks can be dry milled only.

Use of wet-milling capable milling unit is not recommended, due to possible moisture contamination requiring prolonged drying time.

**Material Shrinkage – To reduce the risk of improper fit**

Since the CEREC MTL® Zirconia shrinks by approximately 20-25% during sintering, the shrinkage factor of the respective batch, which is included in the Data Matrix code or on the code on the block (top side), must be read into the software or manually entered.

3.3.3 Manufacturing phase

- Selecting the color

When using zirconia in a CEREC SpeedFire sintering furnace, the color of the block must be selected in advance in the software. This ensures the use of the appropriate parameters.

In the ADMINISTRATION phase, you have selected CEREC MTL® Zirconia.

 1. Select the "Select Color" step.
 2. Select the desired color by choosing the color in the color center.
 3. Press the "OK" button.
- Positioning restorations in the block

The normal positioning of the restoration in the multi transitional layered block is centered.

This captures the entire color gradient. For restorations to have a clearly visible enamel area, they must be positioned in the CAM software as high as possible in the top of the block. It is recommended to place the restoration approximately 1 mm below the top edge of the block to obtain a clearly visible incisal portion.

For a smaller proportion of enamel, we recommend placing the restoration – especially with molars – as far down as possible in the cervical area.

- Starting the production process

Once you have completed the design and assessed the restoration in the preview, you can produce the restoration.
- Enlargement factor

VITA Zahnfabrik determines the enlargement factor in all three spatial dimensions (X-, Y-, Z-direction) and integrates the information as plain text or code in the print on the blank. The barcode is scanned in the device. If this is not possible, the code can also be entered as plain text (e.g., *Z2300EB).

3.4 Sprue removal, adjusting, pre-polishing

- After the milling process and prior to sintering, a diamond or a tungsten carbide milling tool has to be used to reduce marginal edges or carefully re-contouring the fissures as desired. Remove the sprue prior to sintering.
- Bridge restorations may not be separated approximately using a diamond separating disc, since this may result in breaking points on the connectors.
- Do not thin walls beyond the minimum wall thicknesses.

Tips:

- In order to facilitate high-gloss polishing of fully anatomical restorations, it is recommended to smooth the milled restoration before hard sintering with a fine diamond smoothing instrument and/or to pre-polish the restoration with silicone-free polishers.
- For pre-polishing, use a polyurethane-bonded rubber polisher. Silicone polishers may leave residue.

Note:

- No drying is required for dry-milled restorations.
- In the case of high humidity:
 - In an environment with a high level of humidity, the restorations can absorb moisture after drying. For this reason, sintering must take place within a maximum of one hour after milling.

**Airborne particles – To reduce inhalation risk**

- Do not inhale abrasive dusts.
 - Use a vacuum system and wear a mask.
- Adjusting CEREC MTL® Zirconia restorations should always be carried out in the unsintered state.
 - Only use suitable milling tools (e.g., fine-grain diamonds, fine-toothed carbide burs, zirconia stones), low speed ($\leq 20,000$ 1/min), and little pressure. Avoid overheating the ceramic.
 - Make sure that the minimum wall and connector thicknesses are maintained when adjusting/pre-polishing.
 - Remove any traces of milling or polishing residue before sintering to avoid inaccuracy of fit caused by sintered on milling dust. Remove with a brush or oil-free compressed air.
 - Note: The restorations must not be sandblasted or cleaned with a steam jet before sintering!
 - Make sure you do not inhale abrasive dusts. Use a vacuum system and wear a mask.

3.5 Sintering**3.5.1 Preparation for sintering**

After the milling process (and prior to sintering) a tungsten carbide tool shall be used to separate the restoration from the blank. Remove the sprue prior sintering. If further adjustments are needed, it is recommended to do the adjustments before sintering.

**Airborne particles – To reduce inhalation risk**

- Do not inhale abrasive dusts.
- Use a vacuum system and wear a mask.

3.5.2 Sintering the restoration in the CEREC SpeedFire

A CEREC SpeedFire furnace from Dentsply Sirona is recommended for the sintering process.

Final shading of CEREC MTL® Zirconia block restorations depends on correct sintering time

Do not group restorations of different shades/sintering times into one sintering cycle.

Before sintering the restoration, it is recommended to free the restorations from dust using compressed air or a ceramic brush.

Do not inhale abrasive dusts. Use a vacuum system and wear a mask.

NOTE: Restorations made from CEREC MTL® Zirconia must be sintered in dry conditions.

When sintering CEREC MTL® Zirconia restorations in the CEREC SpeedFire, the CEREC software automatically transfers the job to the CEREC SpeedFire if the milling machine and the CEREC SpeedFire are connected.

NOTE: Maximum restoration size.

Observe the maximum furnace chamber size when loading the Furnace:

- Diameter: 38 mm
- Height: 20 mm

The restoration (including the Glazing Support Single/Multi Unit) must not exceed the size of the chamber, length (38 mm) and height (20 mm), otherwise the chamber may be damaged.

For further information please refer to the Operating Instructions of the CEREC SpeedFire.

-  Place the restoration with the occlusal surface facing down directly on the top door insulation.
- Start the process by touching the start icon. The furnace shuts automatically once the process starts.
- The furnace opens automatically following successful heat treatment. The process is not yet complete, as a cooling-off phase occurs in an open condition. There will be a signal tone once the cooling process has completed. When the unit's LED status display is green, the furnace can be unloaded.

Risk of injury

The restoration and parts of the door insulation may still be hot when the unit's LED status display is green. Always use tweezers to remove the restoration from the furnace. Allow the restoration to cool down for another five minutes before picking it up with hands.

- Only unload the furnace using metallic or ceramic tweezers. Plastic tweezers are not suitable for removing restorations, as the restoration is still very hot in this phase. Place the restoration on the fireproof tray for further cooling.

3.6 Adjusting after the sintering process

Surface Condition - To reduce risk of compromised flexural strength

The surface condition of ceramic materials is critical for their flexural strength. Adjustment of the sintered restorations with milling tools, especially in the connector area, must be avoided.

However, if adjustment is necessary, then follow these basic rules:

- Adjustment in the sintered state should be performed with fine diamonds in a high-speed handpiece with water cooling and with low pressure. Adjustments made with diamonds instruments must be followed by polishing.

- As an alternative it is possible to do adjustments with soft diamond rubber polishers and a handpiece at low speed and low pressure. The tool must be applied flat to minimize the chatter.
- Areas that are under tension in clinical use, i.e. primarily the connectors in bridge structures, should not be adjusted.

Tip:

- Work exclusively with polyurethane polishers. Residues of these polishers can be easily removed and burned out without leaving any residue.
- When using silicone-bonded polishers, there is a danger that the abraded silicone cannot be removed without residue.
- This may negatively affect the color/translucency and the bonding area towards the glazing material.

3.7 Polishing, Try-in

CEREC MTL® Zirconia restorations can either be polished, or polished and glazed.

3.7.1 Polishing technique

- The polishing of the occlusal surface, especially of the areas that are in direct contact with the antagonist, is particularly important in monolithic restorations.
- After functional adjusting, re-polish the surfaces of the adjusted occlusal surface very carefully.
- In general, if the surface is polished to a high gloss, it is significantly less abrasive according to laboratory tests. As a result, the high-gloss polish protects the antagonist from unwanted abrasion.

Note:

- With the "polishing technique", the shade effect may differ from that of the shade guide. Depending on the degree of polishing, the shade effect/intensity is increased. If necessary, it is recommended to select a block shade that is one shade brighter than the target shade.
- Polishing all-zirconium oxide restorations results only partially in a mother-of-pearl-type shiny surface, which thus appears to be different.
- The mother-of-pearl effect is the iridescent effect of a polished surface in incident light. The same effect that occurs in sea shells. A similar phenomenon can be observed in well-polished zirconium oxide. Zirconium oxide demonstrates a high refractive index and, when polished well, a high reflection factor. As the shade effect/shade intensity increases depending on the degree of polishing, it is recommended, to select a block shade brighter than the target shade.
- CEREC MTL® Zirconia can be polished with standard polishing agents for zirconia ceramics.

Polishing with VITA SUPRINITY® Polishing Set

- The pre-polishing of ground areas is done with the pink rubber polishers containing diamonds of the VITA SUPRINITY® Polishing Set technical/clinical at a speed of 7,000 to 12,000 rpm.
- High-gloss polishing is then carried out with the diamond-coated, gray rubber polishers at a reduced speed of 4,000 to 8,000 rpm.

For further information please refer Instructions for Use of the VITA SUPRINITY® Polishing Set*.

Polishing with MEISINGER LUSTER® Kits for Zirconia

1. Trimming the contact points if needed (9736H), Recommended rotary Speed: 7,000-12,000 rpm.
2. Smoothing of the occlusal surface (DCA04), Recommended rotary Speed: 7,000-12,000 rpm.
3. Smoothing of the external shape (DCA06), Recommended rotary Speed: 7,000-12,000 rpm.
4. High-shine polishing of the occlusal surfaces without high pressure (DCA10), Recommended rotary Speed: 7,000-12,000 rpm.
5. High-shine polishing of the external shapes (DCA12), Recommended rotary Speed: 7,000-12,000 rpm.

For further information please refer Instructions for Use of the MEISINGER LUSTER® Kits*.

* Not a registered trademark of Dentsply Sirona.

3.7.2 Optional Try-In

- Try-in the restoration for marginal and proximal fit. Make any necessary adjustments as outlined above.

Contamination – To reduce risk of infection

Restorations should be polished, cleaned and disinfected before and after optional try-in. See Hygiene section below.

- Adjustments made during try-in must be re-polished as outlined above, prior to optional stain and glaze application and final delivery.
- The polishing of the occlusal surface, especially of the areas that are in direct contact with the antagonist, is particularly important in monolithic restorations.
- After functional grinding, re-polish the surfaces of the ground occlusal surface very carefully.
- In general, if the surface is polished to a high gloss, it is significantly less or even not abrasive, according to laboratory tests. As a result, the high-gloss polish protects the antagonist from unwanted abrasion.

3.8 Optional Staining and Glazing

If staining and glazing is not required/desired, skip to section 3.9 Pre-cementation surface preparation.

Surface Condition – To reduce risk of excessive wear

- Prior to glazing, ensure restorations are properly sintered and polished.
- Restorations must be clean and dry before application.
- CEREC MTL® Zirconia restorations can be glazed with: Dentsply Sirona Universal Stain and Glaze System.

First, follow the polishing instructions from Section Polishing technique above 3.7.1.

After polishing, clean the restoration using either an ultrasonic cleaner or steam cleaner, prior to Spray Glaze or paint-on application. Ensure the restoration is free of contamination and completely dried with oil-free air prior to paint-on or Spray Glaze application.

3.8.1 Preparation of the restoration for Staining and/or Glazing

Use the Glazing Support Single/Multi Unit holders, and CEREC SpeedPaste as shown in the images below. Fill the restoration with CEREC SpeedPaste, applying the paste evenly up to the restoration margin. Place the Glazing Support Single/Multi Unit holders in the paste and pick up the tweezers, if required. Observe the operating instructions supplied with the CEREC SpeedPaste.

3.8.2 Glazing CEREC MTL® Zirconia with DS Universal Stain and Glaze or DS Universal Spray Glaze Fluo or DS Universal Spray Glaze

Glazing CEREC MTL® Zirconia with DS Universal Stain and Glaze

1. Withdraw desired amount of DS Universal Glaze from the jar and place it on the mixing palette.
NOTE: If the stain or glaze inside the jar has separated, mix thoroughly with a glass, plastic or zirconia spatula.
2. If a thinner consistency is desired, dilute the material with the Dentsply Sirona Stain and Glaze Liquid (REF 601315/15 ml, 601350/50 ml).
3. Apply a thin layer of the glaze material all over the crown surface in the usual manner using a brush. Make sure not to apply the glaze material too thick or too thin. Too thin of a layer will result in a more matte finish. Too thick of a layer can result in puddling and/or pitting of the material.
4. If a more intensive shade effect is desired, DS Universal Stains can be applied on the glazed surface of the restoration. Apply a thin layer of the stain material on the crown surface in the usual manner using a brush. Make sure not to apply the glaze material too thick or too thin.

Follow the Instructions for Use (IFU) for Universal Stain & Glaze.

Glazing CEREC MTL® Zirconia with DS Universal Spray Glaze Fluo or DS Universal Spray Glaze

1. Prior to use, ensure the spray nozzle is not clogged and ensure the nozzle is firmly assembled to the spray head.
2. Protect the intaglio of the restoration and any restoration surface with the CEREC Speed Paste where glaze application is not desired.
3. Vigorously shake the spray can immediately prior to each use. Ensure the mixing balls inside the spray can are free and make an audible noise when shaking the spray can.
4. Keep the spray can in as vertical position as possible while applying spray. And maintain a distance of 6-10 cm (2.5-4.0 in) between the nozzle outlet and the restoration surface.
5. Apply a uniform layer glaze to the surface of the restoration.
6. The applied glaze should be thin and uniform.

Follow the Instructions for Use (IFU) for DS Universal Spray Glaze.

3.8.3 Glazing the CEREC MTL® Zirconia Restoration in the furnace

Glazing CEREC MTL® Zirconia in the CEREC SpeedFire

Position the Support with the restoration centrally on the top door insulation and ensure that the Glazing Support Single/Multi Unit or restoration do not protrude out from the door insulation; otherwise, these may collide with the furnace chamber. Select the “GLAZE” program on the CEREC SpeedFire control panel and allow the unit to cycle.

Glazing CEREC MTL® Zirconia with other furnaces

- For firing the restoration in furnaces from other manufacturers, follow the firing parameters outlined in Table below.

Drying [min]	Closing [min]	Pre-heating temperature [°C]	Pre-heating [min]	Heating rate [°C/min]	Final temperature [°C]	Vacuum [min]	Holding time [min]	Cooling [min]
0	2	400	0	55	760	0	2:00	0

- Position the Support with the restoration centrally on the top door insulation and ensure that the Glazing Support Single/Multi Unit or restoration do not protrude out from the door insulation; otherwise, these may collide with the furnace chamber.
- Observe maximum restoration size.
 - Observe the maximum Furnace chamber size when loading the furnace. Diameter: 38 mm; Height: 20 mm
 - The restoration (including with the Glazing Support Single/Multi Unit) must not exceed this length (38 mm) and height (20 mm), otherwise the chamber may be damaged.
- Only load the furnace with approved materials, which must be assigned uniquely to the heat treatment job in the furnace.
- Use the restoration holders, firing paste (CEREC SpeedPaste) and tweezers supplied in the package for loading the furnace for the glazing process.
- Also observe the operating instructions supplied with the CEREC SpeedPaste.
- Apply as small a portion of the paste as possible to get the firing pin to hold. The restoration must not be completely filled, as undesired expansion may occur, resulting in cracks or fractures.
- Place the tray in the paste and pick up the tweezers, if required.
- After fixing the restoration to the restoration holder, position the holder centrally on the top door insulation and ensure that the holder or restoration do not protrude out from the door insulation; otherwise, these may collide with the furnace chamber.

Warning: Risk of burning

- The restoration and parts of the door insulation may still be hot.
- Always pick up some tweezers to remove the restoration from the furnace.
- Allow the restoration to cool down for another five minutes before picking it up with your hands.

- As the restoration is still very hot in this phase, only unload the furnace using metallic or ceramic tweezers. Plastic tweezers are not suitable for removing restorations. Place the restoration on the fireproof tray (A) for further cooling.

3.9 Pre-Cementation Surface Preparation

- Sand-blast the internal surface of the restoration using 50 µm aluminum oxide at a max pressure of 2.5 bar. (35 PSI)
- Clean and disinfect the restoration as outlined in Hygiene section below.

Contamination – To reduce risk of infection

Restorations should be polished, stained and glazed (optional, if desired) and sandblasted internally before cleaning and disinfection.

Restorations should be cleaned and disinfected immediately prior to delivery.

Follow instructions in Hygiene section below before delivery.

3.10 Cementation

3.10.1 Temporary cementation

- Monolithic CEREC MTL® Zirconia restorations can be cemented temporarily since they exhibit high inherent strength and there is no risk of damaging the veneer when removing the restoration prior to permanent cementation.
- Use hemostats, Baade pliers or other suitable crown removal instrument for removal along the path of insertion. Care should be taken and tension, twisting or rocking must be avoided when removing the restoration
- CEREC MTL® Zirconia restorations are compatible with all temporary cement materials can be used for temporary cementation Follow cement manufacturer's complete Instructions for Use. Non-eugenol cements are recommended if final cementation will be with resin cement.
- Re-sandblasting of internal surfaces followed by cleaning and disinfection (see Hygiene section below) are required prior to definitive cementation.

3.10.2 Final Cementation

Conventional Cements (Full Coverage Crowns and Bridges)

Full coverage crowns and bridges fabricated from CEREC MTL® Zirconia can be conventionally cemented provided that the preparation guidelines are followed. Cement the CEREC MTL® Zirconia restoration with Dentsply Sirona Calibra® Bio Bioceramic Luting Cement following the product Instructions for Use (IFU). Resin modified glass ionomer (RMGI) or glass ionomer (GI) type cements indicated for zirconia ceramic cementation from other manufacturers can be used following their respective instructions for use.

Universal/Self-Adhesive Cements (Full Coverage Crowns and Bridges)

Full coverage crowns and bridges fabricated from CEREC MTL® Zirconia can be cemented using Universal/Self-Adhesive or Adhesive Resin cements provided the preparation guidelines are followed. Cement the CEREC MTL® Zirconia restoration with Dentsply Sirona Calibra® Universal Self-Adhesive Resin Cement or Calibra® Ceram Adhesive Resin Cement following the product Instructions for Use (IFU). Universal or adhesive resin type cements indicated for zirconia ceramic cementation from other manufacturers can be used following their respective instructions for use.

Adhesive Resin Based Cements (Onlays, Inlays, Veneers)

Onlays, inlays, and veneers fabricated from CEREC MTL® Zirconia must be cemented using Adhesive Resin cements provided the preparation guidelines are followed. Cement the CEREC MTL® Zirconia restoration with Dentsply Sirona Calibra® Ceram Adhesive Resin Cement following the product Instructions for Use (IFU). Adhesive resin type cements indicated for zirconia ceramic cementation from other manufacturers can be used following their respective instructions for use.

4. HYGIENE AND DISPOSAL

Cross-contamination

Do not reuse single use products. Dispose of in accordance with local regulations. Finished device should be disinfected per manufacturers' recommendation of disinfection material.

Following materials are considered compatible with CEREC MTL® Zirconia restorations:

- 80% ethanol
- 70% 2-propanol

4.1 Disposal

- Block remains and the block holder do not need to be disposed of separately. They can be disposed of as normal household waste in accordance with local regulations.

Symbol explanations

Medical device	
CE Mark MP Class IIa	
The product may only be sold by a dentist or when prescribed by a dentist (valid only for USA)	
Refer to instructions for use	
Store in a dry location	
Handle with care. Do not throw.	
Refer to information	
Single use	
Manufacturer	
Date of manufacture (YYYY-MM-DD)	
Product number	
Batch description	
Matrix Code	
Recycling	
Distributor	

5. LOT NUMBER, PRODUCTION DATE AND CORRESPONDENCE

- 5.1 For Production date, ISO standard "YYYY-MM-DD" is used.
- 5.2 The following numbers should be quoted in all correspondences:
 - Reorder number (REF)
 - Lot number
 - Production date
- 5.3 Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

Please note: Our products must be used in accordance with the instructions for use. We do not accept any liability for damage resulting from incorrect handling or usage. The user is furthermore obliged to check the product before use with regard to its suitability for the intended area of applications. We cannot accept any liability if the product is used in conjunction with materials or equipment from other manufacturers that are not compatible or not authorized for use with our product and this results in damage.

Date of issue of this information: 2023-05

After the publication of this information for use, any previous versions become obsolete. The current version can be found under <https://www.dentsplysirona.com/ifu>.

VITA Zahnfabrik has been certified and the following product bears the CE mark:

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Distributed by
Degudent GmbH
Rodenbacher Chaussee 4
63457 Hanau
Germany
Tel. +49 6181 59-50

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VITA Zahnfabrik H. Rauter GmbH & Co. KG
Zweigniederlassung Basel c/o Perrig AG
H. Rauter GmbH & Co.KG