IRM® Caps™

Intermediate Restorative Material

Intermediate Restorative Material INSTRUCTIONS FOR USE	ENGLISH	2
Matériau de restauration intermédiaire MODE D'EMPLOI	FRANÇAIS	4
Materiale per restauro provvisorio ISTRUZIONI PER L'USO	ITALIANO	6
Temporäres Füllungsmaterial GEBRAUCHSANWEISUNG	DEUTSCH	8
Material de restauración intermedio INSTRUCCIONES DE USO	ESPAÑOL	10
Material de restauração temporária INSTRUÇÕES DE UTILIZAÇÃO	PORTUGUÊS	12
Provisorisch restauratiemateriaal GEBRUIKSAANWIJZING	NEDERLANDS	14
Provisoriskt restaureringsmaterial BRUKSANVISNING	SVENSKA	16
Provisorisk fyldningsmateriale BRUGSANVISNING	DANSK	18
Midlertidig restaureringsmateriale BRUKSANVISNING	NORSK	20
Väliaikainen paikkamateriaali KÄYTTÖOPAS	SUOMI	22
Tymczasowy materiał wypełniający INSTRUKCJA UŻYCIA	POLSKI	24
Tarpinio restauravimo medžiaga NAUDOJIMO INSTRUKCIJOS	LIETUVIŲ K.	26
Pagaidu plombējamais materiāls LIETOŠANAS INSTRUKCIJA	LATVIEŠU	28
Ajutine restauratsioonimaterjal KASUTUSJUHEND	EESTI	30
Dočasný výplňový materiál NÁVOD K POUŽITÍ	ČESKÝ	32
Material pentru restaurări intermediare INSTRUCȚIUNI DE UTILIZARE	ROMÂNĂ	34
Υλικό μεταβατικών αποκαταστάσεων ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	ΕΛΛΗΝΙΚΑ	36

Manufactured by
Dentsply Caulk
38 West Clarke Avenue
Milford, DE 19963, USA
Tel: 1-302-422-4511
www.dentsplysirona.com

Dentsply DeTrey GmbH
78467 Konstanz
Germany
Tel.: 49-7531-583-0
www.dentsplysirona.com

IRM® Caps™

Intermediate Restorative Material

INSTRUCTIONS FOR USE - ENGLISH

Caution: This is a medical device. For dental use only by dental professionals.
USA: Rx only.

1. PRODUCT DESCRIPTION

IRM® Intermediate Restorative Material is a polymer-reinforced zinc oxide-eugenol composition restorative material designed for intermediate restorations intended to remain in place for no longer than one year.

IRM® Intermediate Restorative Material may also be used as a base under restorative materials and cements that do not contain resin components.

IRM® Material complies with ISO 3107 zinc oxide/eugenol type II: for bases and temporary restorations.

1.1 Indications for Use

IRM® Material is indicated for:

- Intermediate (temporary, less than 1 year) direct restorative filling material
- Base material under non-resin-containing permanent restorations

1.2 Contraindications

- IRM® Material is contraindicated for use with patients who have a known hypersensitivity or severe allergic reaction to eugenol.
- IRM® Material is contraindicated for use with patients who have a known hypersensitivity to polymethylmethacrylate.
- IRM® Material is contraindicated for direct application to dental pulp tissue (direct pulp capping).
- IRM® Material is contraindicated as a base under resin-containing adhesives, restoratives or cements.

1.3 Delivery Forms (Some delivery forms may not be available in all countries)

- Pre-dosed, encapsulated material for machine mixing, minimum deliverable volume 0.15 ml
- Alternatively available as:
- A powder/liquid version for hand mixing

1.4 Composition

Powder: Zinc oxide; Poly-Methyl Methacrylate (PMMA) powder; Zinc Acetate; Pigment
Liquid: Eugenol; Acetic acid

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

Safety alert symbol.

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

1. Liquid contains eugenol. Powder contains polymethacrylate. Components and mixed material may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible individuals.

- **Avoid eye contact** to prevent irritation and possible corneal damage. In case of eye contact, rinse immediately with plenty of water and seek medical attention.
- **Avoid skin contact** to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and alcohol and wash thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
- **Avoid contact with oral soft tissues/mucosa** to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.

2. Improved or reinforced zinc oxide-eugenol products such as IRM® Material are not suitable for use as root canal filling materials as inconsistent, incomplete or delayed resorption may occur.

3. IRM® Material may irritate pulp or soft tissue. Cover dentin close to the pulp (less than 1 mm) with a hard setting calcium hydroxide liner (Dycal® Liner) following manufacturer's instructions.

4. IRM® Material should not be used with patients who have a known hypersensitivity to any of the components.

2.2 Precautions

- This product is intended to be used as specifically outlined in the Instructions for Use. Any use of this product inconsistent with the Instructions for Use is at the discretion and sole responsibility of the practitioner.
- Insufficient data exist to support the use of IRM® material in a root canal as a cement, dressing, medicament or repair material.
- Wear suitable protective eyewear, mask, clothing and gloves. Protective eyewear is recommended for patients.
- Contact with saliva, blood and sulcus fluid during application may cause failure of the restoration. Use adequate isolation such as rubber dam.
- Insufficient data exist to support the use of IRM® Material as a final or temporary cement.
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- Increased temperature and/or humidity reduce available work time and accelerate set time. Material will set faster in the mouth than at the bench. Do not use bench cure rates to judge mouth cure rates.
- Do not open capsule prior to trituration.
- Do not use instruments which have not been properly (re)processed. Before re-use, instruments must be properly cleaned/decontaminated to reduce the risk of cross contamination.
- Slow speed amalgamators or slow speeds on amalgamators are not recommended.
- Capsules are provided ready for use and are designed for single use only. Do not attempt to reprocess for reuse.

- Use only new, clean capsules, stored in covered storage, such as drawers and cabinets, away from potential contamination, and removed from their original packaging using either clean gloved hands or a clean suitable instrument at time of use. To reduce risk of infection, do not use capsules if packaging is damaged or compromised.
- Product is not self-adhesive and may not be used in conjunction with conventional bonding techniques. Cavity preparation must provide for the mechanical retention of the product.
- Numerous factors can affect restoration longevity. Material is intended as an interim material. Check restoration marginal integrity and retention at every periodic examination.

Interactions

Do not use this product in conjunction with (meth)acrylate based dental materials since it may interfere with the hardening of the (meth)acrylate based dental materials.

2.3 Adverse Reactions

1. Product may irritate the eyes and skin.

- **Eye contact:** Irritation and possible corneal damage.
- **Skin contact:** Irritation or possible allergic response. Reddish rashes may be seen on the skin.
- **Contact with mucous membranes:** inflammation. (See Warnings)

2. Product may cause pulpal effects. (See Contraindications).

2.4 Storage Conditions

Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.

- Store in a well ventilated place at temperatures between 10 °C-24 °C / 50 °F-75 °F.
- Allow material to reach room temperature prior to use.
- Do not store with or near materials containing (meth)acrylate resin components.
- Protect from moisture.
- Do not freeze.
- Do not use after expiration date.
- Do not use if the lot number and/or expiration date are missing or illegible.

3. STEP-BY-STEP INSTRUCTIONS

3.1 Cavity preparation

3.1.1 Isolate the field of operation in preferred manner. Rubber dam is recommended.

3.1.2 Cavity preparation must provide for the mechanical retention of the product.

3.1.3 If caries is present, completely excavate using low-speed and/or hand instrumentation.

3.1.4 Place matrix and wedge where indicated.

3.1.5 Gently dry preparation with cotton pellet. Avoid desiccation.

3.2 Pulp Protection

In deep cavities cover the dentin close to the pulp (less than 1 mm) with a hard setting calcium hydroxide liner (Dycal® Liner, see complete Instructions for Use).

To reduce risk of pulp irritation.

Cover dentin close to the pulp (less than 1 mm) with a hard setting calcium hydroxide liner (Dycal® Liner) following manufacturer's instructions.

3.3 Mixing IRM® Caps™ Material

3.3.1 To activate the IRM® Caps™ capsule, hold vertically, grasp bottom half and FIRMLY TIGHTEN the top (see Fig. 1). As top is tightened, you will feel a "snap" as the liquid is released. Continue tightening until you are sure it cannot be seated further.

3.3.2 This complete tightening forces the liquid into the mixing chamber, thereby assuring proper handling and physical characteristics. Without delay, insert capsule (bottom first) into the amalgamator arms (see Fig. 2). Figure 4 gives approximate mixing times for various amalgamators. These are suggested ranges; however, your clinical experience in obtaining a good mix consistency with your unit remains the best guide to mixing time. Use the first five caps to establish the best mix time. (See 3.4 Guidelines for Optimal Mix)


Fig. 1

Mixing Discrepancies – To reduce risk of compromised performance.

- Do not activate capsule until ready to triturate.
- Use of a listed amalgamator is recommended.
- Trial mix to calibrate individual amalgamator performance is recommended
- Slow speed amalgamators or slow speeds on amalgamators are not recommended.


Fig. 2

3.3.3 Remove press cap (see Fig. 3) to gain access to the mixed IRM® Material.

Fig. 4: Approximate Mixing Time (seconds)

AMALGAMATORS	SETTING	SECONDS
ProMix™	High Speed (Rabbit)	12
Vari-Mix®-III	High Speed	12
Silamat® Model C		8
Wig-L-Bug®	High Speed	10
Wig-L-Bug® Model 80	High Speed	10
Pelton Crane*		30
Capmaster*	Not Recommended	


Fig. 3

*Not registered products of Dentsply Sirona

3.4 Guidelines for Optimal Mix

- 3.4.1 Clinical experience with your amalgamator is the best guide to mixing the IRM® Caps™ capsule. Variation in amalgamators (even from the same manufacturer) will occur and the trituration speed of amalgamators will vary with the line voltage on which they are operated. High-speed amalgamators are recommended for best results.
- 3.4.2 If a dry, stiff or rubbery mix is obtained, reduce trituration time, keep speed the same.
- 3.4.3 If a sticky and soft mix is obtained, increase trituration time, keep speed the same. Sticky mixes should be allowed to set a bit before application, or rolled in IRM® powder to eliminate stickiness.
- 3.4.4 If trituration results in a non-coalesced mass (i.e., small balls), increase trituration time, keep speed the same.
- 3.4.5 If the mixed mass of IRM® material does not appear in the press cap (see Fig. 3) then either replace the press cap and mull for 5 seconds or insert a small ended plugger into the cap and remove the mixed mass.

3.5 Placement of Mixed IRM® Material

Mixed material may be carried to the preparation with suitable placement instrument. Adapt, contour and shape with appropriate instruments.

3.5.1 Use as Intermediate Restorative Material

- Trim excess material. Material may be smoothed by lightly burnishing. Initial set occurs approximately 5 minutes from start of mix. Increased temperature, humidity, and powder/liquid ratio may accelerate set.
- Upon set, remove matrix band if present. Check and adjust occlusion as needed. Additional carving or occlusal adjustment may be accomplished with a round bur. Avoid gouging material with sharp burs at high speed. Additional polishing is not recommended.
- Periodically evaluate the intermediate restoration until final restoration is indicated. The intermediate restoration should be removed and final restoration completed within 1 year of placement.

3.5.2 Use as Protective Base

This product is a protective base for restorative material not containing polymeric resin components (e.g. final restoration of amalgam, gold).

Possible Interactions – To reduce risk of restoration failure.

The product should not be used in conjunction with resin-based adhesives, varnishes or restoratives since eugenol may interfere with hardening and cause softening of (meth)acrylate based dental materials.

- Apply mixed material to desired dentin surfaces and depth.
- Following completion of set (approximately 5 minutes from start of mix) cavity refinement may be completed.
- Complete restoration per restorative material manufacturer's Instructions for Use.

4. HYGIENE, PROCESSING AND DISPOSAL

Cross-contamination – To reduce risk of Infection.

- Do not reuse single use products. Dispose of in accordance with local regulations.
- Capsules cannot be reprocessed. Dispose of contaminated capsules in accordance with local regulations.
- When not in use, all supplies and instruments should be in covered storage, such as drawers and cabinets, and away from potential contamination.
- Dispensing materials with clean/disinfected gloves or hands or suitable transfer forceps in a separate room, bringing only what will be used into the operatory is strongly recommended.
- During treatment clinicians with patient contact should only handle items to be used.

4.1 Disposal

- Based on the information in the SDS, the waste generated must not be disposed of together with household garbage. Dispose of material as solid waste in a closed container.
- Dispose of contents/container to authorized hazardous or special waste collection point in accordance with any local regulation.
- Dispose of in accordance with national and local regulations.
- Contaminated containers and accessories that have been in contact with the patient must be disposed of as hazardous waste at risk of biological contamination.
- Dispose of uncontaminated empty used capsules in normal waste stream in accordance with local regulations.
- Always wear gloves when handling packaging and accessories.
- See SDS for further details.

5. LOT NUMBER, EXPIRATION DATE AND CORRESPONDENCE

1. Do not use after expiration date. ISO standard is used: "YYYY-MM-DD."
2. The following numbers should be quoted in all correspondences:
 - Reorder number
 - Lot number
 - Expiration date
3. Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.
4. Device Identification (Basic UDI-DI): ++D002DENTALCEMENTSC6

©2024 Dentsply Sirona. All Rights Reserved. **510201WEB** (R 10/18/24)

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