 **Intensiv SA**
Via al Molino 107
6926 Montagnola
Switzerland
www.intensiv.ch

P.277/09/2014

Instructions for use and safety recommendations

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Packung zu 3, 6

Intensiv Ortho-Strips Opener, Double-Sided, Ref. OS08OP
Intensiv Ortho-Strips Opener, One-Sided, Ref. OS08OP-R, rechts
Intensiv Ortho-Strips Opener, One-Sided, Ref. OS08OP-L, links

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Packung zu 1, 3, 6, 12
Intensiv Ortho-Strips, Ref. OS90
Intensiv Ortho-Strips, Ref. OS60
Intensiv Ortho-Strips, Ref. OS40
Intensiv Ortho-Strips, Ref. OS25
Intensiv Ortho-Strips, Ref. OS15
Intensiv Ortho-Strips Set, Ref. 060A (OS90, OS40, OS25, OS15)
Intensiv Ortho-Strips Set, Ref. 060B (OS60, OS40, OS25, OS15)

Packung zu 3, 6

Intensiv Ortho-Strips Coarse, Double-Sided, Ref. OS80C
Intensiv Ortho-Strips Coarse, One-Sided, Ref. OS80C-R, rechts
Intensiv Ortho-Strips Coarse, One-Sided, Ref. OS80C-L, links

Packung zu 1, 3, 6, 12

Intensiv Ortho-Strips One-Sided, Ref. OS15L, links
Intensiv Ortho-Strips One-Sided, Ref. OS15R, rechts
Intensiv Ortho-Strips One-Sided, Ref. OS25L, links
Intensiv Ortho-Strips One-Sided, Ref. OS25R, rechts
Intensiv Ortho-Strips One-Sided, Ref. OS40L, links
Intensiv Ortho-Strips One-Sided, Ref. OS40R, rechts
Intensiv Ortho-Strips One-Sided, Ref. OS60L, links
Intensiv Ortho-Strips One-Sided, Ref. OS60R, rechts
Intensiv Ortho-Strips One-Sided Set, Ref. 064 (OS15L, OS15R, OS25L, OS25R, OS40L, OS40R)

GB

Intensiv Ortho-Strips Opener

Oscillating instruments for easy and rapid opening of tight interdental contact points.

To be used in combination with the contra-angles:

- Intensiv Swingle, WG-69 LT Ortho, with light
- Intensiv Swingle, WG-69 A Ortho, without light

Alternatively:

- KaVo head 61 LG
- KaVo contra-angle E61C

Product description

- Serrated oscillating diamond-coated metal strip
- Diamond grit-size: 8 µm
- Available in 3 versions:
 - One-Sided Left (L): diamond coating on the left side of the strip
 - One-Sided Right (R): diamond coating on the right side of the strip
 - Double-Sided: diamond coating on both sides of the strip
- Strip height: 3.7 mm
- Length of the working part: 13 mm
- Strip thickness: 0.1±0.03 mm
- Sterilizable and reusable

Indication

Opening of tight interdental contact areas at the initial stage of the orthodontic IPR procedure.

Recommendations for use

- Apply Intensiv Ortho-Strips Opener in interproximal areas at the beginning of the IPR orthodontics protocol for the opening of tight contact points.
- Subsequently proceed with the enamel reduction protocol using Intensiv Ortho-Strips System.
- In case of very tight contact points separate the teeth by placing a wooden wedge cervically.
- Set the contra-angle at a value of 20'000 oscillations per minute max.
- Activate water spray during the entire treatment (minimum 50 ml/min).
- Regulate water flow in case of water excess or shortage
- Insert the strip in a cervical direction while the contra-angle is moving. This prevents the instrument from jamming.
- Operative pressure: bend the strip 1-1.5 mm to achieve the ideal 1-2.5N operating force.
- Do not stop the strip in the interdental space
- Remove the strip from the interdental space whilst the contra-angle is moving

Maintenance and sterilization

- Clean the instruments and remove debris after each use, so as to maintain their abrasive properties.
- Instruments are delivered non sterile. They must be sterilized prior to first use and disinfected, cleaned, and sterilized after each use.
- Use a single-pulsed or fractionated vacuum autoclave and subvacuum drying. Chemiclave sterilizers may also be used. Hot air sterilizers are not suitable for diamond instruments.
- Sterilization must be carried out according to validated procedures.
- Concerning the sterilization process we refer to the ISO standard 17664. We suggest to follow the indications below:
 - Cycles at 134°C:
 - Tmin = 134°C – Tmax = 138°C
 - Pressure = 3.15 bar abs
 - Time = 4 min (raisable)
 - Cycles at 121°C:
 - Tmin = 121°C – Tmax = 125°C
 - Pressure = 2.10 bar abs
 - Time = 16 min (raisable)

Risk warnings

- Avoid bending and levering over the limit: use of excessive force (strong levering effect with the contra-angle) could lead to the breakage of the instrument.
- Avoid jamming or levering actions when oscillating, as this increases the risk of instrument breakage.
- Having passed over the tooth equator line with Intensiv Ortho-Strips Opener, pay attention not to injure the soft tissues.
- Oscillating over 20'000 oscillations per minute could cause breakage of the instrument.
- Avoid temperatures above 180°C which may affect the durability of the instrument.
- Avoid applying forces greater than 3N as this could cause damage to the instrument.
- To ensure traceability of the instruments during their entire application, we recommend keeping the packaging.

Packaging of 3, 6

Intensiv Ortho-Strips Opener, Double-Sided, Ref. OS08OP
Intensiv Ortho-Strips Opener, One-Sided, Ref. OS08OP-R, right
Intensiv Ortho-Strips Opener, One-Sided, Ref. OS08OP-L, left

Additional models

Packaging of 1, 3, 6, 12
Intensiv Ortho-Strips, Ref. OS90
Intensiv Ortho-Strips, Ref. OS60
Intensiv Ortho-Strips, Ref. OS40
Intensiv Ortho-Strips, Ref. OS25
Intensiv Ortho-Strips, Ref. OS15
Intensiv Ortho-Strips Set, Ref. 060A (OS90, OS40, OS25, OS15)
Intensiv Ortho-Strips Set, Ref. 060B (OS60, OS40, OS25, OS15)

Packaging of 3, 6

Intensiv Ortho-Strips Coarse, Double-Sided, Ref. OS80C
Intensiv Ortho-Strips Coarse, One-Sided, Ref. OS80C-R, right
Intensiv Ortho-Strips Coarse, One-Sided, Ref. OS80C-L, left

Packaging of 1, 3, 6, 12

Intensiv Ortho-Strips One-Sided, Ref. OS15L, left
Intensiv Ortho-Strips One-Sided, Ref. OS15R, right
Intensiv Ortho-Strips One-Sided, Ref. OS25L, left
Intensiv Ortho-Strips One-Sided, Ref. OS25R, right
Intensiv Ortho-Strips One-Sided, Ref. OS40L, left
Intensiv Ortho-Strips One-Sided, Ref. OS40R, right
Intensiv Ortho-Strips One-Sided, Ref. OS60L, left
Intensiv Ortho-Strips One-Sided, Ref. OS60R, right
Intensiv Ortho-Strips One-Sided Set, Ref. 064 (OS15L, OS15R, OS25L, OS25R, OS40L, OS40R)

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Intensiv Ortho-Strips Opener, Double-Sided, Ref. OS08OP
Intensiv Ortho-Strips Opener, One-Sided, Ref. OS08OP-R, destra
Intensiv Ortho-Strips Opener, One-Sided, Ref. OS08OP-L, sinistra

Modelli supplementari

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Intensiv Ortho-Strips, Ref. OS90
Intensiv Ortho-Strips, Ref. OS60
Intensiv Ortho-Strips, Ref. OS40
Intensiv Ortho-Strips, Ref. OS25
Intensiv Ortho-Strips, Ref. OS15
Intensiv Ortho-Strips Set, Ref. 060A (OS90, OS40, OS25, OS15)
Intensiv Ortho-Strips Set, Ref. 060B (OS60, OS40, OS25, OS15)

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Intensiv Ortho-Strips Coarse, Double-Sided, Ref. OS80C
Intensiv Ortho-Strips Coarse, One-Sided, Ref. OS80C-R, destra
Intensiv Ortho-Strips Coarse, One-Sided, Ref. OS80C-L, sinistra

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Intensiv Ortho-Strips One-Sided, Ref. OS15L, sinistra
Intensiv Ortho-Strips One-Sided, Ref. OS15R, destra
Intensiv Ortho-Strips One-Sided, Ref. OS25L, sinistra
Intensiv Ortho-Strips One-Sided, Ref. OS25R, destra
Intensiv Ortho-Strips One-Sided, Ref. OS40L, sinistra
Intensiv Ortho-Strips One-Sided, Ref. OS40R, destra
Intensiv Ortho-Strips One-Sided, Ref. OS60L, sinistra
Intensiv Ortho-Strips One-Sided, Ref. OS60R, destra
Intensiv Ortho-Strips One-Sided Set, Ref. 064 (OS15L, OS15R, OS25L, OS25R, OS40L, OS40R)

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Intensiv Ortho-Strips Opener, Double-Sided, Ref. OS08OP
Intensiv Ortho-Strips Opener, One-Sided, Ref. OS08OP-R, derecha
Intensiv Ortho-Strips Opener, One-Sided, Ref. OS08OP-L, izquierda

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Intensiv Ortho-Strips, Ref. OS90
Intensiv Ortho-Strips, Ref. OS60
Intensiv Ortho-Strips, Ref. OS40
Intensiv Ortho-Strips, Ref. OS25
Intensiv Ortho-Strips, Ref. OS15
Intensiv Ortho-Strips Set, Ref. 060A (OS90, OS40, OS25, OS15)
Intensiv Ortho-Strips Set, Ref. 060B (OS60, OS40, OS25, OS15)

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Intensiv Ortho-Strips Coarse, Double-Sided, Ref. OS80C
Intensiv Ortho-Strips Coarse, One-Sided, Ref. OS80C-R, derecha
Intensiv Ortho-Strips Coarse, One-Sided, Ref. OS80C-L, izquierda

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Intensiv Ortho-Strips One-Sided, Ref. OS15L, izquierda
Intensiv Ortho-Strips One-Sided, Ref. OS15R, derecha
Intensiv Ortho-Strips One-Sided, Ref. OS25L, izquierda
Intensiv Ortho-Strips One-Sided, Ref. OS25R, derecha
Intensiv Ortho-Strips One-Sided, Ref. OS40L, izquierda
Intensiv Ortho-Strips One-Sided, Ref. OS40R, derecha
Intensiv Ortho-Strips One-Sided, Ref. OS60L, izquierda
Intensiv Ortho-Strips One-Sided, Ref. OS60R, derecha
Intensiv Ortho-Strips One-Sided Set, Ref. 064 (OS15L, OS15R, OS25L, OS25R, OS40L, OS40R)