- ISTRUZIONI PER L'UTILIZZO DELLE MINI VITI ORTODONTICHE CON SLOT
- INSTRUCTIONS FOR USE ORTHODONTICS MINI SCREW WITH SLOT
- INSTRUCTIONS POUR L'UTILISATION DE MINI-VIS ORTHODONTIQUES AVEC SLOT
- INSTRUCCIONES PARA EL USO DE LOS MINITORNILLOS ORTODONTICOS CON SLOT

Kristal S.r.l. via Edison, 14
20090 Trezzano S/N (MI) - Italia
Tel. +39 024453401 - info@kristalort.net

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Fornitura: Nr. 1 mini vite ortodontica sterile

- INSTRUCTIONS FOR USE OF MINI SCREW WITH SLOT

Description:
The orthodontic mini-screws are made of titanium grade 5 according to ASTM F-136 (Titanium Alloy Ti6Al4V).

The sampling system using screw inserter allows the professional to transform and position the desired orthodontic mini lives in conditions of absolute sterility, avoiding direct manipulation. The orthodontic mini screws is extracted directly from the sterile package using a manual screwdriver, to be placed directly into the surgical site.

All packages bear externally and features: description, batch number, expiry date, code reorganization, and color coding for easy recognition. The Storm orthodontic mini screws are treated for decontamination that guarantees purity surface away from potential contaminants (particulate processing). The packaging is made in the white house and sterilization is done with gamma ray irradiation (25 kGy).

Use:
The orthodontic mini screws are indicated in all cases where it is required orthodontic

mini-screws is recommended when the patient's physical and functional conditions are favorable.

Application:
The operative technique is explained in the surgical protocol. It requires a high degree of accuracy in the preparation of the site through the steps. The thermal trauma from friction of the cutters should be minimized using an appropriate cutter edge, with an appropriate number of turns properly, and adequate healing time. The removal from the package requires the installation of the ideal about the rules of asepsis.

Contraindications:
- Lack of adequate bone support
- Patient in a state of poor health
- Lack of oral hygiene and lack of cooperation by the patient

Warning:
The Storm orthodontic mini screws can only be used by dentists, orthodontists and dentists. The clinical preparation should include case study, detailed history, complete clinical and radiographic analysis. To prevent post-surgical complications is very important to carry out a thorough investigation study before starting treatment:

- Lesions of the vascular nervous bundle in the lower resulting in paralysis leading to numbness of the lip and then the need for device removal.

- Involvement of the sinus and nasal cavity of the naso-palatine nerve causing sinusitis, rhinitis due to perimplant, anesthesia and soreness of the upper lip.

To obtain the replacement of the goods when the quality is not compliant with the manufacturer, or if visibly damaged packaging, you must warn the distributor by telephone within 20 days from receipt of goods, making the return of the product in packaging original document compiling a regular transport of goods. Be careful, keep orthodontic mini screws away from light and heat in a dry place. Not allowed to re-sterilization of orthodontic mini screws in an autoclave, in this case, the manufacturer assumes no responsibility. The device is strictly disposable. Reuse of single use devices creates a potential risk of patient or user infection.

Store the Storm orthodontic mini-screws in the original package until time of placement in a dry place, away from light and heat (store at room temperature).

Mini screws should not be used beyond the expiration date on the package. Do not use if the package is opened or damaged. Do not use if there are some missing labels. The package damage could expose the medical device to the risk of contamination. Contamination of the device may lead to injury, illness, or death of patient.

Modification of the product specific packaging configuration may compromise sterility and functionality of the device. Absence of labeling and instructions at the point of use may lead to usage of product beyond the original shelf life and impair the required traceability.

It should be affixed on the patient's label included showing the main features of the screw inserted. We also recommend prepare sufficient documentation radiological, clinical and photographic as well as completing the informed consent signed by the patient. Known features and technical factors the manufacturer is aware of: Incorrect sterilization

Precaution:
Avoid strenuous physical activity after surgery.

Secondary effects:
Side effects of pure titanium devices are at present unknown. Despite the high degree of success of pure titanium devices, are not excluded risks of failure usually due to technical procedures or biological processes, not the product.

The potential causes of these failures can be represented by faulty surgical techniques and failure contraindications.

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