# Light curing modelling resin

Instructions for use

- DE** Gebrauchsanweisung.
- EN** Instructions for use.
- ES** Instrucciones de uso.
- FR** Mode d'emploi.
- IT** Istruzioni per l'uso.
- PT** Instruções de uso.
- PL** Instrukcja postępowania.
- CS** Návod k použití.
- SK** Návod na použitie.
- SL** Navodila za uporabo.
- ET** Kasutusjuhised.

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DE

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EN

### LIGHT CURING MODELLING RESIN PASTE

Plastic burning without residue, curable by light, suitable for modelling

Curable by light, modelling plastic, ready for use in the form of a gel for multiple indications in dentistry, specifically for pouring mold technology. The application of the workable gel is done from the syringe, the ready to use paste is easy to adapt and customize by its plastic, moldable consistency. This product can be cured in any standard light-curing devices. After the final curing, it is easy to grind both with the micromotor, as well as with milling machines. This product burns without residue during preheating of the mold.

#### APPLICATION AREAS

##### Fixed prosthesis

- Creating impressions for casting crowns and bridges, ceramic frameworks, inlays and onlays.
- Construction of frameworks for copy milling and bridge frameworks in electroplating.

##### Combination technology

- Primary constructions and secondary parts of the telescopic and conical crown technology, individual attachments, grooves/shoulder-pin secondary parts, thrust distributors, turns with interlocks, individual links and bridges etc.

##### Implantology

- Design of impressions for casting implant abutments, implant bridges and implant superstructures and tertiary frameworks. Creation of transmission and fixing aids and the blocking or fixing of single implant parts or segments.

##### Transmission and bite checking

- Transmission bars and key, or insertion aids, bite bars to check the occlusion, stop and checking bites.

#### PROCESSING

##### 1. Preparation

For isolating the gypsum surfaces with conventional alginates. Vaseline is recommended for polyurethane modelling materials. No further isolation is fundamentally required on polished metal or ceramic surfaces. Optionally apply gossamer-thin Vaseline or isolate by a thin oil film.

##### 2. Application

Remove the cap of the syringe. Put application tips on the cartridge and lock. Easily push the plunger for dispensing and modeling. On movement of the tip in the material, it flows easily and can be used, on stopping of the movement, it is firm. In the modeling of large objects with high wall thickness, the material is repeatedly applied in thin layers of 2 mm thickness until it reaches the desired thickness. We recommend carrying out intermediate polymerization for 20 seconds with a hand-held light or for 60 seconds with laboratory device before applying the next layer. Thus incomplete curing or excessive heat or stress in the objects can be avoided. After use, reclose the syringe with the cap. For objects with large volume, it is advisable to coat the impression with a thin layer of casting wax.

##### 3. Light curing

This product can be polymerized in all commercial light-curing devices with a light spectrum of 320 nm -440 nm. The average polymerization time in light curing devices using UV / UV-A as a means of illumination are circa 5 min, for Halogen lamps about 3 minutes, strobe lamps circa 1 min 30 sec and with LED - lamps about 20-40 seconds. Please note that the polymerization time is always dependent on the intensity of illumination from light unit. Consistent curing can only be ensured through regularly maintained equipment. Please pay attention to the instructions of the manufacturer.

##### 4. Important processing notes

- To achieve the best possible casting surfaces, the oxygen inhibition layer on the surface of impression must be completely removed before embedding with isopropanol.
- In milling and polishing, only work at work stations with suction.
- Avoid contact with eyes and mucous membranes. On accidental contact, rinse with running water. Uncured material can, on contact with skin and mucous membrane, cause irritation and sensitization.
- Seal completely after use.
- The product was developed for use in the dental field and may only be used for this purpose by qualified personnel.

##### 5. Technical data

Volumes: 3.0 g injections  
Dosing: Direct application by application tips  
Product colour: blue / red / transparent-clear  
Light-cure: light spectrum 320 nm- 440 nm  
Storage: protect from light

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**PROCLINIC, S.A.U.**
**MEGADENTA Dentalprodukte GmbH**