DE 4DESIGN Dental-Polymer-Block. Gebrauchsanweisung

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EN 4DESIGN dental polymer blank. Instructions for use

Product description

Acrylic polymer block especially suitable for the fabrication of CAD/CAM milled night guards and occlusal splints. Translucent and transparent (clear) colours. 4DESIGN blanks are characterised by their colourfastness, hardness and abrasion resistance.

Composition

Acrylic copolymer based on PMMA > 99% / Inert mineral pigments and other non-toxic substances <1%

Intended use

Night guards and occlusal splints.

Intended users

Dental health professionals, dental technicians and dentists, exclusively.

Indications

Milling of dental structures for therapeutic purposes as occlusal splints or night guards.

Contraindications

Use for fixed and permanent restorations. Permanent implant-supported restorations. Dental splints should not be made with 4DESIGN in case of a known allergic reaction to any of its chemical components.

Instructions for use

Can be used with all open CAD/CAM milling systems (98.5 mm discs) and many special systems (other discs and blocks according to dimensions and contours). Fix the block on the milling machine in the usual way. Carbide (tungsten) cutters can be used for machining. The milled structures can be separated using milling cutters or resin cut-off wheels.

Finishing and polishing: The milled structure can be finished and polished using the usual method for PMMA materials.

Design specifications:

Dental frameworks are designed in CAD software. A 5-axis milling machine is required to reach undercut areas. Minimum wall thickness: 1.5 mm for occlusal sides and 0.9 mm for peripheral sides.

4DESIGN transparent (100% PMMA) can be used for the fabrication of castable dental frameworks and as a test and verification material prior to the fabrication of final final frameworks and for the milling of auxiliary structures such as custom-made surgical guides.

Precautions

The product is not resistant to oxidising chemicals and strong acids. During finishing use Hoover and/or dust mask to avoid breathing dust. The material is not visible with X-rays. Do not reuse the product and do not use expired materials. Do not vary the working parameters established by the manufacturer (mechanical modifications of the material could occur, with consequences for the final manufactured product). Do not contaminate the product during the processing phase. For posterior dental bridges in PMMA, maximum two intermediates.

Interactions

Potential interactions with other products used in the treatment are highly unlikely because the product is a fully polymerised acrylic resin-based material. Before using the product, the dentist should consider known interactions and cross-reactions of the product with other materials already present in the patient's oral environment to ensure that the restoration can be placed safely and durably.

Side effects

In individual cases, rashes, sensitivity or allergies may occur when using the acrylic resin-based splints, in which case do not use it and consult a doctor.

Risks of reuse

The device can be reused as long as there is space available to mill the structures. Clean the device after each use and store it properly identified in its box. Caution; reuse of used dentures and other dental structures on others may lead to infections and transmission of diseases.

Storage conditions The product does not require any special conditions.

Disposal of the product

The machining residues of the blanks do not require specific disposal and can be disposed of in household waste. The accompanying packaging elements such as boxes should be disposed of in the appropriate containers (cardboard). Contaminated products with biological and/or infectious hazards should be disposed of in accordance with the applicable regulations.

Information for the patient

The practitioner shall inform the end-user of the following instructions:

a) Maintain normal oral hygiene and adopt a daily hygiene protocol*, avoid excessive brushing which may cause abrasive wear, and do not immerse the splint in alkalis, acids or other substances that may damage the acrylic materials.

b) The material is not visible to X-rays.

c) The material is relatively tough, but strong bites should be avoided.

d) In case of allergic reactions to the material, consult a doctor.

e) Of the following residual risks: Allergies and irritations caused by the material. Breakage and possible ingestion of the product. Deterioration of the product. Incorrect use of the product. Inadequate hygiene and maintenance of the product.

*Hygiene protocol: patients should brush their splints at least once a day. Brush with a soft bristle brush and a non-abrasive denture cleaner to remove food, plaque and other deposits.

Manufacturer's responsibility

Material for dental use: The manufacturer is not responsible for any damage that may result from use other than in accordance with the instructions for use. Class IIa medical device: Any serious incident involving the device must be reported to the manufacturer and to the competent authority of the Member State in which the incident occurred.

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