**Intended purpose**

Access to the oral cavity of the patient

**Patient target group**

- Patients with permanent teeth
- Patients with deciduous teeth
- Edentulous patients

**Intended users**

- Dentists (clinical procedure)
- Dental assistants (assistance during the clinical procedure)
- Dental hygienists (tooth cleaning procedure in the dental practice)

**Special training**

No further special training is needed.

**Use**

For dental use only.

**Description**

OptraGate® 2 is an auxiliary aid for clinical use which facilitates access to the working field. It provides retraction of lips and cheeks during dental treatment, offers the operator increased visibility and accessibility and facilitates saliva and moisture control in the oral cavity. The device, which features three-dimensional flexibility and elasticity, ensures patient comfort and assists in keeping the patient's mouth open. It is available in sizes Regular and Small.

**Indications**

None

**Areas of application:**

The product is designed to retract the lips and cheeks during diagnostic, preventive, therapeutic and orthodontic dental procedures, e.g.

- initial examination;
- professional teeth cleaning;
- fissure sealing;
- tooth whitening;
- periodontal treatment;
- restorative procedures;
- preparation;
- impression taking;
- cementation of restorations and orthodontic retention elements.

**Contraindications**

The use of the product is contraindicated if the patient is known to be allergic to any of its ingredients.

## OptraGate® 2

[en]	<b>Instructions for Use</b> Flexible aid for full mouth access	[fi]	<b>Käyttöohjeet</b> Joustava aputäline helpottamaan toimenpiteisiin pääsyä suussa
[de]	<b>Gebrauchsinformation</b> Flexibles Hilfsmittel für freien Zugang zum Mundraum	[no]	<b>Bruksanvisning</b> Flexibelt hjelpemiddel for fri tilgang til munnhulen
[fr]	<b>Mode d'emploi</b> Ecarteur flexible permettant un accès facile en bouche	[nl]	<b>Gebruiksaanwijzing</b> Flexibel hulpmiddel voor vrije toegang tot de mondholte
[it]	<b>Istruzioni d'uso</b> Retratore flessibile per l'accesso al cavo orale	[el]	<b>Οδηγίες Χρήσεως</b> Ελαστικό βοηθήμα για πλήρη πρόσβαση στο στόμα
[es]	<b>Instrucciones de uso</b> Ayuda flexible para el acceso completo a la cavidad bucal	[tr]	<b>Kullanma Talimatı</b> Ağız ortamına genişli kolaylaştıran esnek yardımcıdır
[pt]	<b>Instruções de Uso</b> Instrumento flexível para acesso completo da área bucal	[ru]	<b>Инструкция по применению</b> Эластичное вспомогательное средство для свободного доступа к полости рта пациента
[sv]	<b>Bruksanvisning</b> Flexibelt hjälpmedel för fritt arbetsfält	[pl]	<b>Instrukcja stosowania</b>

### Limitations of use

- OptraGate 2 is not sterile and thus not suitable for use in sterile oral surgical or implantological procedures.
- OptraGate 2 is not suitable for use in treatment procedures requiring absolute isolation.
-  For single use only.
- The product is not intended to be reprocessed or reused. The material is not suitable for sterilization. Cleaning and disinfecting agents may have a detrimental effect on the material and lead to undesirable reactions in the patient.

### Side effects

There are no known side effects to date.

### Interactions

There are no known interactions to date.

### Clinical benefit

Easier access to the treatment field.

### Composition

Thermoplastic elastomer (SEBS) and polypropylene.

OptraGate 2 is latex-free.

## 2 Application

The elastic material of OptraGate 2 is stretched by means of two rings. A third ring located on the elastic component supports its anatomic shape. The thicker ring of the device is positioned within the gingivobuccal fold, while the thinner ring, which features two tabs, remains outside the mouth. The elastic component between the two rings thus embraces the lips of the patient and provides retraction due to the restoring force of the rings (see Figure 1).

*Fig 1: OptraGate 2 in place with the tabs pointing downward. It does not interfere with the patient's occlusion.*

1. Select the appropriate size of the OptraGate 2. The selection of the appropriate size should be guided by the distance between the corners of the mouth when the mouth is in a relaxed position (i.e. with the lips closed or slightly open). In order to obtain sufficient retraction of the oral soft tissues, the diameter of the inner ring should ideally be approx. 2 cm wider than the mouth corner distance.

Therefore the following guidelines apply:

- Size Small: Mouth corner distance = 6 cm
- Size Regular: Mouth corner distance = 7 cm

The orientation scale with the respective guide marks imprinted on the foil bag can be used as an aid in the selection of the correct size. It is important not to use a rigid orientation aid for the determination of the mouth corner distance, but to follow the slightly curved course of the lips. In case of doubt, we recommend using the larger size in order to ensure a better fit. The size (R for Regular or S for Small) is indicated on the packaging as well as on the left tab of the OptraGate 2. In order to ensure optimum function, we recommend inserting the largest possible size that fits the patient's mouth.

2. OptraGate 2 may be inserted from both the 12 o'clock position with the patient lying down, or from the 8 o'clock position with the patient sitting upright (or, the 4 o'clock position, if the operator is left-handed).
3. To ease placement, with the tabs of the OptraGate 2 pointing downward, the operator should grasp the thicker intraoral ring by reaching through the thinner extraoral ring. The inner ring should be held between the thumb and middle finger and pressed together slightly (see Figure 2).

*Fig 2: Correct hold on the intraoral ring with thumb and middle finger. The intraoral ring is slightly compressed in the process.*

4. The slightly compressed intraoral ring is then inserted into the buccal corridor on one side, so that the elastic component embraces the corner of the mouth and the outer ring is still located extraorally (see Figure 3).

*Fig. 3: The thicker intraoral ring is positioned between the teeth and the corner of the mouth.*

- As soon as OptraGate 2 is securely in place on one side, the other side is inserted in the same manner by slightly bending it (see Figure 4).

Fig 4: The intraoral ring is positioned behind the two corners of the mouth.

- Subsequently, the intraoral ring is placed behind the lower and upper lips and thus OptraGate 2 assumes its final, stable position (see Figure 5). The intraoral ring can be inserted with more ease if the patient's mouth is relaxed.

Fig 5: Final positioning behind the lower and upper lip.

If, in individual cases, the intraoral ring should show a tendency to slip out of the gingivobuccal fold when the mouth is closed completely, positioning of the intraoral ring deeper inside the vestibule is usually sufficient to solve the problem. Sometimes, the use of another size may also be a solution. The securely positioned OptraGate 2 improves the overall view and accessibility of the working field. Lateral movements of the mandible and opening of the mouth to various degrees may provide additional space and access to the treatment area. The OptraGate 2 may even be left in place when checking the occlusion.

- To remove OptraGate 2, grasp the upper area of the extraoral ring with the help of a paper towel and pull it slightly downwards so that the intraoral ring is loosened from the upper gingivobuccal fold (Figure 6). OptraGate 2 can then be easily removed from the lower gingivobuccal fold and hygienically disposed of using the paper towel.

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Fig 6: Removal of the OptraGate 2

### 3 Safety information

- In the case of serious incidents related to the product, please contact Ivoclar Vivadent AG, Bendererstrasse 2, 9494 Schaan/Liechtenstein, website: [www.ivoclar.com](http://www.ivoclar.com), and your responsible competent authority.
- The current Instructions for Use and explanation of the symbols are available on the website: [www.ivoclar.com](http://www.ivoclar.com).

#### Warnings

- Observe the Safety Data Sheet (SDS) (available at [www.ivoclar.com](http://www.ivoclar.com)).
- During bleaching procedures, care should be taken to prevent any bleaching agent from getting trapped between the OptraGate 2 and the soft tissue in order to prevent irritation of the skin and mucous membrane.

#### Disposal information

Remaining stocks must be disposed of according to the corresponding national legal requirements.

#### Residual risks

Users should be aware that any dental intervention in the oral cavity involves certain risks.

The following clinical residual risks are known:

- Pressure sores

### 4 Shelf life and storage

- Storage temperature 2-28 °C (36-82 °F)
- Keep away from sunlight.
- Do not use the product after the indicated expiry date.
- Expiry date: See note on packaging.
- Before use, visually inspect the packaging and the product for damage. In case of any doubts, please contact Ivoclar Vivadent AG or your local dealer.

### 5 Additional information

Keep material out of the reach of children!

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the materials for their suitability and use for any purpose not explicitly stated in the Instructions.

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