Prior to use, carefully read the instructions for use.

EN

## DENTIN CONDITIONER

POLYACRYLIC ACID CAVITY CONDITIONER

For use only by a dental professional in the recommended indications.

### RECOMMENDED INDICATIONS

DENTIN CONDITIONER is a mild polyacrylic acid solution designed to remove the dentinal smear layer and to condition dentine, thus enhancing the bond between glass ionomer cement and the dentine.

### CONTRAINDICATIONS

In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

### DIRECTIONS FOR USE

- After tooth preparation, apply DENTIN CONDITIONER to the bonding surfaces for 20 seconds using a cotton pellet or sponge.
- Rinse thoroughly with water. Dry by blotting with a cotton pellet or gently blowing with an air syringe. DO NOT DESICCATE. Best results are obtained when prepared surfaces appear moist (glistening).
- Apply GC Fuji BOND LC, GC Fuji II LC or other glass ionomer as a bonding agent, base/liner or restorative material.
- Close the bottle immediately after use.

### STORAGE

Store at temperature of 4-25°C (39.2-77.0°F). (Shelf life : 3 years from date of manufacture).

### PACKAGE

Bottle of 25g (23.8mL).

### CAUTION

- In case of contact with oral tissue or skin, flush immediately with water.
- In case of contact with eyes, flush immediately with water and seek medical attention.

Last revised : 05/2010

**CE** 0086

**'GC.'**

MANUFACTURED by  
**GC CORPORATION**  
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

DISTRIBUTED by  
**GC CORPORATION**  
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

**EU: GC EUROPE N.V.**  
Researchpark Haasrode-Leuven 1240, Interleuvenlaan 33,  
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

**GC AMERICA INC.**  
3737 West 127th Street, Alsip, IL 60803 U.S.A.  
TEL: +1-708-597-0900

**GC ASIA DENTAL PTE. LTD.**  
19 Loyang Way, #06-27 Singapore 508724  
TEL: +65 6546 7588

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