



FICHA TÉCNICA

MASCARILLAS RECTANGULARES BESTDENT		
REFERENCIAS	<p>01502: MASCARILLAS VERDES BESTDENT (50u.)</p> <p>01503: MASCARILLAS AZULES BESTDENT (50u.)</p> <p>15021: MASCARILLAS NEGRAS BESTDENT (50u.)</p> <p>15022: MASCARILLAS BLANCAS BESTDENT (50u.)</p> <p>15023: MASCARILLAS ROSAS BESTDENT (50u.)</p> <p>15024: MASCARILLAS LILAS BESTDENT (50u.)</p>	
FABRICANTE / RESPONSABLE	<p>Hubei V-Medical Products Co., Ltd. No. 67 Xianhong Road West, Xinlirenkou, Zhanggou Town, Xiantao City, 433012 Hubei Province, P.R.China</p>	
CLASIFICACIÓN	Producto Sanitario Clase I, conforme al Reglamento 2017/745/UE	
DESCRIPCIÓN	Mascarillas rectangulares de 3 capas, con elásticos. Tipo IIR, en color verde, azul, negra, blanca rosa o lila. No estériles.	
COMPOSICIÓN	Capa externa	Polipropileno
	Filtro	Polipropileno
	Capa interna	Polipropileno
	Elásticos	Poliéster
ESPECIFICACIONES	<p>Mascarilla tipo IIR</p> <p>Eficacia de la filtración bacteriana (EFB): >98%</p> <p>Presión diferencial: <60 Pa/cm²</p> <p>Presión de resistencia a las salpicaduras: ≥16,0 kPa</p> <p>Medidas: 17,5 x 9,5 cm</p>	
ENVASE	50 mascarillas por envase	
VIDA ÚTIL	3 años a partir de la fecha de fabricación	
FINALIDAD PREVISTA	Se utilizan como barrera para las emisiones respiratorias de las personas (saliva, secreciones...), permitiéndoles respirar sin dificultad. Especialmente diseñadas para situaciones médicas, salas limpias, visitas... Pueden proteger y prevenir de infecciones bacterianas de una manera eficaz.	
INFORMACIÓN SOBRE LA SEGURIDAD	Producto de un solo uso. No usar si el envase está dañado.	
ALMACENAMIENTO Y TRANSPORTE	Almacenar en lugar fresco y seco. Mantener alejado de la luz solar. Proteger de la humedad y evitar el calor excesivo.	
ESTÁNDARES QUE CUMPLEN	EN 14683	

TECHNICAL DATA SHEET

BESTDENT RECTANGULAR FACE MASKS		
REFERENCES	<p>01502: BESTDENT GREEN FACE MASKS (50u.)</p> <p>01503: BESTDENT BLUE FACE MASKS (50u.)</p> <p>15021: BESTDENT BLACK FACE MASKS (50u.)</p> <p>15022: BESTDENT WHITE FACE MASKS (50u.)</p> <p>15023: BESTDENT PINK FACE MASKS (50u.)</p> <p>15024: BESTDENT PURPLE FACE MASKS (50u.)</p>	
MANUFACTURER	<p>Hubei V-Medical Products Co., Ltd. No. 67 Xianhong Road West, Xinlirenkou, Zhanggou Town, Xiantao City, 433012 Hubei Province, P.R.China</p>	
CLASSIFICATION	Medical Device, Class I, according to Regulation 2017/745/UE	
DESCRIPTION	Rectangular three-layered masks, with elastic bands. Type IIR, green, blue black, white, piunk or purple. Non-sterile. For single use.	
COMPOSITION	Outer layer	Polypropylene
	Filter	Polypropylene
	Inner layer	Polypropylene
	Elastic bands	Polyester
SPECIFICATIONS	<p>Type IIR face mask</p> <p>Bacterial Filtration Efficiency (BFE): >98%</p> <p>Differential pressure : <60 Pa/cm²</p> <p>Splash resistance: ≥16,0 kPa</p> <p>Size: 17,5 x 9,5 cm</p>	
PACKAGING	50 masks per box	
SHELF LIFE	3 years from the date of manufacture	
INTENDED USE	Intended to trap the emanations of the respiratory flows (saliva, secretions...) of users, allowing them to breathe without difficulty. They are especially designed for medical situations, clean rooms, and visitors, and can efficiently protect and prevent from bacterial infection.	
SAFETY INFORMATION	For single use only. Do not use if package is damaged.	
HANDLING AND STORAGE CONDITIONS	Store in a cool dry place, away from direct sunlight. Protect from moisture.	
STANDARDS APPLIED	EN 14683	