



FICHA TÉCNICA

PLANCHAS MOUTHGUARD BESTDENT		
REFERENCIAS	21927: 040MB: 1,0 mm 21928: 060MB: 1,5 mm 21929: 120MB: 3,0 mm	
FABRICANTE / RESPONSABLE	J. Ripoll, S.L. Río Tiétar 20, Pol Ind. El Nogal 28110 Algete Madrid (Spain)	
CLASIFICACIÓN	Producto Sanitario Clase I, conforme al Reglamento (UE) 2017/745.	
DESCRIPCIÓN	Planchas de acetato de uso por profesionales dentales. Son compatibles con el uso alimentario y totalmente reciclables. No contienen ftalatos.	
COMPOSICIÓN	Formadas por polímero de etilen-vinil acetato.	
ESPECIFICACIONES	<ul style="list-style-type: none"> • Planchas Mouthguard 1,0 mm: 13x13 cm, 1,0 mm de grosor. • Planchas Mouthguard 1,5 mm: 13x13 cm, 1,5 mm de grosor. • Planchas Mouthguard 3,0 mm: 13x13 cm, 3,0 mm de grosor. 	
ENVASE	Envases conteniendo 25 unidades.	
VIDA ÚTIL	2 años desde la fecha de fabricación.	
FINALIDAD PREVISTA	Las Planchas Mouthguard son productos de uso en odontología para realizar adaptaciones a medida por parte de los laboratorios dentales	
INFORMACIÓN SOBRE LA SEGURIDAD	Solamente para uso profesional dental. Producto no estéril. No se considera un producto peligroso para la salud humana ni para el medio.	
ALMACENAMIENTO Y TRANSPORTE	Almacenar en lugar fresco y seco. Mantener alejado de la luz solar. Proteger de la humedad y evitar el calor excesivo.	

TECHNICAL DATA SHEET

BESTDENT MOUTHGUARD SHEETS		
REFERENCES	<ul style="list-style-type: none"> ▪ 040MB: 1,0 mm ▪ 060MB: 1,5 mm ▪ 120MB: 3,0 mm 	
MANUFACTURER	J. Ripoll, S.L. Río Tiétar 20, Pol Ind. El Nogal 28110 Algete Madrid (Spain)	
CLASSIFICATION	Medical Device, Class I, according to Regulation (EU) 2017/745.	
DESCRIPTION	Acetate sheets for professional dental use only. They are compatible with food contact, totally recyclable and do not contain phthalates.	
COMPOSITION	Made of ethylene-vinyl acetate polymer (EVA).	
SPECIFICATIONS	<ul style="list-style-type: none"> • Mouthguard Sheets 1,0 mm: size: 13x13 cm; thickness: 1,0 mm • Mouthguard Sheets 1,5 mm: size: 13x13 cm, thickness: 1,5 mm • Mouthguard Sheets 2,0 mm: size: 13x13 cm, thickness: 2,0 mm 	
PACKAGING	Packages containing 25 pieces.	
SHELF LIFE	2 years from the date of manufacture.	
INTENDED USE	Mouthguard Sheets are acetate sheets for dental use for the preparation of tailor-made adaptations by dental laboratories.	
SAFETY INFORMATION	For professional dental use only. Non-sterile product. It is not considered a dangerous product for health or the environment.	
HANDLING AND STORAGE CONDITIONS	Store in a cool dry place, away from direct sunlight. Protect from moisture, and avoid excessive heat.	