



	<b>TECHNICAL SHEET</b>	<b>DATE</b>
		<b>14-05-18</b>

**DISPOSABLE DENTAL TRAY 28x18 COLORED - ROUGH BOTTOM**

<b>REFERENCE</b>	BD4182.XX	<b>CLIENTE:</b>	
<b>PACKAGING:</b>	CONVENCIONAL	<b>Customer ref:</b>	
<b>DRAWING:</b>	<input checked="" type="checkbox"/> ADJUNTO		



**MATERIAL :**

P.E.T.

**THICKNESS :** 180 +/- 5% MICRONS

**COLOR :**

DARK BLUE, CYAN, ORANGE, YELLOW

**AVERAGE UNIT WEIGH:** 14,50 grs.      **TECHNIQUE: :** THERMOFORMED


**PACKAGING:**

<b>EMBALAJE :</b>	<b>LOGO :</b>	PROPIO CLIENTE	<b>DISTRIBUTION :</b>
	<b>BOX MODEL:</b>	14 click coin blanca	
	<b>SIZE:</b>	400X300X110 mm	
<b>QUANTITY x BOX :</b>	<b>400 UNITS</b>		2 piles of 200 units each, in a unique polyethylene bag

**PALETIZATION:**

<b>PALLET:</b>	EUROPALET 800 X 1200
<b>MAXIMUM HEIGH:</b>	1.780,00 mm
<b>DISTRIBUTION:</b>	18 boxes x 5 stages + 6 boxes
<b>QUANTITY x PALLET:</b>	<b>38.400 UNITS</b>
<b>FINAL WRAPPING:</b>	Shrink-Wrapping
<b>Packaging observations:</b>	

**REGULATION:**  
 Medical device class 1, European directive 93/42 CEE

QUALITY DEPT.	CLICK COIN, S.L.	CUSTOMER COMPLIANCE	
<b>QUALITY MANAGER</b>	Gemma Bochaca Royo	<b>RESPONSIBLE:</b>	
<b>DATE:</b>	14-05-18	<b>DATE:</b>	
<b>SIGNATURE:</b>		<b>SIGNATURE:</b>	(If needed)