



FICHA TÉCNICA

RETRACTOR DE LABIOS BESTDENT		
REFERENCIAS	25654: Retractor de Labios Adulto+Niño	
FABRICANTE / RESPONSABLE	Huanghua Promisee Dental Co.Ltd Mould Industry Garden, Empolder Section, Huanghua 061100, Hebei, China	
CLASIFICACIÓN	Producto Sanitario Clase I, conforme al Reglamento (UE) 2017/745.	
DESCRIPCIÓN	El Retractor de Labios es un instrumento de uso odontológico que se utiliza para mantener abierta la boca y proporcionar una buena visión de los dientes. No estéril.	
ESPECIFICACIONES	Fabricados en PC. Autoclavable. Color transparente. Tamaños: - 167x55x24 - 145x45x20	
EMBALAJE	Envases con 2 unidades, una de adultos y otra de niños.	
CADUCIDAD	3 años a partir de la fecha de fabricación.	
FINALIDAD PREVISTA	Utilizados para desplazar los tejidos blandos para mejorar la visión y el acceso a la boca, y permitir protección durante procedimientos odontológicos.	
INFORMACIÓN SOBRE LA SEGURIDAD	Solamente para uso profesional dental. Necesitan esterilizarse antes de su uso.	
ALMACENAMIENTO Y TRANSPORTE	Almacenar en lugar fresco y seco. Mantener alejado de la luz solar. Proteger de la humedad y evitar el calor excesivo.	
MSDS	No aplicable	

TECHNICAL DATA SHEET

LIPS RETRACTOR BESTDENT		
REFERENCES	25654: Lips Retractor, Adult + Children	
MANUFACTURER	Huanghua Promisee Dental Co.Ltd Mould Industry Garden, Empolder Section, Huanghua 061100, Hebei, China	
CLASSIFICATION	Medical Device, Class I, according to Regulation (EU) 2017/745.	
DESCRIPTION	Lips Retractor is a hand-held dental instrument used to open the mouth and provide good eyesight of the teeth. Non-sterile.	
SPECIFICATIONS	Made of PC. Transparent and autoclavable. Sizes: - 167x55x24 - 145x45x20.	
PACKAGING	Boxes with 2 units.	
SHELF LIFE	3 years from the date of manufacture.	
INTENDED USE	Used to displace soft tissues to improve vision and access and to afford them protection during surgical procedures.	
SAFETY INFORMATION	For professional dental application only. Non-sterile product. Need to be sterilized before use.	
HANDLING AND STORAGE CONDITIONS	Store in a cool dry place, away from direct sunlight. Protect from moisture and avoid excessive heat.	
MSDS	Not Applicable	