



FICHA TÉCNICA

PINZA PARA DIQUE DE GOMA BESDENT		
REFERENCIAS	Ver Anexo 1	
FABRICANTE / RESPONSABLE	ASA DENTAL SPA Via Valenzana, 60 55054 MASSAROSA – BOZZANO LUCCA - ITALY (IT)	
CLASIFICACIÓN	Producto Sanitario Clase I, según el Reglamento (UE) 2017/745 relativo a productos sanitarios.	
DESCRIPCIÓN	Pinza de acero inoxidable para dique de goma. Instrumento manual no estéril de alta resistencia y elasticidad.	
COMPOSICIÓN	Acero inoxidable adecuado para instrumentos ortopédicos, quirúrgicos y dentales, según definido en el standard UNI EN ISO 7153-1 – <i>Instrumentos Quirúrgicos. Materiales. Parte 1: Metales</i> . Contiene ≤1% de Níquel.	
EMBALAJE	Envase con 1 unidad.	
FINALIDAD PREVISTA	Pinza utilizada por los dentistas como ayuda para asegurar el dique de goma en los dientes.	
LIMPIEZA Y ESTERILIZACIÓN	Todos los instrumentos se deben examinar, limpiar y esterilizar antes de cada uso. Esterilizar en autoclave o con calor seco.	
ALMACENAMIENTO Y TRANSPORTE	Almacenar en lugar fresco y seco. Mantener alejado de la luz solar directa y de fuentes de calor. Proteger de la humedad.	
MSDS	No aplicable	

TECHNICAL DATA SHEET

RUBBER DAM CLAMPS BESDENT		
REFERENCES	See Annex 1	
MANUFACTURER	ASA DENTAL SPA Via Valenzana, 60 55054 MASSAROSA – BOZZANO LUCCA - ITALY (IT)	
CLASSIFICATION	Medical Device class I, according to Regulation (EU) 2017/745 concerning medical devices.	
DESCRIPTION	Stainless steel rubber dam clamp. Non sterile hand-operated device with high strength and elasticity.	
COMPOSITION	Stainless steel adequate for instruments used in general surgery, orthopaedics and dentistry, according to the standard UNI EN ISO 7153-1 – <i>Surgical Instruments. Materials. Part 1: Metals</i> . Contains $\leq 1\%$ Nickel.	
PACKAGING	Individually packed.	
INTENDED USE	Clamp used by dentists as an aid for securing a so-called rubber dam sheet on the teeth.	
STERILISATION	The device must be inspected, cleaned, and sterilized before each use. Methods of sterilisation: autoclave or dry heat.	
HANDLING AND STORAGE CONDITIONS	Store in a cool dry place, away from direct sunlight. Protect from humidity and excessive heat.	
MSDS	Not applicable.	



ANEXO 1 / ANNEX 1

REF	Name
38890	CLAMP N.212 BESTDENT
38891	CLAMP 26N BESTDENT
38892	CLAMP 27N BESTDENT
38893	CLAMP N.W14 BESTDENT
38894	CLAMP N.W2 BESTDENT
38895	CLAMP N.W7 BESTDENT
38896	CLAMP N.W8 BESTDENT
38897	CLAMP N.W8A BESTDENT
38898	CLAMP ALA N.00 BESTDENT
38899	CLAMP ALA N.1 BESTDENT
38900	CLAMP ALA N.12A BESTDENT
38901	CLAMP ALA N.13A BESTDENT
38902	CLAMP ALA N.14 BESTDENT
38903	CLAMP ALA N.14A BESTDENT
38904	CLAMP ALA N.2 BESTDENT
38905	CLAMP ALA N.2A BESTDENT
38906	CLAMP ALA N.4 BESTDENT
38907	CLAMP ALA N.7 BESTDENT
38908	CLAMP ALA N.8 BESTDENT
38909	CLAMP ALA N.8A BESTDENT
38910	CLAMP ALA N.9 BESTDENT