



FICHA TÉCNICA

SUTURAS DE SEDA PROCLINIC		
REFERENCIAS	VER ANEXO 1	
FABRICANTE / RESPONSABLE	ERGON SUTRAMED S.r.l Via Giuseppe Gregoraci, 12 00173 Roma, Italia Tel. 0672677319	
DESCRIPCIÓN	Sutura de seda natural trenzada Proclinic que proporciona un excelente y continuado soporte tisular hasta 3 meses. Revestida de cera y silicona. Alta fuerza tensil, no reabsorbible. Suave con los tejidos blandos, produciendo el mínimo trauma posible. Esterilizadas, en envase hermético para conservar las condiciones de esterilidad. Color negro.	
COMPOSICIÓN	Seda natural, cera y silicona.	
PRESENTACIÓN	12 unidades de 75 cm por caja.	
FINALIDAD	Suturar o ligar tejidos.	
INSTRUCCIONES DE USO	Ver instrucciones en el producto.	
ALMACENAMIENTO Y TRANSPORTE	No se requieren condiciones especiales de almacenamiento y transporte.	

FICHA TÉCNICA

PROCLINIC SILK SUTURES		
REFERENCES	SEE ANEX 1	
MANUFACTURER	<p>ERGON SUTRAMED S.r.l Via Giuseppe Gregoraci, 12 00173 Roma, Italia Tel. 0672677319</p>	
DESCRIPTION	<p>Proclinic braided natural silk suture that provides excellent and continuous tissue support for up to 3 months. Wax and silicone coated. High tensile strength, non-resorbable. Gentle on soft tissues, producing the least possible trauma. Sterilized, in an airtight container to preserve sterility conditions. black colour.</p>	
COMPOSITION	<p>Natural silk, wax and silicone.</p>	
PACKAGING	<p>12 units of 75 cm per box</p>	
INTENDED USE	<p>Suture or ligate tissues.</p>	
INSTRUCTIONS FOR USE	<p>See instructions on the product.</p>	
HANDLING AND STORAGE CONDITIONS	<p>No special storage and transport conditions are required.</p>	

FICHA TÉCNICA

ANEXO 1 / ANEX 1:

ref.	Descripción
49934	SUTURA SEDA 2/0 -(1/2-R.-22mm) 12u. PROC
49935	SUTURA SEDA 2/0 -(1/2-T.-22mm) 12u. PROC
49936	SUTURA SEDA 3/0 -(1/2-R.-20mm) 12u. PROC
49937	SUTURA SEDA 3/0 -(1/2-T.-20mm) 12u. PROC
49938	SUTURA SEDA 3/0 -(1/2-T.-26mm) 12u. PROC
49939	SUTURA SEDA 3/0 -(3/8-T.-19mm) 12u. PROC
49940	SUTURA SEDA 4/0 -(1/2-T.-20mm) 12u. PROC
49941	SUTURA SEDA 4/0 -(1/2-T.-26mm) 12u. PROC
49942	SUTURA SEDA 4/0 -(3/8-T.-19mm) 12u. PROC
49943	SUTURA SEDA 5/0 -(3/8-T.-19mm) 12u. PROC
49968	SUTURA SEDA 2/0 -(3/8-T.-19mm) 12u. PROC
49969	SUTURA SEDA 3/0 -(1/2-R.-22mm) 12u. PROC
49993	SUTURA SEDA 3/0 -(1/2-T.-22mm) 12u. PROC
49994	SUTURA SEDA 3/0 -(3/8-T.-16mm) 12u. PROC
49995	SUTURA SEDA 4/0 -(1/2-R.-22mm) 12u. PROC
49996	SUTURA SEDA 4/0 -(1/2-T.-22mm) 12u. PROC
49997	SUTURA SEDA 4/0 -(3/8-T.-16mm) 12u. PROC
49998	SUTURA SEDA 5/0 -(1/2-T.-22mm) 12u. PROC