



FICHA TÉCNICA

ESPEJOS N.5 C.S RODIO – FRONT SURFACE		
REFERENCIAS	59959	
FABRICANTE / RESPONSABLE	ASA DENTAL SPA Via Valenzana, 60 55054 MASSAROSA – BOZZANO LUCCA - ITALY (IT)	
CLASIFICACIÓN	Producto Sanitario Clase I, bajo reglamento 2017/745 (EU) MDR.	
DESCRIPCIÓN	Espejos fabricados con una capa de rodio para evitar el problema de la doble imagen y a su vez hacerlos más resistentes al rayado y al empañado.	
ESPECIFICACIÓN	Rosca Cone Socket C.S. Nº 5 Diámetro: 24 mm.	
COMPOSICIÓN	Instrumento fabricado en acero inoxidable, espejo con acabado de rodio.	
VIDA ÚTIL	La vida útil del instrumento depende de la frecuencia de uso, el cuidado del usuario y los métodos de reprocesamiento adecuados.	
PRESENTACIÓN	12 unidades.	
FINALIDAD	Para uso por profesionales dentales para una mejor visibilidad dentro del arco dental.	
ESTERILIZACIÓN	El producto se debe esterilizar antes de su uso. Esterilizar con autoclave o con calor seco. No calentar a más de 177 °C.	
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

PROCLINIC RHODIUM MOUTH MIRRORS, No. 5, C.S.		
REFERENCES	59959	
MANUFACTURER	ASA DENTAL SPA Via Valenzana, 60 55054 MASSAROSA – BOZZANO LUCCA - ITALY (IT)	
CLASSIFICATION	Class I Medical Device, under regulation 2017/745 (EU) MDR.	
DESCRIPTION	Mirrors manufactured with a layer of rhodium to avoid the problem of double image and make them more resistant to scratching and fogging. Cone Socket C.S. No. 5.	
SPECIFICATION	Cone Socket C.S. No. 5. Diameter 24 mm.	
COMPOSITION	Instrument made of stainless steel, mirror with rhodium finish.	
SHELF LIFE	Instrument lifetime depends on frequency of use, user care and proper reprocessing methods.	
PACKAGE	12 units.	
FINALITY	For dental professional use only. Used for better visibility within the dental arch.	
STERILIZATION	The product should be sterilised prior to use. Sterilize in an autoclave or by dry heat. Do not heat above 177 °C.	
STORAGE AND TRANSPORT	It does not require special storage and transportation conditions.	