


	FICHA TÉCNICA	TDS_
	PUNTAS PAPEL CONICIDAD 0.04 y 0.06	Fecha: 14/Oct/19
	PROCLINIC	Versión: 01

PRODUCTO	PUNTAS PAPEL TAPER 0.04 y 0.06 Nº 15-40	
REFERENCIAS	81800 .04 N.15 PROCLINIC 81801 .04 N.20 PROCLINIC 81802 .04 N.25 PROCLINIC 81803 .04 N.30 PROCLINIC 81804 .04 N.35 PROCLINIC 81805 .04 N.40 PROCLINIC 81806 .04 N.15-40 PROCL. 81807 .06 N.15 PROCLINIC 81808 .06 N.20 PROCLINIC 81809 .06 N.25 PROCLINIC 81810 .06 N.30 PROCLINIC 81811 .06 N.35 PROCLINIC 81812 .06 N.40 PROCLINIC 81813 .06 N.15-40 PROCL.	
COMPAÑÍA	Pearl Dent Co. Ltd. No3-15B, street No 13, Tan Binh Industrial Park, Hochiminh City, Vietnam.	
CLASSIFICACIÓN	Producto Sanitario clase IIa, Directiva 93/42/CEE	
DESCRIPCIÓN	Las Puntas de Papel son un dispositivo hecho de papel natural altamente absorbente destinado a limpiar el conducto radicular de un diente mediante su inserción y extracción.	
ESPECIFICACIONES	Los puntos de papel tienen una forma redonda y larga y un color blanco, con una cabeza grande y una cabeza pequeña. Tipo de calafateo J & J Type Tamaño estándar ISO Taper 0.04 / Taper 0.06 X - Prototype F1, F2, F3, F4, F5 Tamaño de proceso (R25, R40, R50)	
ENVASE	Embalaje: Caja de plástico compartimentalizada. Cuantificación: La caja se divide en 6 compartimentos. Contiene: 100 puntas / caja; Dentro de la caja: 17 puntas / en 4 compartimentos; 16 puntas / en 2 compartimentos	
USO PREVISTO	Sólo para uso dental. Se usa para secar el canal radicular.	
MÉTODO DE ESTERILIZACIÓN	Producto estéril. No se garantiza la esterilidad después de haber abierto el envase. Método de esterilización: Rayos gamma. No esterilizar en autoclave o por calor.	
CONDICIONES DE ALMACENAMIENTO	Almacenar en lugar fresco y seco. Mantener alejado de la luz solar directa y de fuentes de calor. Proteger de la humedad.	

	TECHNICAL DATA SHEET	TDS_
	PAPER POINTS TAPER 0.04 & 0.06	Date: 14/Oct/19
	PROCLINIC	Version: 01

PRODUCT	PROCLINIC TAPER PAPER POINTS	
REFERENCES	81800 .04 N.15 PROCLINIC 81801 .04 N.20 PROCLINIC 81802 .04 N.25 PROCLINIC 81803 .04 N.30 PROCLINIC 81804 .04 N.35 PROCLINIC 81805 .04 N.40 PROCLINIC 81806 .04 N.15-40 PROCL. 81807 .06 N.15 PROCLINIC 81808 .06 N.20 PROCLINIC 81809 .06 N.25 PROCLINIC 81810 .06 N.30 PROCLINIC 81811 .06 N.35 PROCLINIC 81812 .06 N.40 PROCLINIC 81813 .06 N.15-40 PROCL.	
COMPANY	Pearl Dent Co. Ltd. No3-15B, street No 13, Tan Binh Industrial Park, Hochiminh City, Vietnam.	
CLASSIFICATION	Class IIa medical device, Directive 93/42/EEC	
DESCRIPTION	Absorbent Paper Point is a device made from highly absorbent natural paper intended to clean the root canal of a tooth by insertion and removal.	
SPECIFICATIONS	Paper Points have round long shape and white color, with a large head and a small head. Caulk Type J & J Type ISO Standard Size Taper 0.04/ Taper 0.06 X - Prototype F1, F2, F3, F4, F5 Proc size (R25, R40, R50)	
PACKAGING	Packaging: Box of PS plastic: Slide box. Quantification: Slide box (6 compartments) contain: 100 points / box; Inside of Slide box: 17 points / in 4 compartment; 16 points / in 2 compartment	
INTENDED USE	For dental use only. Used to dry the canal.	
STERILIZATION METHOD	Sterile product. Sterility is not guaranteed after opening the container. Sterilization method: Irradiated gamma ray. The device should not be sterilized in autoclave or heat.	
STORAGE CONDITIONS	Store in a cool dry area, away from direct sunlight. Protect from moisture.	