

Penn Works Owlerton Green Sheffield, S6 2BJ

Telephone: (0114) 234 4231 Fax: (0114) 231 4966

Email: <a href="mailto:ga@swann-morton.com">ga@swann-morton.com</a>
Web: <a href="mailto:www.swann-morton.com">www.swann-morton.com</a>

## STERILE CARBON SURGICAL BLADE SPECIFICATIONS

Supplied in boxes of 100

## **Blades in this range:**

6, 9, 10, 10A, 11, 11P, E/11, Sabre E/11, 12, 12D, 13, 14, 15, Sabre D/15, 15A, 15C, 15T, 16, 40 Fit Handles No.3 Graduated, 3L, B3, B3L, No.7, No.9

18, 19, 20, 21, 22, 22A, 23, 24, 25, 25A, 26, 27, 36 Fit Handles No.4 Graduated, 4L

COMPONENT	MATERIAL	
BLADE	Carbon Steel	
	ANALYSIS:	
	CARBON:	1.2% - 1.3%
	SILICON:	0.10% - 0.35%
	MANGANESE:	0.2% - 0.45%
	CHROMIUM:	0.10% - 0.40%
	SULPHUR:	0.025% max
	PHOSPHORUS:	0.035% max
TOP FOIL	Printed Aluminum with white OPP laminate	
BOTTOM FOIL	Printed Aluminum blade shape with CPP laminate	
RUST INHIBITOR	Waxed VCI impregnated kraft paper	
UNIT CONTAINER BOX	White lined chipboard with litho machine varnish	

Products are manufactured in accordance with BS 2982 and BS EN 27740 (where applicable).

Each box is individually wrapped in cellophane before despatch.

All Components are manufactured in the E.U.

All packets and boxes are marked with the current lot number and expiry dates for traceability; boxes are bar coded where applicable.

Products are sterilised by Cobalt 60 Gamma Irradiation by:-Swann-Morton (Services) Limited, Owlerton Green, Sheffield S6 2BJ