


FICHA TÉCNICA

CUÑAS DE MADERA BESTDENT		
REFERENCIAS	95065: Naranjas 95066: Blancas 95067: Verdes 95068: Amarillas 95069: Rosas 95070: Azules	
FABRICANTE / RESPONSABLE	Huanghua Promisee Dental Co.Ltd Mould Industry Garden, Empolder Section, Huanghua 061100, Hebei, China	
CLASIFICACIÓN	Producto Sanitario Clase I, conforme al Reglamento (UE) 2017/745.	
DESCRIPCIÓN	Las cuñas de madera son productos no estériles de uso dental para colocar en los espacios interdentes. Desechables.	
COMPOSICIÓN	Fabricados en madera de arce.	
ESPECIFICACIONES	Disponible en diferentes colores: naranja, blanco, verde, amarillo, rosa, turquesa.	
ENVASE	Envases con 100 Unidades.	
VIDA ÚTIL	3 años a partir de la fecha de fabricación	
FINALIDAD PREVISTA	Utilizados por su perfecta adaptación a la morfología interdental en las obturaciones.	
INFORMACIÓN SOBRE LA SEGURIDAD	Solamente para uso profesional dental. Producto de un solo uso. No estéril.	
ALMACENAMIENTO Y TRANSPORTE	Almacenar en lugar fresco y seco. Mantener alejado de la luz solar. Proteger de la humedad y evitar el calor excesivo.	

TECHNICAL DATA SHEET

BESTDENT POLY WEDGES		
REFERENCES	95065: Orange 95066: White 95067: Green 95068: Yellow 95069: Pink 95070: Blue	
MANUFACTURER	Huanghua Promisee Dental Co.Ltd Mould Industry Garden, Empolder Section, Huanghua 061100, Hebei, China	
CLASSIFICATION	Medical Device, Class I, according to Regulation (EU) 2017/745.	
DESCRIPTION	Poly wedges are wooden devices for dental use. Non-sterile. Disposable product.	
COMPOSITION	Made of maple wood.	
SPECIFICATIONS	Available in different colours: orange, white, green, yellow, blue, pink	
PACKAGING	Package containing 100 units.	
SHELF LIFE	3 years from the date of manufacture.	
INTENDED USE	Poly Wedges are used to help recover the perfect tooth shape of the adjacent face.	
SAFETY INFORMATION	For single use only. For professional dental application only. Non-sterile product.	
HANDLING AND STORAGE CONDITIONS	Store in a cool dry place, away from direct sunlight. Protect from moisture, and avoid excessive heat.	