

# Fresa para la cavidad

## E0123

LATEX



REF E 0123

### APLICACIONES

- Perforación
- penetración
- trepanado
- excavación
- retención

### CARACTERISTICAS

Fresa redonda de Carburo de Tungsteno

		N°	ISO																
		ISO 001	001	Ø ISO	004	005	006	007	008	009	010	012	014	016	018	021	023	027	
		USA			-	1/4	1/2	-	1	-	2	3	4	5	6	7	8	-	
FG 19 mm						•	•	•	•	•	•	•	•	•	•	•	•	•	•
FG 25 mm							•	•	•	•	•	•	•	•	•	•	•	•	•
RA 22,5 mm							•	•	•	•	•	•	•	•	•	•	•	•	•
HP 44,5 mm					•	•	•	•	•	•	•	•	•	•	•	•	•	•	•
				<input checked="" type="checkbox"/>	5x	5x	5x	5x	5x	5x	5x	5x	5x	5x	5x	5x	5x	5x	5x