



FICHA TÉCNICA

PINZA PARA DIQUE DE GOMA PROCLINIC		
REFERENCIAS	Ver Anexo 1	
FABRICANTE / RESPONSABLE	ASA DENTAL SPA Via Valenzana, 60 55054 MASSAROSA – BOZZANO LUCCA - ITALY (IT)	
CLASIFICACIÓN	Producto Sanitario Clase I, según el Reglamento (UE) 2017/745 relativo a productos sanitarios.	
DESCRIPCIÓN	Pinza de acero inoxidable para dique de goma. Instrumento manual no estéril de alta resistencia y elasticidad.	
COMPOSICIÓN	Acero inoxidable adecuado para instrumentos ortopédicos, quirúrgicos y dentales, según definido en el standard UNI EN ISO 7153-1 – <i>Instrumentos Quirúrgicos. Materiales. Parte 1: Metales</i> . Contiene ≤1% de Níquel.	
EMBALAJE	Envase con 1 unidad.	
FINALIDAD PREVISTA	Pinza utilizada por los dentistas como ayuda para asegurar el dique de goma en los dientes.	
LIMPIEZA Y ESTERILIZACIÓN	Todos los instrumentos se deben examinar, limpiar y esterilizar antes de cada uso. Esterilizar en autoclave o con calor seco.	
ALMACENAMIENTO Y TRANSPORTE	Almacenar en lugar fresco y seco. Mantener alejado de la luz solar directa y de fuentes de calor. Proteger de la humedad.	
MSDS	No aplicable	

TECHNICAL DATA SHEET

RUBBER DAM CLAMPS PROCLINIC		
REFERENCES	See Annex 1	
MANUFACTURER	ASA DENTAL SPA Via Valenzana, 60 55054 MASSAROSA – BOZZANO LUCCA - ITALY (IT)	
CLASSIFICATION	Medical Device class I, according to Regulation (EU) 2017/745 concerning medical devices.	
DESCRIPTION	Stainless steel rubber dam clamp. Non sterile hand-operated device with high strength and elasticity.	
COMPOSITION	Stainless steel adequate for instruments used in general surgery, orthopaedics and dentistry, according to the standard UNI EN ISO 7153-1 – <i>Surgical Instruments. Materials. Part 1: Metals</i> . Contains $\leq 1\%$ Nickel.	
PACKAGING	Individually packed.	
INTENDED USE	Clamp used by dentists as an aid for securing a so-called rubber dam sheet on the teeth.	
STERILISATION	The device must be inspected, cleaned, and sterilized before each use. Methods of sterilisation: autoclave or dry heat.	
HANDLING AND STORAGE CONDITIONS	Store in a cool dry place, away from direct sunlight. Protect from humidity and excessive heat.	
MSDS	Not applicable.	

ANEXO 1 / ANNEX 1

REF	Name
59962	CLAMP 212 PROCLINIC
59963	CLAMP ALETA 00 PROCLINIC
59964	CLAMP ALETA 1 PROCLINIC
59965	CLAMP ALETA 12A PROCLINIC
59966	CLAMP ALETA 13A PROCLINIC
59967	CLAMP ALETA 14 PROCLINIC
59968	CLAMP ALETA 14A PROCLINIC
59969	CLAMP ALETA 2 PROCLINIC
59970	CLAMP ALETA 2A PROCLINIC
59971	CLAMP ALETA 4 PROCLINIC
59972	CLAMP ALETA 7 PROCLINIC
59973	CLAMP ALETA 8 PROCLINIC
59974	CLAMP ALETA 8A PROCLINIC
59975	CLAMP ALETA 9 PROCLINIC
59976	CLAMP ALETA W14 PROCLINIC
59977	CLAMP ALETA W2 PROCLINIC
59978	CLAMP ALETA W7 PROCLINIC
59979	CLAMP ALETA W8 PROCLINIC
59980	CLAMP ALETA W8A PROCLINIC
59981	CLAMP 27 PROCLINIC
4996	CLAMP 26 PROCLINIC
38920	CLAMP B4 PROCLINIC
38911	SET 8 CLAMPS CON SOPORTE PROCLINIC