



FICHA TÉCNICA

CERA ESQUELETICOS EN PLANCHAS LISA PROCLINIC		
REFERENCIAS	H20056: 0.30 MM H20058: 0.40 MM H20060: 0.50 MM H20062: 0.60 MM	
FABRICANTE / RESPONSABLE	Al Dente Dentalprodukte GmbH Borsigstr. 1 - D-38644 Goslar +49 (0) 5321-80031 info@aldente.de Emergency Telephone Number: +49 (0) 171 8637790	
CLASIFICACIÓN	Producto clase protésico.	
DESCRIPCIÓN	Cera calibrada lisa. Tamaño: 75 x 150 mm	
COMPOSICIÓN	Parafina.	
PRESENTACIÓN	15 planchas por envase.	
FINALIDAD	Modelar las placas de esqueléticos completos y parciales.	
INSTRUCCIONES DE USO	Ver instrucciones en el producto.	
ALMACENAMIENTO Y TRANSPORTE	No se requieren condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

SKELETON WAX SMOOTH SHEETS PROCLINIC		
REFERENCES	H20056: 0.30 MM H20058: 0.40 MM H20060: 0.50 MM H20062: 0.60 MM	
MANUFACTURER	Al Dente Dentalprodukte GmbH Borsigstr. 1 - D-38644 Goslar +49 (0) 5321-80031 info@aldente.de Emergency Telephone Number: +49 (0) 171 8637790	
CLASSIFICATION	Prosthetic class product.	
DESCRIPTION	Smooth calibrated wax. Size: 75 x 150mm.	
COMPOSITION	Paraffin.	
PACKAGING	15 plates per-container.	
INTENDED USE	Model complete and partial skeletal plates.	
INSTRUCTIONS FOR USE	See instructions on the product.	
HANDLING AND STORAGE CONDITIONS	No special storage and transport conditions are required.	