



FICHA TÉCNICA

CERA ESQUELETICOS GANCHOS PROCLINIC		
REFERENCIAS	H20078: MOLARES H20080: PREMOLARES H20082: BONYHARD H20084: ANULARES RECTO	
FABRICANTE / RESPONSABLE	Al Dente Dentalprodukte GmbH Borsigstr. 1 - D-38644 Goslar +49 (0) 5321-80031 info@aldente.de Emergency Telephone Number: +49 (0) 171 8637790	
CLASIFICACIÓN	Producto clase protésico.	
DESCRIPCIÓN	Cera calcinable adaptable al modelo sin tensiones. La transparencia de la cera permite ver marcas a través de ella.	
COMPOSICIÓN	Parafina.	
PRESENTACIÓN	200 u por envase.	
FINALIDAD	Agentes auxiliares para la fabricación en la tecnología dental.	
INSTRUCCIONES DE USO	Ver instrucciones en el producto.	
ALMACENAMIENTO Y TRANSPORTE	No se requieren condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

SKELETON WAX HOOKS PROCLINIC		
REFERENCES	H20078: MOLARS H20080: PREMOLARS H20082: BONYHARD H20084: RING STRAIGHT	
MANUFACTURER	Al Dente Dentalprodukte GmbH Borsigstr. 1 - D-38644 Goslar +49 (0) 5321-80031 info@aldente.de Emergency Telephone Number: +49 (0) 171 8637790	
CLASSIFICATION	Prosthetic class product.	
DESCRIPTION	Castable wax adaptable to the model without tensions. The transparency of the wax allows you to see marks through it.	
COMPOSITION	Paraffin.	
PACKAGING	200 u per-container.	
INTENDED USE	Auxiliary agents for manufacturing in dental technology.	
INSTRUCTIONS FOR USE	See instructions on the product.	
HANDLING AND STORAGE CONDITIONS	No special storage and transport conditions are required.	