

# KeyPrint® Compatible Printers



Keystone is committed to expanding our KeyPrint Compatible portfolio to provide you with an unmatched selection of printers to select from. We believe one size printer doesn't fit all, so we extensively validate our resins on leading printer platforms, so you don't have to.

Status: March 19, 2025

## Bio Biocompatible Resins

- Validation complete
- Validation in process

	KeyModel®	KeyPrint® Model®	KeyModel™ Ultra	KeyCast®	KeyMask®	KeyGuide®	KeyTray™	KeyOrtho IBT™	KeyDenture Try-In™	KeySplint Soft®	KeySplint Hard®	KeySplint Soft® Clear	KeySplint Hard® Clear	KeyGuard™	KeyDenture® Base	
Carbon® M1/M2 Series (385 nm)					•	•			•		•	•	•	•	•	□
Carbon® M3 Series (385 nm)						•			•		•	•	•	•	•	
HeyGears A2D (385 nm)										•	•	•	•			
HeyGears A2D 4K (385 nm)										•	•					
HeyGears ChairSide Pro (385 nm)										•	•	□	•			
HeyGears Ultracraft A3D (385 nm)										•	•	•	•			
HeyGears ChairSide (385 nm)										•	•	□	□			
Asiga® MAX (385 nm)	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	□
Asiga® PRO 4k (385 nm)	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	□
Asiga® Ultra (385 nm)										•	•	•	•	•	□	□
Rapid Shape D20+, D30+, D20II, D30II,	•	•	•	•	•	•	•	•	•	•	•	•	•	•	□	
Rapid Shape D70, D90, D100 (385 nm)			•													
Sisma Everes UNO (405 nm)	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	
Microloy Eve (405 nm)	•	•	•	•	•	•	•	•	•	•	•	•	•			
Microloy Versus (385 nm)	•	•	•	•	•	•	•	•	•	•	•	•	•			
W2P SolFlex 135/150/163/170/350/	•	•	•	•	•	•	•	•	•	•	•	•	•			
W2P SolLab 335/370/450 (385 nm)	•	•	•	•	•	•	•	•	•	•	•	•	•			
MiiCraft 125 (405 nm)	•	•	•	•	•	•	•	•	•	•	•	•	•			
MiiCraft Ultra (405 nm)	•	•	•	•	•	•	•	•	•	•	•	•	•			
Dekema trix™ Print (385 nm)	•	•	•	•	•	•	•	•	•	•	•	•	•			
SprintRay Pro (405 nm)				•	•	•	•	•	•	•	•	•	•	•	□	
SprintRay Pro 2 (385 nm)										□	□	□	□	•		□
MoonRay® S100/D75 (385 nm)						•				•						
Flashforge Hunter (405 nm)	•	•		•	•	•	•			•						
Desktop Health Envision One						•		•		•						
Desktop Health Einstein (385nm)					•			•	•	•	•	•	•			□
Ackuretta FreeShape (405nm)	•	•	•	•	•	•	•	•	•	•	•	•	•			
Ackuretta Dentiq (405nm)	•	•	•	•	•	•	•	•	•	•	•	•	•			
Ackuretta SOL Mono (405nm)	•	•	•	•	•	•	•	•	•	•	•	•	•		□	
Ackuretta SOL Plus (405nm)			•							•	•	•	•			
Stratasys Origin One Dental (385 nm)				•	•	•	•	•	•	•	•	•	•	•	□	
Phrozen Sonic XL 4k (405 nm)	•	•	•	•	•	•	•	•	•	•	•	•	•	•	•	
Phrozen Sonic Mini 4k (405 nm)			•	•	•	•	•	•	•	•	•	•	•	•	•	
SHINING 3D® AccuFab-D1/D1s (405nm)	•	•	•	•	•	•	•	•	•	•	•	•	•			
SHINING 3D® Accufab-L4D (405nm)	□	□	•	•	•	•	•	•	•	•	•	•	•	•	•	□
Nexa3D NXD200 (405 nm)			•			•	•	•	•	•	•	•	•			
Nexa3D XiP (405 nm)			•		•	•	•	•	•	•	•	•	•	□		
Uniz Slash 2 (405 nm)			□							•	•					
EPAX DX10 Pro (405 nm)						•			•	•	•					
Planmeca Creo C5 (405nm)			•	□		•	•	•	•	•	•	•				
LuxCreo Lux3+ (405nm)								•		•	•					
LuxCreo iLux Pro (405nm)						•	•	•	•	•	•			•		
NewPro 3D NP1 (405 nm)										•						
RayShape® Shape 1 (405nm)			•	•	□	•	•	•	•	•	•	•	•	•	•	
Zortrax Inkspire 2 Clinical (405nm)			•			•	•	•	•	•	□	•		•		
B9 Creation Dent XL (385nm)						•	•									
Elegoo Mars 3 Pro 4k (385nm)			•	•	•	•	•	•	•	•	•	•	•	•	•	
Anycubic Photon Mono X 6Ks (405 nm)	□	□	•	□	□	□	□	□	□	•	•	•	•	•	□	
Anycubic Photon Mono M5 (405 nm)	□	□	•	□	□	•	•	□	□	•	•	•	•	•	□	□
ProgaPrint PR5 (385 nm)										•		□				
Structo ST-01 (TBD)			•							•						
Roland® 3DXPRINT (TBD)			•	•		•	•	•	•	•	•	•	•	•	•	